

NeXus-CBF Concordance Summary

Herbert J. Bernstein, Aaron S. Brewster, Daniel Paley, Tobias S. Richter, Jonathan Sloan

4 June 2025

Contents

1	Introduction	5
2	How to read this document	5
3	General Mapping Issues	6
4	Mapping from CBF to NeXus	7
4.1	The NeXus Structure Into Which to Map	7
4.2	The NeXus top level	13
4.3	Identifying images	13
4.4	The ARRAY_DATA category	14
4.5	The ARRAY_DATA_EXTERNAL_DATA category	15
4.6	The ARRAY_ELEMENT_SIZE category	15
4.7	The ARRAY_INTENSITIES category	16
4.8	The ARRAY_STRUCTURE category	17
4.9	The ARRAY_STRUCTURE_LIST category	17
4.10	The ARRAY_STRUCTURE_LIST_SECTION category	21
4.11	The ARRAY_STRUCTURE_LIST_AXIS category	22
4.12	The AXIS category	26
	4.12.1 Differences in Coordinate Frames	31
	4.12.2 Mapping Axes	32
4.13	The DIFFRN_DATA_FRAME category	32
4.14	The DIFFRN_DETECTOR category	33
4.15	The DIFFRN_DETECTOR_AXIS category	34
	4.15.1 DIFFRN_DETECTOR_ELEMENT category	34
4.16	The DIFFRN_MEASUREMENT category	34
4.17	The DIFFRN_MEASUREMENT_AXIS category	35
4.18	The DIFFRN_RADIATION category	36
4.19	The DIFFRN_REFLN category	38
	4.19.1 The MAP_SEGMENT category	44
4.20	The VARIANT category	45
5	Mapping from NeXus to CBF/imgCIF	46
5.1	Faithful Mapping of a NeXus Tree	46
5.2	Mapping Fields	47
6	NeXus Base Class Mapping	50
6.1	NXaperture	50
6.2	NXattenuator	50
6.3	NXbeam	51

6.4	NXbeam_stop	53
6.5	NXbending_magnet	54
6.6	NXcapillary	55
6.7	NXcharacterization	56
6.8	NXcollection	57
6.9	NXcollimator	57
6.10	NXcrystal	59
6.11	NXdata	62
6.12	NXdetector	64
6.13	NXdetector_group	70
6.14	NXdisk_chopper	70
6.15	NXentry	71
6.16	NXenvironment	74
6.17	NXevent_data	75
6.18	NXfermi_chopper	76
6.19	NXfilter	77
6.20	NXflipper	79
6.21	NXgeometry	79
6.22	NXguide	80
6.23	NXinsertion_device	82
6.24	NXinstrument	83
6.25	NXlog	85
6.26	NXmirror	86
6.27	NXmoderator	88
6.28	NXmonitor	89
6.29	NXmonochromator	90
6.30	NXnote	91
6.31	NXorientation	92
6.32	NXparameters	92
6.33	NXpolarizer	92
6.34	NXpositioner	93
6.35	NXprocess	94
6.36	NXroot	95
6.37	NXsample	96
6.38	NXsensor	99
6.39	NXshape	100
6.40	NXsource	101
6.41	NXsubentry	103
6.42	NXtransformations	106
6.43	NXtranslation	106
6.44	NXuser	107

6.45	NXvelocity_selector	108
6.46	NXxraylens	109
7	Proposed Pixel Array Detector Application Definitions	110
7.1	NXentry	110
7.2	NXinstrument	112
7.3	NXdetector	113
7.4	detectorSpecific:DECTRIS_detector_specific	117
7.5	detectorModule_nnn:DECTRIS_detector_module	121
7.6	detectorChip_nn:DECTRIS_detector_chip	124
7.7	Consolidated Dectris Eiger Application Definition	124

1 Introduction

This is an updated concordance of the representation of macromolecular crystallographic image data when presented in NeXus/HDF5 or CBF/imgCIF. The original version written by Herbert J. Bernstein, Tobias S. Richter, and Jonathan Sloan was released in August 2013 [3], shortly after release of Dectris Eiger detectors using NeXus/HDF5 and six years after the adoption of CBF/imgCIF for Dectris Pilatus detectors in 2007. As of this writing both image formats are still in widespread use, and conversion between them is a common part of light source workflows.

As of this writing in 2025, the primary reference ontology for coordinate data for macromolecular crystallography is the Protein Data Bank’s PDBX/mmCIF dictionary [9], and the primary reference ontology for raw crystallographic diffraction image data is the CBF/imgCIF dictionary [2] now housed in the DIALS cbflib github repository ¹. The most commonly used container other than CBF/imgCIF for macromolecular crystallographic image data is NeXus/HDF5 [5] using the Dectris Filewriter format as described in the Gold Standard paper [1]. NeXus is tree-oriented and imgCIF/CBF is table-oriented. In this document we summarize an approach to a concordance between them. Both frameworks allow for multiple alternative representations of the same data. HDF5 provides a particularly appropriate format for the management of large numbers of experimental data sets and for the representation of particularly large data sets.

The original concordance in 2013 combined input from Herbert J. Bernstein, Tobias S. Richter and Jonathan Sloan, as well as comments and suggestions by Mark Koennecke in 2010 (see <http://lists.nexusformat.org/pipermail/nexus-developers/attachments/20100201/a9424156/attachment-0001.e1>) and by the members of NIAC.

This draft is by Herbert J. Bernstein, Aaron Brester and Dan Paley, who take complete responsibility for any mistakes and misunderstandings in this version. Please send comments, corrections and suggestions to yayahjb@gmail.com.

The reader is assumed to already have familiarity with imgCIF/CBF, NeXus and HDF5.

2 How to read this document

While the mappings described here are complex and detailed, the approaches to the mappings are simple. If you are familiar with both CBF and NeXus, you should read section 4 “Mapping from CBF to NeXus” on page 7 through “Identifying images” on page 13 and section 5 “Mapping from NeXus to CBF/imgCIF” on page 46 through “Mapping Fields” on page 47. If you are interested on the handling of the Dectris Eiger detector, see section 7 “Proposed Pixel Array Detector Application Definitions” on page 110

¹<http://github.com/dials/cbflib>

3 General Mapping Issues

A CBF is organized as a set of relational tables, each table is called a “category”. The name of a category is essentially the name of a class. Categories are organized into “data blocks” the name of a category must be distinct from the name of all other categories in the same data block. As a relational table, each table has a key. One component of the key is, in general, an “ID”.

A NeXus file is organized as a tree. Each node of the tree has a name, which must be distinct from all the other names of nodes that are children of the same parent, each node of the tree is an HDF5 NeXus class instance or an HDF5 dataset or an HDF5 attribute. Each NeXus class instance has a NeXus class.

In most cases the value of the ID in each row of a CBF table will have to be mapped to the name of a NeXus class instance in the NeXus file, and in most cases the name of a NeXus class instance in a NeXus file will have to be mapped to the value of an ID in some CBF table.

To avoid namespace conflicts, in going from CBF to NeXus, we carry the CBF category name prefixed by “CBF_” along with the ID value in “dotted” notation, using a double underscore, “_”, in place of the period to conform to NeXus naming conventions, but so as to not conflict with the use of the single underscore, “_”, to separate components of names on a single hierarchical level. In going from NeXus to CBF, we strip leading dotted notation name components that match the target category name.

Similar considerations will arise throughout this mapping. The major exception is for the fields in NXdetector, in which the well-established practice in NeXus is to give those fields fixed names, such as “data”. If there is no possibility of a name conflict, we will retain the standard name and use suffixes, rather than prefixes in cases that require disambiguation. In the longer term, it would be best to ensure reliable identification of the function of fields with attributes, rather than by name. In each case, an attribute will be added for that purpose.

Most CBF categories accept identifiers for variants. NeXus does not have an equivalent concept to variants at this time.

CBF carries detailed information about the storage of images, such as the compression used, that is not exposed in NeXus. That information is carried in HDF5.

CBF is designed to organize scans, identified by a scan ID, SCANID, in the context of multiple scans per diffraction experiment, where each diffraction experiment is identified by a diffraction ID, DIFFRNID, in the context of multiple experiments per study of a macromolecule, where each macromolecular study is identified by an entry ID, ENTRYID. This creates a 3-level hierarchy of information supporting the eventual report and structural deposition. At present, NeXus, does not support an equivalent of this hierarchy. This issue should be discussed.

4 Mapping from CBF to NeXus

For the following mapping, where a placement current defined NeXus class instances under NXentry has not yet been specified the mapping is shown as “ \rightarrow ??”, but all CBF items have a mapping under an CBF_cbf NeXus class instance that we propose to place under NXentry, so further mapping could be done entirely in NeXus or HDF5, if desired. Each _CATEGORY.COLUMN value of type CBFTYPE in row NN in data block DATABLOCK is mapped to

```
/entry__ENTRYID:NXentry
  /CBF_cbf:CBF_cbf
    /DATABLOCK:CBF_cbfdb
      /CATEGORY:CBF_cbfcats
        /COLUMN:CBF_cbfcol
          /NN
            /value
```

with the attribute @cbftype=CBFTYPE, so no information will be lost, and all the information in the rest of the NeXus tree will be available in these relational tables. For most traditional NeXus uses of the file, these tables may be ignored or deleted, but they are needed for efficient external management of multiple files in the context of a larger data management system. When very large numbers of datasets have to be handled at a facility, standard practice is to put information on which searches will be done into a relational database. CBF tables are such information. Translating them faithfully into NeXus allows that information to be preserved with the NeXus/HDF5 files, rather than having to deal with two different file formats for this information. In order to index information from a NeXus file into a database, as is done, for example in the iCAT project, first the information needs to be in the NeXus file. Until the CBF to NeXus mapping has become complete and automatic, the extra CBF_cbf class ensures that all the information is available for such database use.

4.1 The NeXus Structure Into Which to Map

The following is the target tree structure. All elements of this mapping should be carefully considered and discussed. The attribute and field names with CBF category name prefixes are potential placeholder for shorter more readable names to be discussed. Names with double asterisks are flagged as particularly worth discussing.

```
/CBF_diffraction_scan__SCANID:NXentry
  /CBF_scan_id="SCANID"
  /CBF_diffraction_id="DIFFRACTIONID"
  /CBF_entry_id="ENTRYID"
  /instrument:NXinstrument
```

```

/ CBF_diffn_detector__DETECTORNAME: NXdetector
/ start_time=STARTDATETIME
/ end_time=ENDDATETIME
/ CBF_diffn_scan_frame__date=DATES
/ CBF_diffn_scan_frame__frame_id=IDS
/ average_count_time=AVGCOUNTTIME
/ @units="sec"
/ average_frame_restart_time=RSTRTIME
/ @units="sec"
/ average_frame_time=TIMEPER
/ @units="sec"
/ count_time=COUNTTIMES
/ @units="sec"
/ frame_time=TIMEPERS
/ @units="sec"
/ frame_restart_time=RSTRTIME
/ @units="sec"
/ frame_start_number=FRAMESTARTNO
/ frame_end_number=FRAMEENDNO
/ distance --> /NXentry/NXinstrument/NXsample/CBF_diffn_measurement__sample_detector_distance
/ data_ARRAYID_BINARYID=DATA
/ @CBF_array_id="ARRAYID"
/ @CBF_binary_id="BINARYID"
/ @CBF_header_contents="HEADER"
/ @CBF_header_convention="HEADERCONVENTION"
/ @x_pixel_size=XPSIZE
/ @y_pixel_size=YPSIZE
/ @CBF_array_intensities__details="DETAILS"
/ @CBF_array_intensities__gain=GAIN
/ @CBF_array_intensities__gain_esd=GAINESD
/ @CBF_array_intensities__linearity="LINEARITY"
/ @CBF_array_intensities__offset=OFFSET
/ @CBF_array_intensities__scaling=SCALING
/ @CBF_array_intensities__overload=OVERLOAD
/ @CBF_array_intensities__undefined_value=UNDEFVAL
/ @CBF_array_intensities__pixel_fast_bin_size=FBINSIZE
/ @CBF_array_intensities__pixel_slow_bin_size=SBINSIZE
/ @CBF_array_intensities__pixel_binning_method="METHOD"
/ deadtime=DTIME
/ description=DESCRIPTION
/ details=DETAILS
/ number_of_axes=NUMDETAXES
/ type=TYPE

/ CBF_array_structure_list__AXISSET1=[]
/ @CBF_array_id="ARRAYID"
/ @CBF_array_structure_list__dimension=DIM1
/ @CBF_array_structure_list__direction="DIR1"
/ @CBF_array_structure_list__index=1
/ @CBF_axis=PRECEDENCE1

```



```

/ CBF_array_structure_list__AXISSET2=[]
/ @CBF_array_id="ARRAYID"
/ CBF_array_structure_list__dimension=DIM2
/ CBF_array_structure_list__direction="DIR2"
/ CBF_array_structure_list__index=2
/ CBF_axis=PRECEDENCE2
/ CBF_array_structure_list_section__SECTIONID=[]
/ @CBF_array_id="ARRAYID"
/ CBF_array_structure_list_section__index=INDEX
/ CBF_array_structure_list_section__end=END
/ CBF_array_structure_list_section__start=START
/ CBF_array_structure_list_section__stride=STRIDE
/ CBF_array_structure_list_axis__AXISID=[] **
/ @CBF_array_structure_list_axis__axis_id="AXISID" **
/ @CBF_array_structure_list_axis__axis_set_id="AXISSETID" **
/ @CBF_array_structure_list_axis__angle=ANGLE **
/ @CBF_array_structure_list_axis__angle_increment=ANGLEINC **
/ @CBF_array_structure_list_axis__displacement=DISP **
/ @CBF_array_structure_list_axis__displacement=FRACDISP **
/ @CBF_array_structure_list_axis__fract_displacement=DISPINC **
/ @CBF_array_structure_list_axis__fract_displacement_increment=FRACINC **
/ @CBF_array_structure_list_axis__angular_pitch=ANGPITCH **
/ @CBF_array_structure_list_axis__radial_pitch=RADPITCH **
/ @CBF_array_structure_list_axis__reference_angle=REFANG **
/ @CBF_array_structure_list_axis__reference_displacement=REFDISP **
/ CBF_diffraction_scan_axis__AXISID=[]
/ @CBF_axis_id="AXISID" **
/ @CBF_diffraction_scan_axis__angle_start=ANGSTART
/ @CBF_diffraction_scan_axis__angle_range=ANGRANGE
/ @CBF_diffraction_scan_axis__angle_increment=ANGINC
/ @CBF_diffraction_scan_axis__angle_rstrt_incr=ANGRSTRT
/ @CBF_diffraction_scan_axis__displacement_start=DISPSTART
/ @CBF_diffraction_scan_axis__displacement_range=DISPRANGE
/ @CBF_diffraction_scan_axis__displacement_increment=DISPINC
/ @CBF_diffraction_scan_axis__displacement_rstrt_incr=DISPRSTRT
/ @CBF_diffraction_scan_axis__reference_angle=ANG
/ @CBF_diffraction_scan_axis__reference_displacement=DISP
/ CBF_diffraction_detector_element__id="ELEMENTID1:ELEMENTID2:..." **
/ CBF_diffraction_detector_element__reference_center_fast=[RCF1,RCF2,...] **
/ CBF_diffraction_detector_element__reference_center_slow=[RCS1,RCS2,...] **
/ CBF_diffraction_detector_element__id="UNITS1:UNITS2:..." **
/ CBF_diffraction_data_frame__section_id=SECTIONIDARRAY **
/ CBF_diffraction_data_frame__binary_id=BINARYIDARRAY **
/ CBF_diffraction_data_frame__center_fast_slow=CENTERARRAY **
/ @units="UNITS"
/ CBF_diffraction_data_frame__details=DETAILSARRAY

/ CBF_diffraction_measurement__GONIOMETER:NXsample
/ CBF_diffraction_measurement__details="DETAILS"
/ CBF_diffraction_measurement__device="DEVICE"

```

```

/ CBF_diffraction_measurement__device_details="DEVDETAILS"
/ CBF_diffraction_measurement__device_type="DEVTYPE"
/ CBF_diffraction_measurement__method="METHOD"
/ number_of_axes=NUMBER
/ CBF_diffraction_measurement__sample_detector_distance=DIST
  /@units="mm"
/ CBF_diffraction_measurement__sample_detector_voffset=VOFFSET
  /@units="mm"
/ CBF_diffraction_measurement__specimen_support="SPECSPRT"
/ CBF_diffraction_radiation__collimation="COLLIMATION"
/ divergence_x=DIVX
  /@units="deg"
/ divergence_y=DIVY
  /@units="deg"
/ CBF_diffraction_radiation__div_x_y_source=DIVXY
  /@units="deg^2"
/ CBF_diffraction_radiation__filter_edge=ABSEGE
  /@units="angstroms"
/ CBF_diffraction_radiation__inhomogeneity=HWIDTH
  /@units="mm"

/monochromator:NXmonochromator
  /wavelength=WAVELENGTH
/ CBF_diffraction_radiation__monochromator="MONOCHROMATOR"
/ CBF_diffraction_radiation__polarization_norm=POLNANG
  /@units="deg"
/ CBF_diffraction_radiation__polarization_ratio=POLRAT
/ CBF_diffraction_radiation__polarization_source_norm=POLSNANG
  /@units="deg"
/ CBF_diffraction_radiation__polarization_source_ratio=POLSRAT
/ CBF_diffraction_radiation__probe="RADIATION"
/ CBF_diffraction_radiation__type="SIEGBAHNTYPE"
/ CBF_diffraction_radiation__xray_symbol="IUPACXRAYSYMB"
/ CBF_diffraction_scan_SCANID:NXscan
/ CBF_diffraction_scan_frame_monitor__DETECTORNAME:NXmonitor
  /data=MONITORVALUES
  /count_time=INTEGRATIONTIMES
  /@units="sec"

```

In the years since the original proposal, NeXus has chosen to gather the information about the various transformations based rotation, translation, and general axes in one or more NXtransformations groups.

This is a fragmentary example of the axis definitions as proposed for NeXus taken from the CBF fragment afterwards.

```

/instrument:NXinstrument
  / CBF_diffraction_detector__DETECTOR:NXdetector
    CBF_axis__DETECTOR_PITCH=[0.]
      @units="deg"
    @CBF_location="image_1.axis.vector.10"

```

```

        @depends_on="axis__DETECTOR_Y"
        @transformation_type="rotation"
        @vector=[-1, 0, 0]
CBF_axis__DETECTOR_Y=[0.]
        @units="mm"
        @CBF_location="image_1.axis.vector.9"
        @depends_on="axis__DETECTOR_Z"
        @transformation_type="translation"
        @vector=[0, -1, 0]
CBF_axis__DETECTOR_Z=[250.]
        @units="mm"
        @CBF_location="image_1.axis.vector.8"
        depends_on="."
        @transformation_type="translation"
        @vector=[0, 0, 1]
CBF_axis__ELEMENT_X=[0.]
        @units="mm"
        @offset_units="mm"
        @CBF_location="image_1.axis.vector.11"
        @depends_on="axis__DETECTOR_PITCH"
        @transformation_type="translation"
        @vector=[-1, 0, 0]
        @offset=[-211.818, -217.322, 0]
CBF_axis__ELEMENT_Y=[0.]
        @units="mm"
        @CBF_location="image_1.axis.vector.12"
        @depends_on="axis__ELEMENT_X"
        @transformation_type="translation"
        @vector=[0, 1, 0]
/CBF_coordinate_system:NXcoordinate_system
/CBF_axis__SLS_X=[]
        @CBF_location="image_1.axis.vector.0"
        @depends_on="."
        @vector=[-1, 0, 0]
/CBF_axis__SLS_Y=[]
        @CBF_location="image_1.axis.vector.1"
        @depends_on="."
        @vector=[0, -1, 0]
/CBF_axis__SLS_Z=[]
        @CBF_location="image_1.axis.vector.2"
        @depends_on="."
        @vector=[0, 0, 1]
/CBF_axis__GRAVITY=[]
        @cbf_location="image_1.axis.vector.7"
        @depends_on="."
        @vector=[0, -1, 0]
/CBF_axis__BEAM=[]
        @cbf_location="image_1.axis.vector.6"
        @depends_on="."
        @vector=[0, 0, 1]

```

```

/CCP4_diffraction_measurement__GONIOMETER:NXsample (Note: Changed from NXgoniometer, 8Aug13)
/CCP4_axis__GONIOMETER_KAPPA=[0]
  @units="deg"
  @CCP4_location="image_1.axis.vector.4"
  @depends_on="axis__GONIOMETER_OMEGA"
  @transformation_type="rotation"
  @vector= [-0.64279, 0.76604, 0]
/CCP4_axis__GONIOMETER_OMEGA=[0]
  @units="deg"
  @CCP4_location=image_1.axis.vector.3
  @depends_on="."
  @transformation_type="rotation"
  @vector= [1, 0, 0]
/CCP4_axis__GONIOMETER_PHI=[0]
  @units="deg"
  @CCP4_location=image_1.axis.vector.5
  @depends_on="axis__GONIOMETER_KAPPA"
  @transformation_type="rotation"
  @vector= [1, 0, 0]

loop_
_axis.id          #__
_axis.type        #__\_____
_axis.equipment   #__|_____|_____
_axis.depends_on  #__|_____|_____
_axis.vector[1]   #__|_____|_____
_axis.vector[2]   #__|_____|_____
_axis.vector[3]   #__|_____|_____
_axis.offset[1]   #__|_____|_____
_axis.offset[2]   #__|_____|_____
_axis.offset[3]   #__|_____|_____
#                # | | | | | | | | | |
#                # | | | | | | | | | |

#####|#####|#####|#####|##|##|###|##|##|
# The SLS Beamline axis convention is similar to the imgCIF convention, but |
# Y and Z run the other way | | | | | | | | |
#####|#####|#####|#####|##|##|###|##|##|
# | | | | | | | | | |
SLS_X      general      general      .      1 0 0 0 0 0
SLS_Y      general      general      .      0 -1 0 0 0 0
SLS_Z      general      general      .      0 0 -1 0 0 0
# | | | | | | | | | |
#####|#####|#####|#####|##|##|###|##|##|
# We define a kappa geometry with a left-handed omega and phi and a right- |
# handed kappa. The kappa axis arm is at the top when omega is zero | |
#####|#####|#####|#####|##|##|###|##|##|
# | | | | | | | | | |
GONIOMETER_OMEGA  rotation  goniometer      .      -1 0 0 . . .
GONIOMETER_KAPPA  rotation  goniometer      GONIOMETER_OMEGA
                                0.64279 0.76604
                                0 . . .

```

```

GONIOMETER_PHI    rotation    goniometer    GONIOMETER_KAPPA
                  -1  0  0  .  .  .
#####|#####|#####|#####|##|##|###|##|##|
BEAM              general    source    .    0  0 -1
.  .  .
GRAVITY           general    gravity    .    0 -1  0
.  .  .
#####|#####|#####|#####|##|##|###|##|##|
#              |          |          |          | | | | | | |
#####|#####|#####|#####|##|##|###|##|##|
# The detector is assumed to be mounted on an arm parallel to the beam | |
# with a DETECTOR_Y vertical translation and a pitch axis              | |
#####|#####|#####|#####|##|##|###|##|##|
#              |          |          |          | | | | | | |
DETECTOR_Z        translation detector    .    0  0 -1  0  0  0
DETECTOR_Y        translation detector    DETECTOR_Z    0 -1  0  0  0  0
DETECTOR_PITCH    rotation    detector    DETECTOR_Y    1  0  0  0  0  0
#####|#####|#####|#####|##|##|###|##|##|
# This detector is assumed to have the 0,0 corner at +X and -Y | | | | |
# we assume a 2463 x 2527 detector on a 0.172 mm pixel pitch | | | | |
#####|#####|#####|#####|##|##|###|##|##|
#              |          |          |          | | | | | | |
ELEMENT_X         translation detector    DETECTOR_PITCH 1  0  0
                  211.818 -217.322  0
ELEMENT_Y         translation detector    ELEMENT_X    0  1  0  0  0  0

```

4.2 The NeXus top level

The top level presented on the NeXus side of this mapping is

```
/CBF_diffn_scan__SCANID:NXentry
```

These scans can then be NeXus class instanceed into diffraction experiments and then into macromolecular studies by careful organization of files in tree of directories, but as the number of files and directories become large and individual scans get moved, there is a serious risk of loss of critical information if the identifiers for higher levels of the hierarchy are not recoverable from individual scans. For this reason, the fields

```

/CBF_scan_id="SCANID"
/CBF_diffn_id="DIFFRNID"
/CBF_entry_id="ENTRYID"

```

are given at the top level, if available, even if the same information is carried deeper in the hierarchy.

4.3 Identifying images

In many experimental disciplines, including X-ray crystallography, multiple related images are generated a single experiment. In some cases these can be gathered into a single array, but there are also many cases in which multiple data arrays are needed. Therefore, it is not sufficient to reserve the name “data” for the data. In imgCIF and array of data is identified by two identifiers, an array ID, specifying the structure

of the array, and a binary ID identifying a specific array. We name such data arrays by composing the commonly used name “data” these two identifiers separated by an underscores “_”.

When only one array ID and only one binary ID are present, then it will be sufficient to just use the field name “data”.

It is important to note that an image from an single detector may consist of only a slice (or “section”) taken from a larger array of data, and that a multi-element detector may be stored either as a single slice from one array common to all the elements, or as separate slices from different arrays, one array (or even more than one array) per detector element.

4.4 The ARRAY_DATA category

Data items in the ARRAY_DATA category are the containers for the array data items described in the category ARRAY_STRUCTURE.

It is recognized that the data in this category needs to be used in two distinct ways. During a data collection the lack of ancillary data and timing constraints in processing data may dictate the need to make a ‘miniCBF’ nothing more than an essential minimum of information to record the results of the data collection. In that case it is proper to use the ARRAY_DATA category as a container for just a single image and a compacted, beam-line dependent list of data collection parameter values. In such a case, only the tags ‘_array_data.header_convention’, ‘_array_data.header_contents’ and ‘_array_data.data’ need be populated.

For full processing and archiving, most of the tags in this dictionary will need to be populated.

The mapping from CBF into NeXus for the ARRAY_DATA category is:

```
_array_data.array_id ARRAYID
_array_data.binary_id BINARYID
_array_data.data DATAARRAY
_array_data.header_contents HEADER
_array_data.header_convention HEADERCONVENTION
```

→

```
/entry:NXentry
/instrument:NXinstrument
/DETECTORNAME:NXdetector_group
/DETECTORELEMENTNAME:NXdetector
/data_ARRAYID_BINARYID -->
    /entry/data_ARRAYID_BINARYID/data_ARRAYID_BINARYID
    @CBF_array_id="ARRAYID"
    @CBF_binary_id="BINARYID"
    @CBF_header_contents="HEADER"
    @CBF_header_convention="HEADERCONVENTION"
```

4.5 The ARRAY_DATA_EXTERNAL_DATA category

Data items in the ARRAY_DATA_EXTERNAL_DATA category optionally record the location and essential characteristics of arrays of data for use in ARRAY_DATA that are found external to the cif_img file.

The mapping from CBF into NeXus for the ARRAY_DATA_EXTERNAL_DATA category is:

```
_array_data.array_id ARRAYID
_array_data.binary_id BINARYID
_array_data.external_data_id EXTERNAL_DATA_ID
_array_data.header_contents HEADER
_array_data.header_convention HEADERCONVENTION

_array_data_external_data.archive_format ARCHIVEFORMAT
_array_data_external_data.archive_path ARCHIVEPATH
_array_data_external_data.id EXTERNAL_DATA_ID
_array_data_external_data.format HDF5
_array_data_external_data.uri DATAURI
_array_data_external_data.path DATAPATH
_array_data_external_data.frame DATAFRAME
```

→

```
/entry:NXentry
/data_ARRAYID_BINARYID:NXdata
  data_ARRAYID_BINARYID --> DATAURI//DATAPATH
/instrument:NXinstrument
  /DETECTORNAME:NXdetector_group
  /DETECTORELEMENTNAME:NXdetector
    /data_ARRAYID_BINARYID -->
      /entry/data_ARRAYID_BINARYID/data_ARRAYID_BINARYID
        @CBF_array_id="ARRAYID"
        @CBF_binary_id="BINARYID"
        @CBF_header_contents="HEADER"
        @CBF_header_convention="HEADERCONVENTION"
        @CBF_external_data_archive_format="ARCHIVEFORMAT"
        @CBF_external_data_archive_path="ARCHIVEPATH"
        @CBF_external_data_uri="DATAURI"
        @CBF_external_data_path="DATAPATH"
        @CBF_external_data_frame="DATAFRAME"
```

4.6 The ARRAY_ELEMENT_SIZE category

Data items in the ARRAY_ELEMENT_SIZE category record the physical size of array elements along each array dimension.

The mapping from CBF into NeXus for the ARRAY_ELEMENT_SIZE category is:

There is no general equivalent to ARRAY_ELEMENT_SIZE in NeXus, but for NXdetector_module, the pixel size is given as values in the array fast_pixel_direction and slow_pixel_direction.

4.7 The ARRAY_INTENSITIES category

Data items in the ARRAY_INTENSITIES category record the information required to recover the intensity data from the set of data values stored in the ARRAY_DATA category.

The detector may have a complex relationship between the raw intensity values and the number of incident photons. In most cases, the number stored in the final array will have a simple linear relationship to the actual number of incident photons, given by `_array_intensities.gain`. If raw, uncorrected values are presented (e.g. for calibration experiments), the value of `_array_intensities.linearity` will be 'raw' and `_array_intensities.gain` will not be used.

The mapping from CBF into NeXus for the ARRAY_DATA_EXTERNAL_DATA category is:

```
_array_intensities.array_id ARRAYID
_array_intensities.binary_id BINARYID
_array_intensities.details DETAILS
_array_intensities.gain GAIN
_array_intensities.gain_esd GAINESD
_array_intensities.linearity LINEARITY
_array_intensities.offset OFFSET
_array_intensities.scaling SCALING
_array_intensities.overload OVERLOAD
_array_intensities.undefval UNDEFVAL
_array_intensities.underload UNDERLOAD
_array_intensities.pixel_fast_bin_size FBINSIZE
_array_intensities.pixel_slow_bin_size SBINSIZE
_array_intensities.pixel_binning_method METHOD
```

→

```
/entry:NXentry
/data_ARRAYID_BINARYID:NXdata
/data_ARRAYID_BINARYID
  @CBF_array_id="ARRAYID"
  @CBF_binary_id="BINARYID"
  @details="DETAILS"
  @gain=[GAIN]
  @gain_esd=[GAINESD]
  @linearity="LINEARITY"
  @offset=[OFFSET]
  @saturation_value=[OVERLOAD]
  @scaling_factor=[SCALING]
  @undefval=[UNDEFVAL]
  @underload_value=[UNDERLOAD]
  @CBF_array_intensities__pixel_fast_bin_size=[FBINSIZE]
```



```

        @CBF_array_intensities__pixel_slow_bin_size=[SBINSIZE]
        @CBF_array_intensities__pixel_binning_method="METHOD"
/instrument:NXinstrument
/DETECTORNAME:NXdetector_group
/DETECTORELEMENTNAME:NXdetector
/data_ARRAYID_BINARYID -->
    /entry/data_ARRAYID_BINARYID/data_ARRAYID_BINARYID

```

The argument has been made that these attributes are not needed because NeXus files are supposed to have 'true values' stored. In many cases that is true and then none of these attributes are needed. However, with some detectors and some experiments there are good technical and scientific reasons to bring in values that will need processing later to derive 'true values', and in those case some or all of these attributes will be needed. They are provided for such cases.

The same attributes could be used as fields in the case of a single data array, but in that case links for all the fields would be needed from NXdata to NXdetector, so it is preferable to use attributes even in the case of a single data array. The reverse mapping will support both uses.

4.8 The ARRAY_STRUCTURE category

Data items in the ARRAY_STRUCTURE category record the organization and encoding of array data that may be stored in the ARRAY_DATA category.

Note that this is essentially a type that may apply to multiple binary images, and corresponds to some of the detailed HDF5 information about an array. The following mapping is a placeholder for the names given for future reference, if needed.

The information in this category is the byte order, the compression information, and the encoding, which is carried in and retrievable from the HDF5 types, properties lists, etc.

At present NeXus does not expose this information. This should be discussed.

4.9 The ARRAY_STRUCTURE_LIST category

Data items in the ARRAY_STRUCTURE_LIST category record the size and organization of each array dimension.

The relationship to physical axes may be given.

The mapping from CBF into NeXus for the ARRAY_STRUCTURE_LIST category is:

```
_array_structure_list.axis_set_id AXISSETID
```

```

_array_structure_list.array_id ARRAYID
_array_structure_list.array_section_id ARRAYSECTIONID
_array_structure_list.dimension DIM
_array_structure_list.direction DIR
_array_structure_list.index INDEX
_array_structure_list.precedence PRECEDENCE

_array_structure_list_section.array_id ARRAYID -->
_array_structure_list_section.id ARRAYSECTIONID -->
_array_structure_list_section.index PRECEDENCE-->
_array_structure_list_section.end END -->
_array_structure_list_section.start START -->
_array_structure_list_section.stride STRIDE -->

loop_
_array_structure_list_axis.axis_id
_array_structure_list_axis.axis_set_id
_array_structure_list_axis.angle
_array_structure_list_axis.angle_increment
_array_structure_list_axis.displacement
_array_structure_list_axis.fract_displacement
_array_structure_list_axis.displacement_increment
_array_structure_list_axis.fract_displacement_increment
_array_structure_list_axis.angular_pitch
_array_structure_list_axis.reference_angle
_array_structure_list_axis.reference_displacement REFDISP
AXISID1 AXISSETID ANGLE1 ANGLEINC1 DISP1 FRACTDISP1
    DISPINC1 FRACTINC1 ANGPITCH1 REFANG1
AXISID2 AXISSETID ANGLE2 ANGLEINC2 DISP2 FRACTDISP2
    DISPINC2 FRACTINC2 ANGPITCH2 REFANG2
AXISID3 AXISSETID ANGLE3 ANGLEINC3 DISP3 FRACTDISP3
    DISPINC3 FRACTINC3 ANGPITCH3 REFANG3

_diffn_data_frame.array_id ARRAYID
_diffn_data_frame.binary_id BINARYID
_diffn_data_frame.center_fast CENF
_diffn_data_frame.center_slow CENS
_diffn_data_frame.center_derived CENDERIVED
_diffn_data_frame.center_units UNITS
_diffn_data_frame.detector_element_id ELEMENTID
_diffn_data_frame.id FRAMEID
_diffn_data_frame.details DETAILS

→

...

/entry:NXentry
/data_ARRAYID_BINARYID:NXdata
    @signal="data_ARRAYID_BINARY_ID"

```

```

/data_ARRAYID_BINARYID[ the data for the array ARRAYID,
                        binary BINARYID, all sections,
                        all FRAMES]
@axes=[...,AXISID1,...] with AXISID1 inserted at PRECEDENCE-1
@AXISID1_indices=[PRECEDENCE-1]
@AXISID2_indices=[PRECEDENCE-1]
@AXISID3_indices=[PRECEDENCE-1]
@AXISID1_origins=[origin1] (default 0)
@AXISID2_origins=[origin2] (default 0)
@AXISID3_origins=[origin3] (default 0)
@AXISID1_sizes=[size1]
@AXISID2_sizes=[size2]
@AXISID3_sizes=[size3]
@AXISID1_strides=[stride1]
@AXISID2_strides=[stride2]
@AXISID3_strides=[stride3]

...
/AXISID1 -->
    /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID1
/AXISID2 -->
    /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID2
/AXISID2 -->
    /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID3
...

/instrument:NXinstrument
/DETECTORNAME:NXdetector_group
/DETECTORELEMENTNAME:NXdetector
    /data_ARRAYID_BINARYID -->
        /entry/data_ARRAYID_BINARYID/data_ARRAYID_BINARYID
    /AXISID1 -->
        /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID1
    /AXISID2 -->
        /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID2
    /AXISID2 -->
        /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID3
/ARRAYSECTIONID:NXdetector_module
    /data_origin=[...] -- the 0-based origins indices of ARRAYSECTIONID
    /data_size=[...] the sizes in pixels of ARRAYSECTIONID
    /data_stride=[...] the strides of ARRAYSECTIONID
    ..
/transformations:NXtransformations
    /AXISID1=[DISP1,DISP1+DISPINC1,...]
        (or using angles where appropriate)
    @depends_on=... determined from AXIS definitions
    @equipment="detector"
    @offset=[...] determined from AXIS definitions
    @offset_units="mm"
    @transformation_type="..." from AXIS definitions

```

```

@units="mm"
@vector=[...] determined from AXIS definitions
@CBF_array_structure_list_axis__axis_id="AXISID1"
@CBF_array_structure_list_axis__axis_set_id="AXISSETID"
@CBF_array_structure_list_axis__angle=ANGLE1
@CBF_array_structure_list_axis__angle_increment=ANGLEINC1
@CBF_array_structure_list_axis__displacement=DISP1
@CBF_array_structure_list_axis__displacement=FRACDISP1
@CBF_array_structure_list_axis__displacement_increment=DISPINC1
@CBF_array_structure_list_axis__fract_displacement_increment
    =FRACTINC1
@CBF_array_structure_list_axis__angular_pitch=ANGPITCH1
@CBF_array_structure_list_axis__radial_pitch=RADPITCH1
@CBF_array_structure_list_axis__reference_angle=REFANG1
@CBF_array_structure_list_axis__reference_displacement=REFDISP1
/AXISID2=[DISP2,DISP2+DISPINC2,...]
    (or using angles where appropriate)
@depends_on=... determined from AXIS definitions
@equipment="detector"
@offset=[...] determined from AXIS definitions
@offset_units="mm"
@transformation_type="..." from AXIS definitions
@units="mm"
@vector=[...] determined from AXIS definitions
@CBF_array_structure_list_axis__axis_id="AXISID2"
@CBF_array_structure_list_axis__axis_set_id="AXISSETID"
@CBF_array_structure_list_axis__angle=ANGLE2
@CBF_array_structure_list_axis__angle_increment=ANGLEINC2
@CBF_array_structure_list_axis__displacement=DISP2
@CBF_array_structure_list_axis__displacement=FRACDISP2
@CBF_array_structure_list_axis__displacement_increment=DISPINC2
@CBF_array_structure_list_axis__fract_displacement_increment
    =FRACTINC2
@CBF_array_structure_list_axis__angular_pitch=ANGPITCH2
@CBF_array_structure_list_axis__radial_pitch=RADPITCH2
@CBF_array_structure_list_axis__reference_angle=REFANG2
@CBF_array_structure_list_axis__reference_displacement=REFDISP2
/AXISID3=[DISP3,DISP3+DISPINC3,...]
    (or using angles where appropriate)
@depends_on=... determined from AXIS definitions
@equipment="detector"
@offset=[...] determined from AXIS definitions
@offset_units="mm"
@transformation_type="..." from AXIS definitions
@units="mm"
@vector=[...] determined from AXIS definitions
@CBF_array_structure_list_axis__axis_id="AXISID3"
@CBF_array_structure_list_axis__axis_set_id="AXISSETID"
@CBF_array_structure_list_axis__angle=ANGLE3
@CBF_array_structure_list_axis__angle_increment=ANGLEINC3

```

```

@CBF_array_structure_list_axis__displacement=DISP3
@CBF_array_structure_list_axis__displacement=FRACDISP3
@CBF_array_structure_list_axis__displacement_increment=DISPINC3
@CBF_array_structure_list_axis__fract_displacement_increment
    =FRACTINC3
@CBF_array_structure_list_axis__angular_pitch=ANGPITCH3
@CBF_array_structure_list_axis__radial_pitch=RADPITCH3
@CBF_array_structure_list_axis__reference_angle=REFANG3
@CBF_array_structure_list_axis__reference_displacement=REFDISP3

```

The same axis AXISIDn may appear in multiple axis sets for different values of PRECEDENCE of the data array, in which case the values in AXISIDn_indices will be the sorted list of PRECEDENCE-1 values and the array section information will be organized by the same ordering.

This mapping of precedence of axes and directions may seem redundant with the array storage parameters maintained by HDF5, but changed information is used by applications to remap images to undo rotations and flips and needs to be accessible without rewriting the image, which may be impractical. Failure to expose at least this information in the NeXus will greatly restrict the range of applications that will be able to use the NeXus version of these files.

4.10 The ARRAY_STRUCTURE_LIST_SECTION category

Data items in the ARRAY_STRUCTURE_LIST_SECTION category identify the dimension-by-dimension start, end and stride of each section of an array that is to be referenced.

For any array with identifier ARRAYID, array section ids of the form ARRAYID(start1:end1:stride1,start2:end2:stride2, ...) are defined by default.

For the given index, the elements in the section are of indices:

```

_array_structure_list_section.start,
_array_structure_list_section.start + _array_structure_list_section.stride,
_array_structure_list_section.start + 2*_array_structure_list_section.stride,
...

```

stopping either when the indices leave the limits of the indices of that dimension or

```

[min(_array_structure_list_section.start, _array_structure_list_section.end),
 max(_array_structure_list_section.start, _array_structure_list_section.end)].

```

The ordering of these elements is determined by the overall ordering of _array_structure_list_section.array_id and not by the ordering implied

by the stride.

The mapping from CBF into NeXus for the ARRAY_STRUCTURE_LIST_SECTION category is:

```
_array_structure_list_section.array_id ARRAYID -->
_array_structure_list_section.id SECTIONID -->
_array_structure_list_section.index INDEX-->
_array_structure_list_section.end END -->
_array_structure_list_section.start START -->
_array_structure_list_section.stride STRIDE -->
/instrument:NXinstrument
  /DETECTORNAME:NXdetector_group
  /DETECTORELEMENTNAME:NXdetector
    /data_ARRAYID_BINARYID -->
      /entry/data_ARRAYID_BINARYID/data_ARRAYID_BINARYID
    /AXISID1 -->
      /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID1
    /AXISID2 -->
      /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID2
    /AXISID2 -->
      /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID3
  /ARRAYSECTIONID:NXdetector_module
    /data_origin=[...] -- the 0-based origins indices of ARRAYSECTIONID
    /data_size=[...] the sizes in pixels of ARRAYSECTIONID
    /data_stride=[...] the strides of ARRAYSECTIONID
```

4.11 The ARRAY_STRUCTURE_LIST_AXIS category

Data items in the ARRAY_STRUCTURE_LIST_AXIS category describe the physical settings of sets of axes for the centres of pixels that correspond to data points described in the ARRAY_STRUCTURE_LIST category.

In the simplest cases, the physical increments of a single axis correspond to the increments of a single array index. More complex organizations, e.g. spiral scans, may require coupled motions along multiple axes.

Note that a spiral scan uses two coupled axes: one for the angular direction and one for the radial direction. This differs from a cylindrical scan for which the two axes are not coupled into one set.

Axes may be specified either for an entire array or for just a section of an array.

The mapping from CBF into NeXus for the ARRAY_STRUCTURE_LIST_AXIS category is:

```
_array_structure_list.axis_set_id AXISSETID
_array_structure_list.array_id ARRAYID
```

```

_array_structure_list.array_section_id ARRAYSECTIONID
_array_structure_list.dimension DIM
_array_structure_list.direction DIR
_array_structure_list.index INDEX
_array_structure_list.precedence PRECEDENCE

_array_structure_list_section.array_id ARRAYID -->
_array_structure_list_section.id ARRAYSECTIONID -->
_array_structure_list_section.index PRECEDENCE-->
_array_structure_list_section.end END -->
_array_structure_list_section.start START -->
_array_structure_list_section.stride STRIDE -->

loop_
_array_structure_list_axis.axis_id
_array_structure_list_axis.axis_set_id
_array_structure_list_axis.angle
_array_structure_list_axis.angle_increment
_array_structure_list_axis.displacement
_array_structure_list_axis.fract_displacement
_array_structure_list_axis.displacement_increment
_array_structure_list_axis.fract_displacement_increment
_array_structure_list_axis.angular_pitch
_array_structure_list_axis.reference_angle
_array_structure_list_axis.reference_displacement REFDISP
AXISID1 AXISSETID ANGLE1 ANGLEINC1 DISP1 FRACTDISP1
    DISPINC1 FRACTINC1 ANGPITCH1 REFANG1
AXISID2 AXISSETID ANGLE2 ANGLEINC2 DISP2 FRACTDISP2
    DISPINC2 FRACTINC2 ANGPITCH2 REFANG2
AXISID3 AXISSETID ANGLE3 ANGLEINC3 DISP3 FRACTDISP3
    DISPINC3 FRACTINC3 ANGPITCH3 REFANG3

_diffraction_data_frame.array_id ARRAYID
_diffraction_data_frame.binary_id BINARYID
_diffraction_data_frame.center_fast CENF
_diffraction_data_frame.center_slow CENS
_diffraction_data_frame.center_units UNITS
_diffraction_data_frame.detector_element_id ELEMENTID
_diffraction_data_frame.id FRAMEID
_diffraction_data_frame.details DETAILS

→

...

/entry:NXentry
/data_ARRAYID_BINARYID:NXdata
    @signal="data_ARRAYID_BINARY_ID"
    /data_ARRAYID_BINARYID[ the data for the array ARRAYID,

```

```

        binary BINARYID, all sections,
        all FRAMES]
@axes=[...,AXISID1,...] with AXISID1 inserted at PRECEDENCE-1
@AXISID1_indices=[PRECEDENCE-1]
@AXISID2_indices=[PRECEDENCE-1]
@AXISID3_indices=[PRECEDENCE-1]
@AXISID1_origins=[origin1] (default 0)
@AXISID2_origins=[origin2] (default 0)
@AXISID3_origins=[origin3] (default 0)
@AXISID1_sizes=[size1]
@AXISID2_sizes=[size2]
@AXISID3_sizes=[size3]
@AXISID1_strides=[stride1]
@AXISID2_strides=[stride2]
@AXISID3_strides=[stride3]

...
/AXISID1 -->
    /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID1
/AXISID2 -->
    /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID2
/AXISID2 -->
    /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID3
...

/instrument:NXinstrument
/DETECTORNAME:NXdetector_group
/DETECTORELEMENTNAME:NXdetector
    /data_ARRAYID_BINARYID -->
        /entry/data_ARRAYID_BINARYID/data_ARRAYID_BINARYID
    /AXISID1 -->
        /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID1
    /AXISID2 -->
        /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID2
    /AXISID2 -->
        /entry/instrument/DETECTORELEMENTNAME/transformations/AXISID3
/ARRAYSECTIONID:NXdetector_module
    /data_origin=[...] -- the 0-based origins indices of ARRAYSECTIONID
    /data_size=[...] the sizes in pixels of ARRAYSECTIONID
    /data_stride=[...] the strides of ARRAYSECTIONID
    ..
/transformations:NXtransformations
    /AXISID1=[DISP1,DISP1+DISPINC1,...] (or using angles where appropriate)
        @depends_on=... determined from AXIS definitions
        @equipment="detector"
        @offset=[...] determined from AXIS definitions
        @offset_units="mm"
        @transformation_type="..." from AXIS definitions
        @units="mm"
        @vector=[...] determined from AXIS definitions

```



```

@CBF_array_structure_list_axis__axis_id="AXISID1"
@CBF_array_structure_list_axis__axis_set_id="AXISSETID"
@CBF_array_structure_list_axis__angle=ANGLE1
@CBF_array_structure_list_axis__angle_increment=ANGLEINC1
@CBF_array_structure_list_axis__displacement=DISP1
@CBF_array_structure_list_axis__displacement=FRACDISP1
@CBF_array_structure_list_axis__displacement_increment=DISPINC1
@CBF_array_structure_list_axis__fract_displacement_increment=FRACTINC1
@CBF_array_structure_list_axis__angular_pitch=ANGPITCH1
@CBF_array_structure_list_axis__radial_pitch=RADPITCH1
@CBF_array_structure_list_axis__reference_angle=REFANG1
@CBF_array_structure_list_axis__reference_displacement=REFDISP1
/AXISID2=[DISP2,DISP2+DISPINC2,...] (or using angles where appropriate)
@depends_on=... determined from AXIS definitions
@equipment="detector"
@offset=[...] determined from AXIS definitions
@offset_units="mm"
@transformation_type="..." from AXIS definitions
@units="mm"
@vector=[...] determined from AXIS definitions
@CBF_array_structure_list_axis__axis_id="AXISID2"
@CBF_array_structure_list_axis__axis_set_id="AXISSETID"
@CBF_array_structure_list_axis__angle=ANGLE2
@CBF_array_structure_list_axis__angle_increment=ANGLEINC2
@CBF_array_structure_list_axis__displacement=DISP2
@CBF_array_structure_list_axis__displacement=FRACDISP2
@CBF_array_structure_list_axis__displacement_increment=DISPINC2
@CBF_array_structure_list_axis__fract_displacement_increment=FRACTINC2
@CBF_array_structure_list_axis__angular_pitch=ANGPITCH2
@CBF_array_structure_list_axis__radial_pitch=RADPITCH2
@CBF_array_structure_list_axis__reference_angle=REFANG2
@CBF_array_structure_list_axis__reference_displacement=REFDISP2
/AXISID3=[DISP3,DISP3+DISPINC3,...] (or using angles where appropriate)
@depends_on=... determined from AXIS definitions
@equipment="detector"
@offset=[...] determined from AXIS definitions
@offset_units="mm"
@transformation_type="..." from AXIS definitions
@units="mm"
@vector=[...] determined from AXIS definitions
@CBF_array_structure_list_axis__axis_id="AXISID3"
@CBF_array_structure_list_axis__axis_set_id="AXISSETID"
@CBF_array_structure_list_axis__angle=ANGLE3
@CBF_array_structure_list_axis__angle_increment=ANGLEINC3
@CBF_array_structure_list_axis__displacement=DISP3
@CBF_array_structure_list_axis__displacement=FRACDISP3
@CBF_array_structure_list_axis__displacement_increment=DISPINC3
@CBF_array_structure_list_axis__fract_displacement_increment=FRACTINC3
@CBF_array_structure_list_axis__angular_pitch=ANGPITCH3
@CBF_array_structure_list_axis__radial_pitch=RADPITCH3

```

```
@CBF_array_structure_list_axis__reference_angle=REFANG3
@CBF_array_structure_list_axis__reference_displacement=REFDISP3
```

The same axis AXISIDn may appear in multiple axis sets for different values of PRECEDENCE of the data array, in which case the values in AXISIDn_indices will be the sorted list of PRECEDENCE-1 values and the array section information will be organized by the same ordering.

4.12 The AXIS category

Data items in the AXIS category record the information required to describe the various goniometer, detector, source and other axes needed to specify a data collection or the axes defining the coordinate system of an image.

The location of each axis is specified by two vectors: the axis itself, given by a unit vector in the direction of the axis, and an offset to the base of the unit vector.

The vectors defining an axis are referenced to an appropriate coordinate system. The axis vector, itself, is a dimensionless unit vector. Where meaningful, the offset vector is given in millimetres. In coordinate systems not measured in metres. the offset is not specified and is taken as zero.

The available coordinate systems are:

- The imgCIF standard laboratory coordinate system
- The direct lattice (fractional atomic coordinates)
- The orthogonal Cartesian coordinate system (real space)
- The reciprocal lattice
- An abstract orthogonal Cartesian coordinate frame

For consistency in this discussion, we call the three coordinate system axes X, Y and Z. This is appropriate for the imgCIF standard laboratory coordinate system, and last two Cartesian coordinate systems, but for the direct lattice, X corresponds to a, Y to b and Z to c, while for the reciprocal lattice, X corresponds to a*, Y to b* and Z to c*.


For purposes of visualization, all the coordinate systems are taken as right-handed, i.e., using the convention that the extended thumb of a right hand could point along the first (X) axis, the straightened pointer finger could point along the second (Y) axis and the middle finger folded inward could point along the third (Z) axis. |

THE IMGCIF STANDARD LABORATORY COORDINATE SYSTEM

The imgCIF standard laboratory coordinate system is a right-handed orthogonal coordinate similar to the MOSFLM coordinate system, but imgCIF puts Z along the X-ray beam, rather than putting X along the X-ray beam as in MOSFLM.

The vectors for the imgCIF standard laboratory coordinate system form a right-handed Cartesian coordinate system with its origin in the sample or specimen. The origin of the axis system should, if possible, be defined in terms of mechanically stable axes to be both in the sample and in the beam. If the sample goniometer or other sample positioner has two axes the intersection of which defines a unique point at which the sample should be mounted to be bathed by the beam, that will be the origin of the axis system. If no such point is defined, then the midpoint of

the line of intersection between the sample and the center of the beam will define the origin. For this definition the sample positioning system will be set at its initial reference position for the experiment.



Axis 1 (X): The X-axis is aligned to the mechanical axis pointing from the sample or specimen along the principal axis of the goniometer or sample positioning system if the sample positioning system has an axis that intersects the origin and which form an angle of more than 22.5 degrees with the beam axis.

Axis 2 (Y): The Y-axis completes an orthogonal right-handed system defined by the X-axis and the Z-axis (see below).

Axis 3 (Z): The Z-axis is derived from the source axis which goes from the sample to the source. The Z-axis is the component of the source axis in the direction of the source orthogonal to the X-axis in the plane defined by the X-axis and the source axis.

If the conditions for the X-axis can be met, the coordinate system will be based on the goniometer or other sample positioning system and the beam and not on the orientation of the detector, gravity etc. The vectors necessary to specify all other axes are given by sets of three components in the order (X, Y, Z). If the axis involved is a rotation axis, it is right-handed, i.e. as one views the object to be rotated from the origin (the tail) of the unit vector, the rotation is clockwise. If a translation axis is specified, the direction of the unit vector specifies the sense of positive translation.

Note: This choice of coordinate system is similar to but significantly different from the choice in MOSFLM [6]. In MOSFLM, X is along the X-ray beam (the CBF/imgCIF Z axis) and Z is along the rotation axis.

In some experimental techniques, there is no goniometer or the principal axis of the goniometer is at a small acute angle with respect to the source axis. In such cases, other reference axes are needed to define a useful coordinate system. The order of priority in defining directions in such cases is to use the detector, then gravity, then north.

If the X-axis cannot be defined as above, then the direction (not the origin) of the X-axis should

be parallel to the axis of the primary detector element corresponding to the most rapidly varying dimension of that detector element's data array, with its positive sense corresponding to increasing values of the index for that dimension. If the detector is such that such a direction cannot be defined (as with a point detector) or that direction forms an angle of less than 22.5 degrees with respect to the source axis, then the X-axis should be chosen so that if the Y-axis is chosen in the direction of gravity, and the Z-axis is chosen to be along the source axis, a right-handed orthogonal coordinate system is chosen. In the case of a vertical source axis, as a last resort, the X-axis should be chosen to point North.

All rotations are given in degrees and all translations are given in mm.

Axes may be dependent on one another. The X-axis is the only goniometer axis the direction of which is strictly connected to the hardware. All other axes are specified by the positions they would assume when the axes upon which they depend are at their zero points.

When specifying detector axes, the axis is given to the beam centre. The location of the beam centre on the detector should be given in the DIFFRN_DETECTOR category in distortion-corrected millimetres from the (0,0) corner of the detector.

It should be noted that many different origins arise in the definition of an experiment. In particular, as noted above, it is necessary to specify the location of the beam centre on the detector in terms of the origin of the detector, which is, of course, not coincident with the centre of the sample.

The unit cell, reciprocal cell and crystallographic orthogonal Cartesian coordinate system are defined by the CELL and the matrices in the ATOM_SITES category.

THE DIRECT LATTICE (FRACTIONAL COORDINATES)

The direct lattice coordinate system is a system of fractional coordinates aligned to the crystal, rather than to the laboratory. This is a natural coordinate system for maps and atomic coordinates. It is the simplest coordinate system in which to apply symmetry. The axes are determined by the cell edges, and are not necessarily orthogonal. This coordinate system is not uniquely defined and depends on the cell parameters in the CELL category and the settings chosen to index the crystal.

Molecules in a crystal studied by X-ray diffraction are organized into a repeating regular array of unit cells. Each unit cell is defined by three vectors, a, b and c. To quote from [4],

“The choice of the unit cell is not unique and therefore, guidelines have been established for selecting the standard basis vectors and the origin. They are based on symmetry and metric considerations:

- “(1) The axial system should be right handed.
- (2) The basis vectors should coincide as much as possible with directions of highest symmetry.
- (3) The cell taken should be the smallest one that satisfies|
- (4) Of all the lattice vectors, none is shorter than a.
- (5) Of those not directed along a, none is shorter than b.
- (6) Of those not lying in the ab plane, none is shorter than c.
- (7) The three angles between the basis vectors a, b and c are either all acute (< 90 degrees)

or all obtuse (greater than equal to 90 degrees).''

These rules do not produce a unique result that is stable under the assumption of experimental errors, and the the resulting cell may not be primitive.

In this coordinate system, the vector (.5, .5, .5) is in the middle of the given unit cell.

Grid coordinates are an important variation on fractional coordinates used when working with maps. In imgCIF, the conversion from fractional to grid coordinates is implicit in the array indexing specified by `_array_structure_list.dimension`. Note that this implicit grid-coordinate scheme is 1-based, not zero-based, *i.e.* the origin of the cell for axes along the cell edges with no specified `_array_structure_list.axis.displacement` will have grid coordinates of (1,1,1), *i.e.* array indices of (1,1,1).

THE ORTHOGONAL CARTESIAN COORDINATE SYSTEM (REAL SPACE)

The orthogonal Cartesian coordinate system is a transformation of the direct lattice to the actual physical coordinates of atoms in space. It is similar to the laboratory coordinate system, but is anchored to and moves with the crystal, rather than being anchored to the laboratory. The transformation from fractional to orthogonal cartesian coordinates is given by the `_atom_sites.Cartn_transf_matrix[i][j]` and `_atom_sites.Cartn_transf_vector[i]` tags. A common choice for the matrix of the transformation is given in the 1992 PDB format document

$$\begin{vmatrix} a & b \cos(\gamma) & c \cos(\beta) \\ 0 & b \sin(\gamma) & c (\cos(\alpha) - \cos(\beta) \cos(\gamma)) / \sin(\gamma) \\ 0 & 0 & V / (a b \sin(\gamma)) \end{vmatrix}$$

This is a convenient coordinate system in which to do fitting of models to maps and in which to understand the chemistry of a molecule.

THE RECIPROCAL LATTICE

The reciprocal lattice coordinate system is used for diffraction intensities. It is based on the reciprocal cell, the dual of the cell, in which reciprocal cell edges are derived from direct cell faces:

$$\begin{aligned} a^* &= bc \sin(\alpha) / V & b^* &= ac \sin(\beta) / V & c^* &= ab \sin(\gamma) / V \\ \cos(\alpha^*) &= (\cos(\beta) \cos(\gamma) - \cos(\alpha)) / (\sin(\beta) \sin(\gamma)) \\ \cos(\beta^*) &= (\cos(\alpha) \cos(\gamma) - \cos(\beta)) / (\sin(\alpha) \sin(\gamma)) \\ \cos(\gamma^*) &= (\cos(\alpha) \cos(\beta) - \cos(\gamma)) / (\sin(\alpha) \sin(\beta)) \\ V &= abc \sqrt{(1 - \cos(\alpha)^2 - \cos(\beta)^2 - \cos(\gamma)^2 + 2\cos(\alpha)\cos(\beta)\cos(\gamma))} \end{aligned}$$

In this form the dimensions of the reciprocal lattice are in reciprocal Angstroms (\AA^{-1}). A dimensionless form can be obtained by multiplying by the wavelength. Reflections are commonly indexed against this coordinate system as (h, k, l) triples. See [8] and [7].

The mapping from CBF into NeXus for the AXIS category is::

<code>_axis.id</code>	AXISID
<code>_axis.type</code>	AXISTYPE

_axis.equipment	AXISEQUIPMENT
_axis.equipment_component	AXISEQUIPCOMP
_axis.depends_on	AXISDEPENDSON
_axis.rotation_axis	AXISROTAXIS
_axis.rotation	AXISROTATION
_axis.vector[1]	AXISV1
_axis.vector[2]	AXISV2
_axis.vector[3]	AXISV3
_axis.offset[1]	AXISO1
_axis.offset[2]	AXISO2
_axis.offset[3]	AXISO3
_axis.system	AXISSYSTEM

→

```
{
  /entry:NXentry
    /instrument:NXinstrument
      /DETECTORELEMENTNAME:NXdetector
for AXISEQUIPMENT=="detector"}
{
  /entry:NXentry
    /sample:NXsample
for AXISEQUIPMENT=="goniometer"}
{
  /entry:NXentry
for AXISEQUIPMENT=="general"}
  /transformations:NXtransformations
    /AXISID=[]
      @units="mm" if AXISTYPE=="translation"
      or @units="degrees" if AXISTYPE=="rotation"
      @transformation_type="AXISTYPE"
      @equipment_component="AXISEQUIPCOMP"
      @depends_on="AXISDEPENDSON"
      @rotation_axis="AXISROTAXIS"
      @rotation=AXISROTATION
      @rotation_units="degrees"
      @offset=offsetxform([01,02,03])
      @offset_inits="mm"
      @vector=coorndxform([V1,V2,V3])
```

- _axis.system →
Only a laboratory coordinate system is handled in NeXus at present. The addition of more coordinate systems in NeXus is pending. See the discussion of transformation to the McStas coordinate system below.
- _axis.vector[1] V1 →
/CBF_axis__AXISID=[]
@vector=coorndxform([V1,V2,V3])
- _axis.vector[2] V2 →
/CBF_axis__AXISID=[]
@vector=coorndxform([V1,V2,V3])

- `_axis.vector[3] V3 →`
`/CBF_axis_AXISID=[]`
`@vector=cooridxform([V1,V2,V3])`
- `_axis.variant → ??` NeXus does not handle variants at this time

4.12.1 Differences in Coordinate Frames

The standard coordinate frame in NeXus is the McStas coordinate frame, in which the Z-axis points in the direction of the incident beam, the X-axis is orthogonal to the Z-axis in the horizontal plane and pointing left as seen from the source and the Y-axis points upwards. The origin is in the sample.

The standard coordinate frame in imgCIF/CBF aligns the X-axis to the principal goniometer axis, chooses the Z-axis to point from the sample into the beam. If the beam is not orthogonal to the X-axis, the Z-axis is the component of the vector points into the beam orthogonal to the X-axis. The Y-axis is chosen to complete a right-handed axis system.

Let us call the NeXus coordinate axes, X_{nx} , Y_{nx} and Z_{nx} and the imgCIF/CBF coordinate axes, X_{cbf} , Y_{cbf} and Z_{cbf} and the direction of gravity, $Gravity$. In order to translate a vector $v_{nx} = (x, y, z)$ from the NeXus coordinate system to the imgCIF coordinate system, we also need two additional axes, as unit vectors, $Gravity_{cbf}$, the downwards direction, and $Beam_{cbf}$, the direction of the beam (e.g. $(0, 0, -1)$).

In practice, the beam is not necessarily perfectly horizontal, so Y_{nx} is not necessarily perfectly vertical. Therefore, in order to generate X_{nx} , Y_{nx} and Z_{nx} some care is needed. The cross product between two vectors \vec{a} and \vec{b} is a new vector \vec{c} orthogonal to both \vec{a} and \vec{b} , chosen so that $\vec{a}, \vec{b}, \vec{c}$ is a right handed system. If \vec{a} and \vec{b} are orthogonal unit vectors, this right-handed system is an orthonormal coordinate system.

In the CBF coordinate frame, Z_{nx} is aligned to $Beam_{cbf}$:

$$Z_{nx} = Beam_{cbf}$$

X_{nx} is defined as being horizontal at right angles to the beam, pointing to the left when seen from the source. Assuming the beam is not vertical, we can compute X_{nx} as the normalized cross product of the beam and the gravity:

$$X_{nx} = (Beam_{cbf} \times Gravity_{cbf}) / ||Beam_{cbf} \times Gravity_{cbf}||$$

To see that this satisfies the constraint of being horizontal and pointing to the left, consider the case of $Beam = (0, 0, -1)$ and $Gravity = (0, 0, 1)$ then we would have $X_{nx} = (1, 0, 0)$ from the cross product above. The normalization is only necessary if the beam is not horizontal.

Finally Y_{nx} is computed as the cross product of the beam and X_{nx} , completing a orthonormal right-handed system with Y_{nx} pointing upwards:

$$Y_{nx} = Beam_{cbf} \times X_{nx}$$

Then we know that in the imgCIF/CBF coordinate frame

$$v_{nx} = x \cdot X_{nx} + y \cdot Y_{nx} + z \cdot Z_{nx}$$

Thus, given the imgCIF/CBF vectors for the true direction of the beam and the true direction of gravity, we have a linear transformation from the NeXus coordinate frame to the imgCIF/CBF coordinate frame. The origins of the two frames agree. The inverse linear transformation will transform a vector in the imgCIF/CBF coordinate frame into the NeXus coordinate frame.

In the common case in which the beam is orthogonal to the principal goniometer axis so that $Beam_{cbf} = (0, 0, -1)$ and the imgCIF/CBF Y-axis points upwards, the transformation inverts the X

and Z axes. In the other common case in which the beam is orthogonal to the principal goniometer axis and the imgCIF/CBF Y-axis points downwards, the transformation inverts the Y and Z axes.

4.12.2 Mapping Axes

There are two transformations needed: `coordxform(v)` which takes a vector, `v`, in the CBF imgCIF Standard Laboratory Coordinate System and returns the equivalent McStas coordinate vector, and `offsetxform(o)` which takes an offset, `o`, in the CBF imgCIF Standard Laboratory Coordinate System and returns the equivalent NeXus offset. As of this writing, it has not been decided as to whether the NeXus offset should also be relative (in which case `offsetxform = coordxform`) or whether the NeXus offset should be absolute.

In imgCIF/CBF all the information about all axes other than their settings are gathered in one `AXIS` category. The closest equivalent container in NeXus is the `NXinstrument` class. We put the information about detector axes into an `detector:NXdetector` NeXus class instance, information about the goniometer into an `goniometer:NXsample` NeXus class instance, etc. Additionally, in view of the general nature of some axes, such as the coordinate frame axes and gravity, we add a `coordinate.system:NXcoordinate_system` NeXus class instance with `axis_gravity`, `axis_beam` and other axes not tied to specific equipment.

We have applied the coordinate frame transformation changing the CBF laboratory coordinates into McStas coordinates. Notice that X and Z have changed direction, but Y has not. In other experimental setup, other transformations may occur. The offsets for dependent axes are given relative to the total offset of axes on which that axis is dependent. Note that the axis settings do not enter into this calculation, because the offsets of dependent axes are given with all axes at their zero settings.

The `cbf_location` attribute gives a mapping back into the CBF `AXIS` category in dotted notation. The first component is the data block. The second component is `'axis'`. The third component is either `'vector'` or `'offset'` for information drawn from the `AXIS.VECTOR[...]` or `AXIS.OFFSET[...]` respectively. The last component is the CBF row number to facilitate recovering the original CBF layout.

4.13 The DIFFRN_DATA_FRAME category

Data items in the `DIFFRN_DATA_FRAME` category record the details about each frame of data.

The items in this category were once in a `DIFFRN_FRAME_DATA` category, which is now deprecated. The items from the old category are provided as aliases but should not be used for new work.

The mapping from CBF into NeXus for the `DIFFRN_DATA_FRAME` category is:

```
_diffrn_data_frame.array_id ARRAYID
_diffrn_data_frame.array_section_id SECTIONID
_diffrn_data_frame.binary_id BINID
_diffrn_data_frame.center_fast CENF
_diffrn_data_frame.center_slow CENS
_diffrn_data_frame.center_units UNITS
_diffrn_data_frame.detector_element_id ELEMENTID
_diffrn_data_frame.id FRAMEID
_diffrn_data_frame.details DETAILS
```

→


```

/entry:NXentry
/instrument:NXinstrument
/DETECTORNAME:NXdetector_group
/DETECTORELEMENTNAME:NXdetector
  /CBF_diffraction_data_frame__section_id=[SECTIONIDARRAY]
  /CBF_diffraction_data_frame__binary_id=[BINARYIDARRAY]
  /CBF_diffraction_data_frame__center_fast_slow=[CENTERARRAY]
    @units="UNITS"
  /CBF_diffraction_data_frame__details=["DETAILSARRAY"]

```

inserts either ARRAYID (if SECTIONID is not specified or SECTIONID into the element of SECTIONIDARRAY for this frame and for this detector element (see below);

inserts BINID into the element of BINARYIDARRAY for this frame and for this detector element (see below);

inserts CENF into the element of CENTERARRAY for this frame, for this detector element and for the fast centre (see below);

inserts CENS into the element of CENTERARRAY for this frame, for this detector element and for the slow centre (see below);

only one CENTERARRAY unit is provided. If there is variation, the values in CENTERARRAY should be rescaled to uniform units.

_diffraction_data_frame.detector_element_id ELEMENTID -->
ELEMENTID used to index into the arrays of this category by the ordinal of the matching ELEMENTID in DIFFRACTION_DETECTOR_ELEMENT__id for the fast index;

FRAMEID used to index into the arrays of this category by the ordinal of the matching ELEMENTID in DIFFRACTION_DETECTOR_ELEMENT__id for the slow index by matching FRAMEID against _diffraction_scan_frame.frame_id and using _diffraction_scan_frame.frame_number from the same row.

inserts DETAILS into the element of DETAILSARRAY for this frame and for this detector element (see below);

The arrays created in the mapping have a slow index of the number of frames and a fast index of the number of detector elements. There is a middle index for CENTERARRAY in the order fast and then slow.

4.14 The DIFFRACTION_DETECTOR category

Data items in the DIFFRACTION_DETECTOR category describe the detector used to measure the scattered radiation, including any analyser and post-sample collimation.

The mapping from CBF into NeXus for the DIFFRACTION_DETECTOR category is:

```

_diffraction_detector.diffraction_id DIFFRACTIONID
_diffraction_detector.id DETECTORNAME

```

```

_diffrn_detector.details DETAILS
_diffrn_detector.detector DETECTOR
_diffrn_detector.dtime DTIME
_diffrn_detector.gain_setting GAINSETTING
_diffrn_detector.number_of_axes NAXES
_diffrn_detector.type DETTYPE

```

→

```

/entry:NXentry
  /CBF_scan_id="SCANID"
  /CBF_diffrn_id="DIFFRNID"
  /instrument:NXinstrument
    /DETECTORNAME:NXdetector_group
    /DETECTORELEMENTNAME:NXdetector
      /details="DETAILS"
      /type="DETECTOR"
      /deadtime=DTIME
      /number_of_axes=NAXES
      /description="DETTYPER"
      /gain_setting="GAINSETTING"

```

4.15 The DIFFRN_DETECTOR_AXIS category

Data items in the DIFFRN_DETECTOR_AXIS category associate axes with detectors.

The mapping from CBF into NeXus for the DIFFRN_DETECTOR_AXIS category is:

```

_diffrn_detector_axis.axis_id AXISID -->
_diffrn_detector_axis.detector_id DETECTORNAME -->
  /instrument:NXinstrument
    /DETECTORNAME:NXdetector_group
    /DETECTORELEMENTNAME:NXdetector
      /transformations:NXtransformations
      /AXISID=[]

```

This information normally will duplicate information obtained from the ARRAY_STRUCTURE_LIST_AXIS.

4.15.1 DIFFRN_DETECTOR_ELEMENT category

Data items in the DIFFRN_DETECTOR_ELEMENT category record the details about spatial layout and other characteristics of each element of a detector which may have multiple elements.

In most cases, giving more detailed information in ARRAY_STRUCTURE_LIST and ARRAY_STRUCTURE_LIST_AXIS is preferable to simply providing the centre of the detector element.

4.16 The DIFFRN_MEASUREMENT category

Data items in the DIFFRN_MEASUREMENT category record details about the device used to orient and/or position the crystal during data measurement and the manner in which the diffraction data

were measured.

The mapping from CBF into NeXus for the DIFFRN_MEASUREMENT category is:

```
_diffrn_measurement.diffrn_id DIFFRNID
_diffrn_measurement.details DETAILS
_diffrn_measurement.device DEVICE
_diffrn_measurement.device_details DEVDETAILS
_diffrn_measurement.device_type DEVTYPE
_diffrn_measurement.id GONIOMETER
_diffrn_measurement.method METHOD
_diffrn_measurement.number_of_axes NUMBER
_diffrn_measurement.sample_detector_distance DIST
_diffrn_measurement.sample_detector_distance_derived DISTDERIVED
_diffrn_measurement.sample_detector_voffset VOFST
_diffrn_measurement.specimen_support SPECSPT
```

→

```
/entry:NXentry
/ CBF_scan_id="SCANID"
/ CBF_diffrn_id="DIFFRNID"
/ instrument:NXinstrument
/ CBF_diffrn_measurement__GONIOMETER:NXgoniometer
/ details="DETAILS"
/ local_name="DEVICE"
/ description="DEVDETAILS"
/ type="DEVTYPE"
/ CBF_diffrn_measurement__method="METHOD"
/ number_of_axes=NUMBER
/ CBF_diffrn_measurement__specimen_support="SPECSPT"
/ CBF_diffrn_detector__DETECTORNAME:NXdetector
/ distance=DIST
@units="mm"
/ distance_derived=DISTDERIVED
/ CBF_diffrn_measurement__sample_detector_voffset=VOFST
@units="mm"
```

4.17 The DIFFRN_MEASUREMENT_AXIS category

Data items in the DIFFRN_MEASUREMENT_AXIS category associate axes with goniometers.

The mapping from CBF into NeXus for the DIFFRN_MEASUREMENT_AXIS category is:

```
_diffrn_measurement_axis.axis_id AXISID
_diffrn_measurement_axis.measurement_device DEVICE
_diffrn_measurement_axis.measurement_id GONIOMETER
```

→

```

/entry:NXentry
  /CBF_scan_id="SCANID"
  /CBF_diffn_id="DIFFRNID"
  /sample:NXsample
    /transformations:NXtransformations
      /AXISID=[]
  /instrument:NXinstrument
    /CBF_diffn_measurement__GONIOMETER:NXgoniometer
      /CBF_diffn_measurement__device="DEVICE"

```

4.18 The DIFFRN_RADIATION category

Data items in the DIFFRN_RADIATION category describe the radiation used for measuring diffraction intensities, its collimation and monochromatization before the sample.

Post-sample treatment of the beam is described by data items in the DIFFRN_DETECTOR category.

The mapping from CBF into NeXus for the DIFFRN_RADIATION category is:

_diffn_radiation.collimation	COLLIMATION
_diffn_radiation.diffn_id	DIFFRNID
_diffn_radiation.div_x_source	DIVX
_diffn_radiation.div_y_source	DIVY
_diffn_radiation.div_x_y_source	DIVXY
_diffn_radiation.filter_edge'	ABSEEDGE
_diffn_radiation.inhomogeneity	HWIDTH
_diffn_radiation.monochromator	MONOCHROMATOR
_diffn_radiation.polarisn_norm	POLNANG
_diffn_radiation.polarisn_ratio	POLRAT
_diffn_radiation.polarizn_source_norm	POLSNANG
_diffn_radiation.polarizn_source_ratio	POLSRAT
_diffn_radiation.polarizn_Stokes_I	SVECI
_diffn_radiation.polarizn_Stokes_Q	SVECQ
_diffn_radiation.polarizn_Stokes_U	SVECU
_diffn_radiation.polarizn_Stokes_V	SVECV
_diffn_radiation.polarisn_norm_esd	POLNANGESD
_diffn_radiation.polarisn_ratio_esd	POLRATESD
_diffn_radiation.polarizn_source_norm_esd	POLSNANGESD
_diffn_radiation.polarizn_source_ratio_esd	POLSRATESD
_diffn_radiation.polarizn_Stokes_I_esd	SVECIESD
_diffn_radiation.polarizn_Stokes_Q_esd	SVECQESD
_diffn_radiation.polarizn_Stokes_U_esd	SVECUESD
_diffn_radiation.polarizn_Stokes_V_esd	SVECVESD
_diffn_radiation.probe	RADIATION
_diffn_radiation.type	SIEGBAHNTYPE
_diffn_radiation.xray_symbol	IUPACXRAYSYMB
_diffn_radiation.wavelength_id	ID
_diffn_radiation_wavelength.id	ID
_diffn_radiation_wavelength.wavelength	WAVELENGTH
_diffn_radiation_wavelength.wt	WEIGHT

```

_diffn_scan_frame.polarizn_Stokes_I STOKESI
_diffn_scan_frame.polarizn_Stokes_Q STOKESQ
_diffn_scan_frame.polarizn_Stokes_U STOKESU
_diffn_scan_frame.polarizn_Stokes_V STOKESV
_diffn_scan_frame.polarizn_Stokes_I_esd STOKESIESD
_diffn_scan_frame.polarizn_Stokes_Q_esd STOKESQESD
_diffn_scan_frame.polarizn_Stokes_U_esd STOKESUESD
_diffn_scan_frame.polarizn_Stokes_V_esd STOKESVESD
_diffn_scan_frame.frame_number FRAMENO

```

→

```

/entry:NXentry
  /CBF_scan_id="SCANID"
  /sample:NXsample
  /beam:NXbeam
    /incident_divergence_x=DIVX
      @units="degrees"
    /incident_divergence_y=DIVY
      @units="degrees"
    /incident_divergence_xy=DIXXY
      @units="degrees^2"
    /CBF_diffn_radiation_wavelength__wavelength_id=[WAVELENGTH_ID]
    /incident_wavelength=[WAVELENGTH]
      @units="A"
    /weight=[WEIGHT]
    /incident_polarisation_stokes_average=[SVECI,SVECQ,SVECU,SVECV]
    /incident_polarisation_stokes_average_uncertainty=[SVECIESD,SVECQESD,SVEQUESD,SVECVESD]
    /incident_polarisation_stokes=[STOKESI,STOKESQ,STOKESU,STOKESV]
    /incident_polarisation_stokes_uncertainty=[STOKESIESD,STOKESQESD,STOKESUESD,STOKESVESD]
      @units="Watts/meter^2"
    /CBF_diffn_radiation__polarisn_norm=POLNANG
      @units="deg"
    /CBF_diffn_radiation__polarisn_ratio=POLRAT
    /CBF_diffn_radiation__polarisn_norm_uncertainty=POLNANGESD
      @units="deg"
    /CBF_diffn_radiation__polarisn_ratio_uncertainty=POLRATESD
    /CBF_diffn_radiation__polarisn_source_norm=POLSNANG
      @units="deg"
    /CBF_diffn_radiation__polarizn_source_ratio=POLSRAT
    /CBF_diffn_radiation__polarisn_source_norm_uncertainty=POLSNANGESD
      @units="deg"
    /CBF_diffn_radiation__polarizn_source_ratio_uncertainty=POLSRATESD
    /CBF_diffn_radiation__filter_edge=ABSEGE
      @units="angstroms"
    /CBF_diffn_radiation__inhomogeneity=HWIDTH
      @units="mm"
  /instrument:NXinstrument
    /monochromator:NXmonochromator
      /description="MONOCHROMATOR"

```

```

/source:NXsource
/probe="RADIATION"
/CBF_diffraction_radiation__type="SIEGBAHNTYPE"
/CBF_diffraction_radiation__xray_symbol="IUPACXRAYSYMB"

```

With the incident_polarisation_stokes array indexed by FRAMENO

```

\end{itemize}
}

```

```

\subsection{The DIFFRN\_RADIATION\_WAVELENGTH category}

```

```

\footnotesize{\begin{verbatim}
The DIFFRN_RADIATION_WAVELENGTH category is defined in
the Protein Data Bank's PDBx/mmCIF dictionary
\cite{westbrook2022pdbx}

```

Data items in the DIFFRN_RADIATION_WAVELENGTH category describe the wavelength of the radiation used to measure the diffraction intensities. Items may be looped to identify and assign weights to distinct components of a polychromatic beam.

- `_diffraction_radiation.wavelength.id` ID →
`/instrument:NXinstrument`
`/CBF_diffraction_radiation_wavelength_id=["ID"]`
- `_diffraction_radiation.wavelength.wavelength_id` WAVELENGTH_ID →
`/instrument:NXinstrument`
`/CBF_diffraction_radiation_wavelength_wavelength_id=["WAVELENGTH_ID"]`
- `_diffraction_radiation.wavelength.wavelength` WAVELENGTH →
`/instrument:NXinstrument`
`/CBF_diffraction_radiation_wavelength_wavelength=[WAVELENGTH]`
- `_diffraction_radiation.wavelength.wt` WEIGHT → `/instrument:NXinstrument`
`/CBF_diffraction_radiation_wavelength_wt=[WEIGHT]`
- `_diffraction_radiation.wavelength.variant` → ?? (*NeXus does not handle variants at this time*)

4.19 The DIFFRN_REFLN category

This category redefinition has been added to extend the key of the standard DIFFRN_REFLN category from mmCIF.

Data items in the DIFFRN_REFLN category record details about the intensities in the diffraction data set identified by `_diffraction_refl.diffraction_id`.

The DIFFRN_REFLN data items refer to individual intensity measurements and must be included in looped lists.

The DIFFRN_REFLNS data items specify the parameters that apply to all intensity measurements in the particular diffraction data set identified by `_diffnr_reflms.diffnr_id` and `_diffnr_reflms.frame_id`

There are proposals on handling reflections in NeXus in NXreflections. When the discussions are resolved this section will be updated.

```
\begin{itemize}
```

```
\item{\_diffnr\_refln.frame\_id $\rightarrow$ ??}
```

```
\item{\_diffnr\_refln.variant $\rightarrow$ ??}
```

```
\end{itemize}
```

```
}
```

```
\subsection{The DIFFRN\_SCAN category}
```

```
\footnotesize{\begin{verbatim}
```

```
Data items in the DIFFRN_SCAN category describe the parameters of one or more scans, relating axis positions to frames.
```

The mapping from CBF into NeXus for the DIFFRN_SCAN category is:

```
_diffnr_scan.id SCANID
_diffnr_scan.date_end ENDDATETIME
_diffnr_scan.date_end_estimated ENDDATETIMEEST
_diffnr_scan.date_start STARTDATETIME
_diffnr_scan.integration_time AVGCOUNTTIME
_diffnr_scan.frame_id_start FRAMESTARTID
_diffnr_scan.frame_id_end FRAMEENDID
_diffnr_scan.frames FRAMES
_diffnr_scan.time_period TIMEPER
_diffnr_scan.time_rstrt_incr RSTRTTIME
```

→

```
/entry:NXentry
/entry_scan_id="SCANID"
/end_time=ENDDATETIME
/end_time_estimated=ENDDATETIMEEST
/start_time=STARTDATETIME
/average_count_time=AVGCOUNTTIME
@units="sec"
/average_frame_time=TIMEPER
```

```

        @units="sec"
        /average_frame_restart_time=RSTRTIME
        @units="sec"
        /instrument:NXinstrument
        /DETECTORNAME:NXdetector_group
        /DETECTORELEMENTNAME:NXdetector
        /frame_start_number=FRAMESTARTNO
        /frame_end_number=FRAMEENDNO

FRAMESTARTNO is the value of _diffrn_scan_frame.frame_number
for which the value of _diffrn_scan_frame.frame_id equals FRAMESTARTID
FRAMEENDNO is the value of _diffrn_scan_frame.frame_number
for which the value of _diffrn_scan_frame.frame_id equals FRAMEENDID
\end{itemize}
}

```

\subsection{The DIFFRN_SCAN_AXIS category}

\footnotesize{\begin{verbatim}

Data items in the DIFFRN_SCAN_AXIS category describe the settings of axes for particular scans. Unspecified axes are assumed to be at their zero points.

The mapping from CBF into NeXus for the DIFFRN_SCAN_AXIS category is:

```

_diffrn_scan_axis.axis_id AXISID-->
_diffrn_scan_axis.angle_start ANGSTART
_diffrn_scan_axis.angle_range ANGRANGE
_diffrn_scan_axis.angle_increment ANGINC
_diffrn_scan_axis.angle_rstrt_incr ANGRSTRT
_diffrn_scan_axis.displacement_start DISPSTART
_diffrn_scan_axis.displacement_range DISPRANGE
_diffrn_scan_axis.displacement_increment DISPINC
_diffrn_scan_axis.displacement_increment DISPINC
_diffrn_scan_axis.displacement_rstrt_incr DISPRSTRT
_diffrn_scan_axis.reference_angle ANG
_diffrn_scan_axis.reference_displacement DISP
_diffrn_scan_axis.scan_id SCANID

```

→

```

{ /entry:NXentry
  /CBF_scan_id="SCANID"
  /instrument:NXinstrument
    /DETECTORNAME:NXdetector_group
    /DETECTORELEMENTNAME:NXdetector
  for AXISEQUIPMENT=="detector"}
{ /CBF_diffrn_scan__SCANID:NXentry
  /sample:NXsample

```



```

for AXISEQUIPMENT=="goniometer"}
{ /CBF_diffraction_scan__SCANID:NXentry
for AXISEQUIPMENT=="general"}
  /transformations:NXtransformations
  /AXISID=[]
    @diffraction_scan_axis__angle_start=ANGSTART
    @diffraction_scan_axis__angle_range=ANGRANGE
    @diffraction_scan_axis__angle_increment=ANGINC
    @diffraction_scan_axis__angle_rstrt_incr=ANGRSTRT
    @diffraction_scan_axis__displacement_start=DISPSTART
    @diffraction_scan_axis__displacement_range=DISPRANGE
    @diffraction_scan_axis__displacement_increment=DISPINC
    @diffraction_scan_axis__displacement_rstrt_incr=DISPRSTRT
    @diffraction_scan_axis__reference_angle=ANG
    @diffraction_scan_axis__reference_displacement=DISP
\end{itemize}
}

```

`\footnotesize{\begin{verbatim}`
 Data items in the DIFFRACTION_SCAN_FRAME category describe
 the relationships of particular frames to scans.

The mapping from CBF into NeXus for the DIFFRACTION_SCAN_FRAME category is:

```

_diffraction_scan_frame.date DATETIME
_diffraction_scan_frame.frame_id ID
_diffraction_scan_frame.frame_number FRAMENUMBER
_diffraction_scan_frame.integration_time COUNTTIME
_diffraction_scan_frame.polarization_Stokes_I STOKESI
_diffraction_scan_frame.polarization_Stokes_Q STOKESQ
_diffraction_scan_frame.polarization_Stokes_U STOKESU
_diffraction_scan_frame.polarization_Stokes_V STOKESV
_diffraction_scan_frame.scan_id SCANID
_diffraction_scan_frame.time_period FRAMETIME
_diffraction_scan_frame.time_rstrt_incr RSTRTTIME

```

→

```

/entry:NXentry
  /CBF_scan_id="SCANID"
  /instrument:NXinstrument
    /DETECTORNAME:NXdetector_group
    /DETECTORELEMENTNAME:NXdetector
      /CBF_diffraction_scan_frame__date=["DATETIME"]
      /CBF_diffraction_scan_frame__frame_id=["ID"]
      /count_time=[COUNTTIME]
      /frame_time=[FRAMETIME]
      /frame_restart_time=[RSTRTTIME]
  /sample:NXsample

```

```

        /beam:NXbeam
        /incident_polarisation_stokes=[STOKESI,STOKESQ,STOKESU,STOKESV]
        @units="Watts/meter^2"
        where each array element is inserted at index FRAMENUMBER
    \end{itemize}
}

```

```

\subsection{The DIFFRN_SCAN_FRAME_AXIS category}

```

```

\footnotesize{\begin{verbatim}
Data items in the DIFFRN_SCAN_FRAME_AXIS category describe the
settings of axes for particular frames. Unspecified axes are
assumed to be at their zero points. If, for any given frame,
nonzero values apply for any of the data items in this category,
those values should be given explicitly in this category and not
simply inferred from values in DIFFRN_SCAN_AXIS.

```

The mapping from CBF into NeXus for the DIFFRN_SCAN_FRAME_AXIS category is:

```

_diffrn_scan_frame_axis.axis_id AXISID
_diffrn_scan_frame_axis.angle ANGLE
_diffrn_scan_frame_axis.angle_increment ANGLEINCREMENT
_diffrn_scan_frame_axis.angle_rstrt_incr ANGLERSTRTINCREMENT
_diffrn_scan_frame_axis.displacement DISP
_diffrn_scan_frame_axis.displacement_increment DISPINCREMENT
_diffrn_scan_frame_axis.displacement_rstrt_incr DISPRSTRTINCREMENT
_diffrn_scan_frame_axis.reference_angle REFANGLE
_diffrn_scan_frame_axis.reference_displacement REFDISP

```

→

```

{ /entry:NXentry
  /CBF_scan_id="SCANID"
  /instrument:NXinstrument
  /DETECTORNAME:NXdetector_group
  /DECTORELEMENTNAME:NXdetector
for AXISEQUIPMENT=="detector"}
{ /entry:NXentry
  /sample:NXsample
for AXISEQUIPMENT=="goniometer"}
{ /entry:NXentry
for AXISEQUIPMENT=="general"}
  /transformations:NXtransformations
  /AXISID=[]
    @diffrn_scan_frame_axis__angle_start=[ANGSTART]
    @diffrn_scan_frame_axis__angle_range=[ANGRANGE]
    @diffrn_scan_frame_axis__angle_increment=[ANGINC]
    @diffrn_scan_frame_axis__angle_rstrt_incr=[ANGRSTRT]
    @diffrn_scan_frame_axis__displacement_start=[DISPSTART]

```

```

@diffn_scan_frame_axis__displacement_range=[DISPRANGE]
@diffn_scan_frame_axis__displacement_increment=[DISPINC]
@diffn_scan_frame_axis__displacement_rstrt_incr=[DISPRSTRT]
@diffn_scan_frame_axis__reference_angle=[ANG]
@diffn_scan_frame_axis__reference_displacement=[DISP]

```

note that @units="mm" or @units="deg" should also be specified.

The dimensions of the array depend on np (the number of frames = the value of _diffn_scan.frames)

either DISP OR ANGLE is inserted as the i-th element counting from 1 in AXISID where i is the value of _diffn_scan.frame.frame_number for which the value of _diffn_scan.frame.frame_id agrees with the value of _diffn_scan.frame.axis.frame_id

```

The remaining tags similarly populate the attribute arrays
\end{itemize}
}

```

```

\subsection{The DIFFRN\_SCAN\_FRAME\_MONITOR category}

```

```

\footnotesize{\begin{verbatim}

```

Data items in the DIFFRN_SCAN_FRAME_MONITOR category record the values and details about each monitor for each frame of data during a scan.

Each monitor value is uniquely identified by the combination of the scan id given by _diffn_scan.frame.scan_id, the frame id given by _diffn_scan.frame_monitor.frame_id, the monitor's detector id given by _diffn_scan.frame_monitor.detector_id, and a 1-based ordinal given by _diffn_scan.frame_monitor.id.

If there is only one frame for the scan, the value of _diffn_scan.frame_monitor.frame_id may be omitted.

A single frame may have more than one monitor value, and each monitor value may be the result of integration over the entire frame integration time given by the value of _diffn_scan.frame.integration_time, or many monitor values may be reported over shorter times given by the value of _diffn_scan.frame_monitor.integration_time. If only one monitor value for a given monitor is collected during the integration time of the frame, the value of _diffn_scan.frame_monitor.id may be omitted.

The mapping from CBF into NeXus for the DIFFRN_SCAN_FRAME_MONITOT category is:

```

_diffrn_scan_frame_monitor.id MONID -->
_diffrn_scan_frame_monitor.detector_id DETECTORNAME -->
_diffrn_scan_frame_monitor.scan_id SCANID -->
_diffrn_scan_frame_monitor.frame_id FRAMEID -->
_diffrn_scan_frame_monitor.integration_time INTEGRATIONTIME -->
_diffrn_scan_frame_monitor.monitor_value MONITORVALUE -->

→

/entry:NXentry
/  CBF_scan_id="SCANID"
/    instrument:NXinstrument
/      CBF_diffrn_scan_frame_monitor__DETECTORNAME_MONID:NXmonitor
/        @CBF_detector_id="DETECTORNAME"
/        @CBF_diffrn_scan_frame_monitor__id="MONID"
/        /data=[MONITORVALUE]
/        /count_time=[INTEGRATIONTIME]
\end{itemize}
}

\subsection{The MAP category}

{\footnotesize\begin{verbatim}
Data items in the MAP category record
the details of a maps. Maps record values of parameters,
such as density, that are functions of position within
a cell or are functions of orthogonal coordinates in
three space.
```

A map may be composed of one or more map segments specified in the MAP_SEGMENT category.

- `_map.details` → ??
- `_map.diffrn_id` → ??
- `_map.entry_id` → ??
- `_map.id` → ??
- `_map.variant category` → ??

To be mapped to NeXus in the future.

4.19.1 The MAP_SEGMENT category

Data items in the MAP_SEGMENT category record the details about each segment (section or brick) of a map.

- `_map_segment.array_id` → ??
- `_map_segment.array_section_id` → ??

- `_map.segment.binary_id` → ??
- `_map.segment.mask.array_id` → ??
- `_map.segment.mask.array.section_id` → ??
- `_map.segment.mask.binary_id` → ??
- `_map.segment.id` → ??
- `_map.segment.map_id` → ??
- `_map.segment.details` → ??
- `_map.segment.variant` category → ??

To be mapped to NeXus in the future.

4.20 The VARIANT category

Data items in the VARIANT category record the details about sets of VARIANTS of data items.

There is sometimes a need to allow for multiple versions of the same data items in order to allow for refinements and corrections to earlier assumptions, observations and calculations. In order to allow data sets to contain more than one VARIANT of the same information, an optional `...variant` data item as a pointer to `_variant.variant` has been added to the key of every category, as an implicit data item with a null (empty) default value.

All rows in a category with the same VARIANT value are considered to be related to one another and to all rows in other categories with the same VARIANT value. For a given VARIANT, all such rows are also considered to be related to all rows with a null VARIANT value, except that a row with a null VARIANT value is for which all other components of its key are identical to those entries in another row with a non-null VARIANT value is not related the the rows with that non-null VARIANT value. This behavior is similar to the convention for identifying alternate conformers in an atom list.

An optional role may be specified for a VARIANT as the value of `_variant.role`. Possible roles are null, "preferred", "raw data", "unsuccessful trial".

VARIANTS may carry an optional timestamp as the value of `_variant.timestamp`.

VARIANTS may be related to other VARIANTS from which they were derived by the value of `_variant.variant_of`

Further details about the VARIANT may be specified as the value of `_variant.details`.

In order to allow VARIANT information from multiple datasets to be combined, `_variant.diffraction_id` and/or `_variant.entry_id` may be used.

- `_variant.details` → ??
- `_variant.role` → ??
- `_variant.timestamp` → ??
- `_variant.variant` → ??
- `_variant.variant_of` → ??

To be mapped to NeXus in the future.

5 Mapping from NeXus to CBF/imgCIF

Mapping from NeXus to CBF/imgCIF is matter of establishing appropriate tables and columns for each of the NeXus classes. In general, a NeXus class will correspond to a category, while the specific name will be part of the key of that category, usually the id of the category. The terms in NeXus may be defined in base classes (see http://download.nexusformat.org/doc/html/classes/base_classes/index.html) or in application definitions (see <http://download.nexusformat.org/doc/html/classes/applications/index.html>). We will consider them one at a time.

As with the mapping from CBF to NeXus, where we do both a complete faithful mapping of a CBF to a NeXus tree in CBF.cbf, and a more structured tag-by-tag mapping throughout the NeXus tree, in this direction we will address a complete faithful mapping of all features of a NeXus tree into a single CBF data block as well as a more structured item-by-item mapping.

5.1 Faithful Mapping of a NeXus Tree

A NeXus tree is a tree of HDF5 groups, each with an associated NeXus class. For the full faithful mapping in to a single CBF datablock, which will be given the name `NeXus_Tree`, each NeXus class will be mapped to a CBF category of the same name, and each instance of that class will be mapped to a single row in that table. The unique identifier of a row, in the column `NX_tree.path` will be the rooted path to the particular NeXus class instance in the tree. Each path component will consist of the NeXus class composed with the actual name of the NeXus class instance in dotted notation, with an underscore before each NeXus class name. Each subgroup of a NeXus class instance (also a NeXus class instance) will be entered into a column named with name of the subgroup, with a value equal to the rooted path to that subgroup.

For example, in the NeXus tree

```
\entry:NXentry
  \instrument:NXinstrument
    \detector:NXdetector
```

the mapping would be to

```
_datablock_NeXus_Tree

loop_
  _NXentry.NX_tree_path  "/_NXentry.entry"
  _NXentry.instrument
```

```

"/_NXentry.entry"    "/_NXentry.entry/_NXinstrument.instrument"

loop_
_NXinstrument.NX_tree_path
_NXinstrument.detector
"/_NXentry.entry/_NXinstrument.instrument" "/_NXentry.entry/_NXinstrument.instrument/_NXdetector__detector"

_NXdetector.detector "/_NXentry.entry/_NXinstrument.instrument/_NXdetector__detector"

```

An attribute of a NeXus class other than the HDF5 attribute for the NeXus class name itself will be given a column name composed from the prefix ‘NX.class_attribute_’ and the name of the attribute.

5.2 Mapping Fields

Fields share some of the characteristics of classes and some of the characteristics of attributes. When a rooted path to a field is needed in this discussion, the rooted path be composed of the NeXus class composed with the name of the NeXus class instance composed with the name of the field in dotted notation, with an underscore before each NeXus class name.

If a NeXus class instance is given in the NeXus tree by a link, rather than directly an extra tag having the name of the column with the suffix ‘__link’ will be used to carry the link path to the target in addition to giving the path to the origin of the link.

The mapping of fields in NeXus class instances will be handled in one of two ways:

- If instances of a NeXus class are each permitted to contain multiple instances of a field, each field instance will be handled like a NeXus class instance, with the modification to rooted paths noted above. Each such field will be assigned a column name that is generic to the instances, such as ‘data’ or ‘axis_poise’, and a category name beginning with ‘NX_’ will be generated to hold the field instance value and attribute values. The value in the column will be a list of the rooted paths to the field instances.
- if each instance of a NeXus class may only contain one unique instance of a field, then the field name will be used as the column name. The value of the column for a field will be the verbatim value of the field, using the CBF binary data type when needed, or CIF 2 lists and tuples when needed. Attributes of fields are handled by creating a column beginning with the name of the column of the field and appending the attribute name in dotted notation, converting the dots to double underscores.

For example, in the NeXus tree

```

\entry:NXentry
  \instrument:NXinstrument
    \detector:NXdetector
      \data=BINDATA
      \@signal=1

```

the mapping would be to

```
_datablock_NeXus_Tree
```

```
_loop_
```

```

_NXentry.NX_tree_path  "/_NXentry.entry"
_NXentry.instrument
"/_NXentry.entry"    "/_NXentry.entry/_NXinstrument.instrument"

loop_
_NXinstrument.NX_tree_path
_NXinstrument.detector
"/_NXentry.entry/_NXinstrument.instrument" "/_NXentry.entry/_NXinstrument.instrument/_NXdetector__detector"

loop_
_NXdetector.NX_tree_path
_NXdetector.data      BINDATA
_NXdetector.data__signal

"/_NXentry.entry/_NXinstrument.instrument/_NXdetector__detector"  BINDATA  1

```

While, in the NeXus tree

```

/entry:NXentry
/instrument:NXinstrument
/CBF\_diffraction\_detector__detector:NXdetector
  CBF\_axis\_DETECTOR\_PITCH=[0.]
    @units="deg"
    @cbf_location="image_1.axis.vector.10"
    @depends_on="CBF\_axis\_DETECTOR\_Y"
    @transformation_type="rotation"
    @vector=[-1, 0, 0]
  CBF\_axis\_DETECTOR\_Y=[0.]
    @units="mm"
    @cbf_location="image_1.axis.vector.9"
    @depends_on="CBF\_axis\_DETECTOR\_Z"
    @transformation_type="translation"
    @vector=[0, -1, 0]
  CBF\_axis\_DETECTOR\_Z=[250.]
    @units="mm"
    @cbf_location="image_1.axis.vector.8"
    depends_on="."
    @transformation_type="translation"
    @vector=[0, 0, 1]
  CBF\_axis\_ELEMENT\_X=[0.]
    @units="mm"
    @offset_units="mm"
    @cbf_location="image_1.axis.vector.11"
    @depends_on="CBF\_axis\_DETECTOR\_PITCH"
    @transformation_type="translation"
    @vector=[-1, 0, 0]
    @offset=[-211.818, -217.322, 0]
  CBF\_axis\_ELEMENT\_Y=[0.]
    @units="mm"

```



```

        @cbf_location=image_1.axis.vector.12
        @depends_on="CBF_axis__ELEMENT_X"
        @transformation_type="translation"
        @vector= [0, 1, 0]

the mapping would be to

_datablock_NeXus_Tree

_loop_
_NXentry.NX_tree_path  "/"_NXentry.entry"
_NXentry.instrument
"/_NXentry.entry"      "/"_NXentry.entry/_NXinstrument.instrument"

loop_
_NXinstrument.NX_tree_path
_NXinstrument.detector
"/_NXentry.entry/_NXinstrument.instrument" "/"_NXentry.entry/_NXinstrument.instrument/_NXdetector__diffrn_detect

loop_
_NXdetector.NX_tree_path
_NXdetector.axis_poise

"/_NXentry.entry/_NXinstrument.instrument/_NXdetector__diffrn_detector__detector__axis_poise"
["CBF_axis__DETECTOR_PITCH", "CBF_axis__DETECTOR_Y", "CBF_axis__DETECTOR_Z", "CBF_axis__ELEMENT_X", "CBF_axis__ELEME

loop_
_NX_axis_poise.NX_tree_path
_NX_axis_poise.value
_NX_axis_poise.units
_NX_axis_poise.offset_units
_NX_axis_poise.cbf_location
_NX_axis_poise.depends_on
_NX_axis_poise.transformation_type
_NX_axis_poise.vector
_NX_axis_poise.offset

"/_NXentry.entry/_NXinstrument.instrument/_NXdetector__diffrn_detector__detector__axis_poise__axis__DETECTOR_PITCH"
0. "deg" . "image_1.axis.vector.10" "axis__DETECTOR_Y" "rotation" [-1,0,0] .
"/_NXentry.entry/_NXinstrument.instrument/_NXdetector__diffrn_detector__detector__axis_poise__axis__DETECTOR_Y"
0. "mm" . "image_1.axis.vector.9" "axis__DETECTOR_Z" "translation" [0,-1,0] .
"/_NXentry.entry/_NXinstrument.instrument/_NXdetector__diffrn_detector__detector__axis_poise__axis__DETECTOR_Z"
250. "mm" . "image_1.axis.vector.8" . "translation" [0,0,1] .
"/_NXentry.entry/_NXinstrument.instrument/_NXdetector__diffrn_detector__detector__axis_poise__axis__ELEMENT_X"
0. "mm" "mm" "image_1.axis.vector.11" "axis__DETECTOR_PITCH" "translation" [-1,0,0] [-211.818,-217.322,0]
"/_NXentry.entry/_NXinstrument.instrument/_NXdetector__diffrn_detector__detector__axis_poise__axis__ELEMENT_Y"
0. "mm" . "image_1.axis.vector.12" "axis__ELEMENT_X" "translation" [0,1,0] .

```

While the full faithful mapping of the NeXus tree into CBF is primarily for development and debugging, portions of it will also become essentially the final detailed mapping. Notice, for

example, that the `NX.axis.poise` category in the faithful mapping has all the information needed to populate the standard CBF axis category with minimal transformations. In the final detailed mapping, in addition to the `NX.tree.path`, each category will have a ‘`NX_id`’ column to provide a unique identifier when a CBF originates the data and no path is available to use as the key. Normally this will be the last component of the path.

6 NeXus Base Class Mapping

The details of the NeXus base classes are provided at
http://download.nexusformat.org/doc/html/classes/base_classes/index.html
 which should be consulted for getails.

6.1 NXaperture

Template of a beamline aperture.

NXaperture (base class, version 1.0)

```
description:NX_CHAR
material:NX_CHAR
NXgeometry
NXgeometry
NXnote
```

- APERTURE:NXaperture →
 - _NXaperture.NX_tree_path NEXUSTREEPATH
 - _NXaperture.NX_id APERTURE
 - _NXaperture.NX_scan_id SCANID
 - _NXaperture.NX_diffn_id DIFFRNID
 - _NXaperture.NX_entry_id ENTRYID
- description:NX_CHAR=DESCRIPTION →
 - _NXaperture.description DESCRIPTION
- material:NX_CHAR=MATERIAL →
 - _NXaperture.material MATERIAL
- geometry_aperture:NXgeometry →
 - _NXaperture.NXgeometry_id geometry_aperture
- geometry_blades:NXgeometry →
 - _NXaperture.NXgeometry_id geometry.blades
- note:NXnote →
 - _NXaperture.NXnote_id note

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘`/NXaperture__APERTURE`’ where APERTURE is the name of this group, typically ‘`aperture`’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

6.2 NXattenuator

Description of a device that reduces the intensity of a beam by attenuation. If uncertain whether to use `NXfilter` (band-pass filter) or `NXattenuator` (reduces beam intensity), then choose `NXattenuator`.

NXattenuator (base class, version 1.0)

```
absorption_cross_section:NX_FLOAT
attenuator_transmission:NX_FLOAT
distance:NX_FLOAT
scattering_cross_section:NX_FLOAT
status:NX_CHAR
    @time
thickness:NX_FLOAT
type:NX_CHAR
```

- ATTENUATOR:NXattenuator →
 _NXattenuator.NX_tree_path NEXUSTREEPATH
 _NXattenuator.NX_id ATTENUATOR
 _NXattenuator.NX_scan_id SCANID
 _NXattenuator.NX_diffn_id DIFFRNID
 _NXattenuator.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXattenuator_’ where ATTENUATOR is the name of this group, typically ‘attenuator’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- absorption_cross_section:NX_FLOAT=ABSORPTION_CROSS_SECTION →
 _NXattenuator.absorption_cross_section ABSORPTION_CROSS_SECTION
- attenuator_transmission:NX_FLOAT=ATTENUATOR_TRANSMISSION →
 _NXattenuator.attenuator_transmission ATTENUATOR_TRANSMISSION
- distance:NX_FLOAT=DISTANCE →
 _NXattenuator.distance DISTANCE
- scattering_cross_section:NX_FLOAT=SCATTERING_CROSS_SECTION →
 _NXattenuator.scattering_cross_section SCATTERING_CROSS_SECTION
- status:NX_CHAR=STATUS →
 _NXattenuator.status STATUS
- @time=TIME →
 _NXattenuator.status_time TIME
- thickness:NX_FLOAT=THICKNESS →
 _NXattenuator.thickness THICKNESS
- type:NX_CHAR=TYPE →
 _NXattenuator.type TYPE

6.3 NXbeam

Template of the state of the neutron or X-ray beam at any location. It will be referenced by beamline component groups within the NXinstrument group or by the NXsample group. Note that variables such as the incident energy could be scalar values or arrays. This group is especially valuable in storing the results of instrument simulations in which it is useful to specify the beam profile, time distribution *etc.* at each beamline component. Otherwise, its most likely use is in the NXsample group in which it defines the results of the neutron scattering by the sample, *e.g.*, energy transfer, polarizations.

Note: there has been a recent NIAC discussion on this class, in which Tobias Richter pointed out that the polarization is unspecified, and the Stokes parameters would be 4 array position, not 2. In addition, there was a ‘j’ index floating around and the 2-dimensional arrays appear to have been transposed. The version has added the 2-parameter Denzo polarization parameters in the order norm, then ratio. In addition we have added the Stokes parameters in the order I, Q, U, V. We propose that the unspecified polarization be deprecated.

```
NXbeam (base class, version 1.0)
  distance:NX_FLOAT
  energy_transfer:NX_FLOAT[i]
  final_beam_divergence:NX_FLOAT[i,2]
  final_energy:NX_FLOAT[i]
  final_polarization:NX_FLOAT[i,2]
  final_polarization_Denzo:NX_FLOAT[i,2]
  final_polarization_Stokes:NX_FLOAT[i,4]
  final_wavelength:NX_FLOAT[i]
  final_wavelength_spread:NX_FLOAT[i]
  flux:NX_FLOAT[i]
  incident_beam_divergence:NX_FLOAT[i,2]
  incident_energy:NX_FLOAT[i]
  incident_polarization:NX_FLOAT[i,2]
  incident_polarization_Denzo:NX_FLOAT[i,2]
  incident_polarization_Stokes:NX_FLOAT[i,4]
  incident_wavelength:NX_FLOAT[i]
  incident_wavelength_spread:NX_FLOAT[i]
NXdata
```

- BEAM:NXbeam →
 _NXbeam.NX_tree_path NEXUSTREEPATH
 _NXbeam.NX_id BEAM
 _NXbeam.NX_scan_id SCANID
 _NXbeam.NX_diffraction_id DIFFRACTIONID
 _NXbeam.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXbeam_BEAM’ where BEAM is the name of this group, typically ‘beam’. The SCANID, DIFFRACTIONID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- distance:NX_FLOAT=DISTANCE →
 _NXbeam.distance DISTANCE
 is an alias for
 _diffraction_measurement.sample_detector_distance DISTANCE
- energy_transfer:NX_FLOAT[i]=ENERGY_TRANSFER →
 _NXbeam.energy_transfer ENERGY_TRANSFER
- final_beam_divergence:NX_FLOAT[i,2]=FINAL_BEAM_DIVERGENCE →
 _NXbeam.final_beam_divergence FINAL_BEAM_DIVERGENCE
- final_energy:NX_FLOAT[i]=FINAL_ENERGY →
 _NXbeam.final_energy FINAL_ENERGY
- final_polarization:NX_FLOAT[i,2]=FINAL_POLARIZATION →
 _NXbeam.final_polarization FINAL_POLARIZATION

- `final_polarization_Denzo:NX_FLOAT[i,2]=FINAL_POLARIZATION_DENZO →`
`_NXbeam.final_polarization FINAL_POLARIZATION_DENZO`
- `final_polarization_Stokes:NX_FLOAT[i,4]=FINAL_POLARIZATION_STOKES →`
`_NXbeam.final_polarization FINAL_POLARIZATION_STOKES`
- `final_wavelength:NX_FLOAT[i]=FINAL_WAVELENGTH →`
`_NXbeam.final_wavelength FINAL_WAVELENGTH`
- `final_wavelength_spread:NX_FLOAT[i]=FINAL_WAVELENGTH_SPREAD →`
`_NXbeam.final_wavelength_spread FINAL_WAVELENGTH_SPREAD`
- `flux:NX_FLOAT[i]=FLUX →`
`_NXbeam.flux FLUX`
- `incident_beam_divergence:NX_FLOAT[i,2]=INCIDENT_BEAM_DIVERGENCE →`
`_NXbeam.incident_beam_divergence INCIDENT_BEAM_DIVERGENCE`
- `incident_energy:NX_FLOAT[i]=INCIDENT_ENERGY →`
`_NXbeam.incident_energy INCIDENT_ENERGY`
- `incident_polarization:NX_FLOAT[i,2]=INCIDENT_POLARIZATION →`
`_NXbeam.incident_polarization INCIDENT_POLARIZATION`
- `incident_polarization_Denzo:NX_FLOAT[i,4]=INCIDENT_POLARIZATION_DENZO →`
`_NXbeam.incident_polarization_Stokes INCIDENT_POLARIZATION_DENZO`
- `incident_polarization_Stokes:NX_FLOAT[i,4]=INCIDENT_POLARIZATION_STOKES →`
`_NXbeam.incident_polarization_Stokes INCIDENT_POLARIZATION_STOKES`
- `incident_wavelength:NX_FLOAT[i]=INCIDENT_WAVELENGTH →`
`_NXbeam.incident_wavelength INCIDENT_WAVELENGTH`
- `incident_wavelength_spread:NX_FLOAT[i]=INCIDENT_WAVELENGTH_SPREAD →`
`_NXbeam.incident_wavelength_spread INCIDENT_WAVELENGTH_SPREAD`
- `data:NXdata →`
`_NXbeam.NXdata_id data`

The final Denzo polarization from the beam component immediately prior to the beam being incident on the sample should agree with the values of the CBF tags `_diffrn_radiation.polarizn_source_norm` and `_diffrn_radiation.polarizn_source_ratio`.

6.4 NXbeam_stop

A class for a beamstop. Beamstops and their positions are important for SANS and SAXS experiments.

`NXbeam_stop` (base class, version 1.0)

```
description:NX_CHAR
distance_to_detector:NX_FLOAT
size:NX_FLOAT
status:NX_CHAR
x:NX_FLOAT
y:NX_FLOAT
NXgeometry
```

- `BEAM_STOP:NXbeam_stop →`
`_NXbeam_stop.NX_tree_path NEXUSTREEPATH`

```
_NXbeam_stop.NX_id BEAM_STOP
_NXbeam_stop.NX_scan_id SCANID
_NXbeam_stop.NX_diffraction_id DIFFRANID
_NXbeam_stop.NX_entry_id ENTRYID
```

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘‘/NXbeam_stop_BEAM_STOP’’ where BEAM_STOP is the name of this group, typically ‘‘beam_stop’’. The SCANID, DIFFRANID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- description:NX_CHAR=DESCRIPTION →
_NXbeam_stop.description DESCRIPTION
- distance_to_detector:NX_FLOAT=DISTANCE_TO_DETECTOR →
_NXbeam_stop.distance_to_detector DISTANCE_TO_DETECTOR
- size:NX_FLOAT=SIZE →
_NXbeam_stop.size SIZE
- status:NX_CHAR=STATUS →
_NXbeam_stop.status STATUS
- x:NX_FLOAT=X →
_NXbeam_stop.x X
- y:NX_FLOAT=Y →
_NXbeam_stop.y Y
- geometry1:NXgeometry →
_NXbeam_stop.NXgeometry_id geometry1

6.5 NXbending magnet

description for a bending magnet

```
NXbending_magnet (base class, version 1.0)
accepted_photon_beam_divergence:NX_FLOAT
bending_radius:NX_FLOAT
critical_energy:NX_FLOAT
divergence_x_minus:NX_FLOAT
divergence_x_plus:NX_FLOAT
divergence_y_minus:NX_FLOAT
divergence_y_plus:NX_FLOAT
magnetic_field:NX_FLOAT
source_distance_x:NX_FLOAT
source_distance_y:NX_FLOAT
spectrum:NXdata
NXgeometry
```

- BENDING_MAGNET:NXbending_magnet →
_NXbending_magnet.NX_tree_path NEXUSTREEPATH
_NXbending_magnet.NX_id BENDING_MAGNET
_NXbending_magnet.NX_scan_id SCANID

- ```

_NXbending_magnet.NX.diffraction_id DIFFRACTION_ID
_NXbending_magnet.NX.entry_id ENTRY_ID

```
- where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘‘/NXbending\_magnet\_BENDING\_MAGNET’’ where BENDING\_MAGNET is the name of this group, typically ‘‘bending\_magnet’’. The SCANID, DIFFRACTION\_ID and ENTRY\_ID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.
- `accepted_photon_beam_divergence:NX_FLOAT=ACCEPTED_PHOTON_BEAM_DIVERGENCE →`  
`_NXbending_magnet.accepted_photon_beam_divergence ACCEPTED_PHOTON_BEAM_DIVERGENCE`
  - `bending_radius:NX_FLOAT=BENDING_RADIUS →`  
`_NXbending_magnet.bending_radius BENDING_RADIUS`
  - `critical_energy:NX_FLOAT=CRITICAL_ENERGY →`  
`_NXbending_magnet.critical_energy CRITICAL_ENERGY`
  - `divergence_x_minus:NX_FLOAT=DIVERGENCE_X_MINUS →`  
`_NXbending_magnet.divergence_x_minus DIVERGENCE_X_MINUS`
  - `divergence_x_plus:NX_FLOAT=DIVERGENCE_X_PLUS →`  
`_NXbending_magnet.divergence_x_plus DIVERGENCE_X_PLUS`
  - `divergence_y_minus:NX_FLOAT=DIVERGENCE_Y_MINUS →`  
`_NXbending_magnet.divergence_y_minus DIVERGENCE_Y_MINUS`
  - `divergence_y_plus:NX_FLOAT=DIVERGENCE_Y_PLUS →`  
`_NXbending_magnet.divergence_y_plus DIVERGENCE_Y_PLUS`
  - `magnetic_field:NX_FLOAT=MAGNETIC_FIELD →`  
`_NXbending_magnet.magnetic_field MAGNETIC_FIELD`
  - `source_distance_x:NX_FLOAT=SOURCE_DISTANCE_X →`  
`_NXbending_magnet.source_distance_x SOURCE_DISTANCE_X`
  - `source_distance_y:NX_FLOAT=SOURCE_DISTANCE_Y →`  
`_NXbending_magnet.source_distance_y SOURCE_DISTANCE_Y`
  - `spectrum:NXdata →`  
`_NXbending_magnet.NXdata_id spectrum`
  - `geometry1:NXgeometry →`  
`_NXbending_magnet.NXgeometry_id geometry1`

## 6.6 NXcapillary

This is a dictionary of field names to use for describing a capillary as used in X-ray beamlines. Based on information provided by Gerd Wellenreuther.

```

NXcapillary (base class, version 1.0)
accepting_aperture:NX_FLOAT
focal_size:NX_FLOAT
manufacturer:NX_CHAR
maximum_incident_angle:NX_FLOAT
type:NX_CHAR
working_distance:NX_FLOAT
gain:NXdata
transmission:NXdata

```

- CAPILLARY:NXcapillary →  
     \_NXcapillary.NX\_tree\_path NEXUSTREEPATH  
     \_NXcapillary.NX\_id CAPILLARY  
     \_NXcapillary.NX\_scan\_id SCANID  
     \_NXcapillary.NX\_diffn\_id DIFFRNID  
     \_NXcapillary.NX\_entry\_id ENTRYID  
     where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘‘/NXcapillary\_CAPILLARY’’ where CAPILLARY is the name of this group, typically ‘‘capillary’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.
- accepting\_aperture:NX\_FLOAT=ACCEPTING\_APERTURE →  
     \_NXcapillary.accepting\_aperture ACCEPTING\_APERTURE
- focal\_size:NX\_FLOAT=FOCAL\_SIZE →  
     \_NXcapillary.focal\_size FOCAL\_SIZE
- manufacturer:NX\_CHAR=MANUFACTURER →  
     \_NXcapillary.manufacturer MANUFACTURER
- maximum\_incident\_angle:NX\_FLOAT=MAXIMUM\_INCIDENT\_ANGLE →  
     \_NXcapillary.maximum\_incident\_angle MAXIMUM\_INCIDENT\_ANGLE
- type:NX\_CHAR=TYPE →  
     \_NXcapillary.type TYPE
- working\_distance:NX\_FLOAT=WORKING\_DISTANCE →  
     \_NXcapillary.working\_distance WORKING\_DISTANCE
- gain:NXdata →  
     \_NXcapillary.NXdata\_id gain
- transmission:NXdata →  
     \_NXcapillary.NXdata\_id transmission

## 6.7 NXcharacterization

note: This base class may be removed in future releases of NXDL. If you have a use for this base class, please provide a description of your intended use to the NIAC ([nexus-committee@nexusformat.org](mailto:nexus-committee@nexusformat.org)).

NXcharacterization (base class, version 1.0)

```
@source
@location
@mime_type
definition:NX_CHAR
@version
@URL
```

- CHARACTERIZATION:NXcharacterization →  
     \_NXcharacterization.NX\_tree\_path NEXUSTREEPATH  
     \_NXcharacterization.NX\_id CHARACTERIZATION  
     \_NXcharacterization.NX\_scan\_id SCANID  
     \_NXcharacterization.NX\_diffn\_id DIFFRNID  
     \_NXcharacterization.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under



score and, finally, the name of the group, ending with ‘‘/NXcharacterization\_\_CHARACTERIZATION’’ where CHARACTERIZATION is the name of this group, typically ‘‘characterization’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- @source=SOURCE →  
\_NXcharacterization.NX\_class\_attribute\_\_source SOURCE
- @location=LOCATION →  
\_NXcharacterization.NX\_class\_attribute\_\_location LOCATION
- @mime\_type=MIME.TYPE →  
\_NXcharacterization.NX\_class\_attribute\_\_mime\_type MIME.TYPE
- definition:NX.CHAR=DEFINITION →  
\_NXcharacterization.definition DEFINITION
- @version=VERSION →  
\_NXcharacterization.definition\_\_version VERSION
- @URL=URL →  
\_NXcharacterization.definition\_\_URL URL

## 6.8 NXcollection

Use NXcollection to gather together any set of terms. The original suggestion is to use this as a container class for the description of a beamline.

For NeXus validation, NXcollection will always generate a warning since it is always an optional group. Anything (groups, fields, or attributes) placed in an NXcollection group will not be validated.

NXcollection (contributed definition, version 1.0)  
beamline:NX\_CHAR

- COLLECTION:NXcollection →  
\_NXcollection.NX\_tree\_path NEXUSTREEPATH  
\_NXcollection.NX\_id COLLECTION  
\_NXcollection.NX\_scan\_id SCANID  
\_NXcollection.NX\_diffn\_id DIFFRNID  
\_NXcollection.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘‘/NXcollection\_\_COLLECTION’’ where COLLECTION is the name of this group, typically ‘‘collection’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- beamline:NX\_CHAR=BEAMLINE →  
\_NXcollection.beamline BEAMLINE

## 6.9 NXcollimator

Template of a beamline collimator.

NXcollimator (base class, version 1.0)  
absorbing\_material:NX\_CHAR

```

blade_spacing:NX_FLOAT
blade_thickness:NX_FLOAT
divergence_x:NX_FLOAT
divergence_y:NX_FLOAT
frequency:NX_FLOAT
soller_angle:NX_FLOAT
transmitting_material:NX_CHAR
type:NX_CHAR
NXgeometry
frequency_log:NXlog

```

- COLLIMATOR:NXcollimator →  
    \_NXcollimator.NX\_tree\_path NEXUSTREEPATH  
    \_NXcollimator.NX\_id COLLIMATOR  
    \_NXcollimator.NX\_scan\_id SCANID  
    \_NXcollimator.NX\_diffn\_id DIFFRNID  
    \_NXcollimator.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXcollimator\_COLLIMATOR’ where COLLIMATOR is the name of this group, typically ‘collimator’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- absorbing\_material:NX\_CHAR=ABSORBING\_MATERIAL →  
    \_NXcollimator.absorbing\_material ABSORBING\_MATERIAL
- blade\_spacing:NX\_FLOAT=BLADE\_SPACING →  
    \_NXcollimator.blade\_spacing BLADE\_SPACING
- blade\_thickness:NX\_FLOAT=BLADE\_THICKNESS →  
    \_NXcollimator.blade\_thickness BLADE\_THICKNESS
- divergence\_x:NX\_FLOAT=DIVERGENCE\_X →  
    \_NXcollimator.divergence\_x DIVERGENCE\_X
- divergence\_y:NX\_FLOAT=DIVERGENCE\_Y →  
    \_NXcollimator.divergence\_y DIVERGENCE\_Y
- frequency:NX\_FLOAT=FREQUENCY →  
    \_NXcollimator.frequency FREQUENCY
- soller\_angle:NX\_FLOAT=SOLLER\_ANGLE →  
    \_NXcollimator.soller\_angle SOLLER\_ANGLE
- transmitting\_material:NX\_CHAR=TRANSMITTING\_MATERIAL →  
    \_NXcollimator.transmitting\_material TRANSMITTING\_MATERIAL
- type:NX\_CHAR=TYPE →  
    \_NXcollimator.type TYPE
- geometry1:NXgeometry →  
    \_NXcollimator.NXgeometry\_id geometry1
- frequency\_log:NXlog →  
    \_NXcollimator.NXlog\_id frequency\_log

See also \_diffn\_radiation.div\_x\_source and \_diffn\_radiation.div\_y\_source

## 6.10 NXcrystal

Template of a crystal monochromator or analyzer. Permits double bent monochromator comprised of multiple segments with anisotropic Gaussian mosaic.

If curvatures are set to zero or are absent, array is considered to be flat.

Scattering vector is perpendicular to surface. Crystal is oriented parallel to beam incident on crystal before rotation, and lies in vertical plane.

```
NXcrystal (base class, version 1.0)
 azimuthal_angle:NX_FLOAT[i]
 bragg_angle:NX_FLOAT[i]
 chemical_formula:NX_CHAR
 curvature_horizontal:NX_FLOAT
 curvature_vertical:NX_FLOAT
 cut_angle:NX_FLOAT
 cylindrical_orientation_angle:NX_NUMBER
 d_spacing:NX_FLOAT
 density:NX_NUMBER
 is_cylindrical:NX_BOOLEAN
 mosaic_horizontal:NX_FLOAT
 mosaic_vertical:NX_FLOAT
 order_no:NX_INT
 orientation_matrix:NX_FLOAT[3,3]
 polar_angle:NX_FLOAT[i]
 reflection:NX_INT[3]
 scattering_vector:NX_FLOAT
 segment_columns:NX_FLOAT
 segment_gap:NX_FLOAT
 segment_height:NX_FLOAT
 segment_rows:NX_FLOAT
 segment_thickness:NX_FLOAT
 segment_width:NX_FLOAT
 space_group:NX_CHAR
 temperature:NX_FLOAT
 temperature_coefficient:NX_FLOAT
 thickness:NX_FLOAT
 type:NX_CHAR
 unit_cell:NX_FLOAT[n_comp,6]
 unit_cell_a:NX_FLOAT
 unit_cell_alpha:NX_FLOAT
 unit_cell_b:NX_FLOAT
 unit_cell_beta:NX_FLOAT
 unit_cell_c:NX_FLOAT
 unit_cell_gamma:NX_FLOAT
 unit_cell_volume:NX_FLOAT
 usage:NX_CHAR
 wavelength:NX_FLOAT[i]
 reflectivity:NXdata
 transmission:NXdata
 NXgeometry
 temperature_log:NXlog
```

shape:NXshape

- CRYSTAL:NXcrystal →  
\_NXcrystal.NX\_tree\_path NEXUSTREEPATH  
\_NXcrystal.NX\_id CRYSTAL  
\_NXcrystal.NX\_scan\_id SCANID  
\_NXcrystal.NX\_diffn\_id DIFFRNID  
\_NXcrystal.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXcrystal\_\_CRYSTAL’ where CRYSTAL is the name of this group, typically ‘crystal’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- azimuthal\_angle:NX\_FLOAT[i]=AZIMUTHAL\_ANGLE →  
\_NXcrystal.azimuthal\_angle AZIMUTHAL\_ANGLE
- bragg\_angle:NX\_FLOAT[i]=BRAGG\_ANGLE →  
\_NXcrystal.bragg\_angle BRAGG\_ANGLE
- chemical\_formula:NX\_CHAR=CHEMICAL\_FORMULA →  
\_NXcrystal.chemical\_formula CHEMICAL\_FORMULA
- curvature\_horizontal:NX\_FLOAT=CURVATURE\_HORIZONTAL →  
\_NXcrystal.curvature\_horizontal CURVATURE\_HORIZONTAL
- curvature\_vertical:NX\_FLOAT=CURVATURE\_VERTICAL →  
\_NXcrystal.curvature\_vertical CURVATURE\_VERTICAL
- cut\_angle:NX\_FLOAT=CUT\_ANGLE →  
\_NXcrystal.cut\_angle CUT\_ANGLE
- cylindrical\_orientation\_angle:NX\_NUMBER=CYLINDRICAL\_ORIENTATION\_ANGLE →  
\_NXcrystal.cylindrical\_orientation\_angle CYLINDRICAL\_ORIENTATION\_ANGLE
- d\_spacing:NX\_FLOAT=D\_SPACING →  
\_NXcrystal.d\_spacing D\_SPACING
- density:NX\_NUMBER=DENSITY →  
\_NXcrystal.density DENSITY
- is\_cylindrical:NX\_BOOLEAN=IS\_CYLINDRICAL →  
\_NXcrystal.is\_cylindrical IS\_CYLINDRICAL
- mosaic\_horizontal:NX\_FLOAT=MOSAIC\_HORIZONTAL →  
\_NXcrystal.mosaic\_horizontal MOSAIC\_HORIZONTAL
- mosaic\_vertical:NX\_FLOAT=MOSAIC\_VERTICAL →  
\_NXcrystal.mosaic\_vertical MOSAIC\_VERTICAL
- order\_no:NX\_INT=ORDER\_NO →  
\_NXcrystal.order\_no ORDER\_NO
- orientation\_matrix:NX\_FLOAT[3,3]=ORIENTATION\_MATRIX →  
\_NXcrystal.orientation\_matrix ORIENTATION\_MATRIX
- polar\_angle:NX\_FLOAT[i]=POLAR\_ANGLE →  
\_NXcrystal.polar\_angle POLAR\_ANGLE
- reflection:NX\_INT[3]=REFLECTION →  
\_NXcrystal.reflection REFLECTION

- scattering\_vector:NX\_FLOAT=SCATTERING\_VECTOR →  
\_NXcrystal.scattering\_vector SCATTERING\_VECTOR
- segment\_columns:NX\_FLOAT=SEGMENT\_COLUMNS →  
\_NXcrystal.segment\_columns SEGMENT\_COLUMNS
- segment\_gap:NX\_FLOAT=SEGMENT\_GAP →  
\_NXcrystal.segment\_gap SEGMENT\_GAP
- segment\_height:NX\_FLOAT=SEGMENT\_HEIGHT →  
\_NXcrystal.segment\_height SEGMENT\_HEIGHT
- segment\_rows:NX\_FLOAT=SEGMENT\_ROWS →  
\_NXcrystal.segment\_rows SEGMENT\_ROWS
- segment\_thickness:NX\_FLOAT=SEGMENT\_THICKNESS →  
\_NXcrystal.segment\_thickness SEGMENT\_THICKNESS
- segment\_width:NX\_FLOAT=SEGMENT\_WIDTH →  
\_NXcrystal.segment\_width SEGMENT\_WIDTH
- space\_group:NX\_CHAR=SPACE\_GROUP →  
\_NXcrystal.space\_group SPACE\_GROUP
- temperature:NX\_FLOAT=TEMPERATURE →  
\_NXcrystal.temperature TEMPERATURE
- temperature\_coefficient:NX\_FLOAT=TEMPERATURE\_COEFFICIENT →  
\_NXcrystal.temperature\_coefficient TEMPERATURE\_COEFFICIENT
- thickness:NX\_FLOAT=THICKNESS →  
\_NXcrystal.thickness THICKNESS
- type:NX\_CHAR=TYPE →  
\_NXcrystal.type TYPE
- unit\_cell:NX\_FLOAT[n\_comp,6]=UNIT\_CELL →  
\_NXcrystal.unit\_cell UNIT\_CELL
- unit\_cell\_a:NX\_FLOAT=UNIT\_CELL\_A →  
\_NXcrystal.unit\_cell.a UNIT\_CELL\_A
- unit\_cell\_alpha:NX\_FLOAT=UNIT\_CELL\_ALPHA →  
\_NXcrystal.unit\_cell.alpha UNIT\_CELL\_ALPHA
- unit\_cell\_b:NX\_FLOAT=UNIT\_CELL\_B →  
\_NXcrystal.unit\_cell.b UNIT\_CELL\_B
- unit\_cell\_beta:NX\_FLOAT=UNIT\_CELL\_BETA →  
\_NXcrystal.unit\_cell.beta UNIT\_CELL\_BETA
- unit\_cell\_c:NX\_FLOAT=UNIT\_CELL\_C →  
\_NXcrystal.unit\_cell.c UNIT\_CELL\_C
- unit\_cell\_gamma:NX\_FLOAT=UNIT\_CELL\_GAMMA →  
\_NXcrystal.unit\_cell.gamma UNIT\_CELL\_GAMMA
- unit\_cell\_volume:NX\_FLOAT=UNIT\_CELL\_VOLUME →  
\_NXcrystal.unit\_cell.volume UNIT\_CELL\_VOLUME
- usage:NX\_CHAR=USAGE →  
\_NXcrystal.usage USAGE
- wavelength:NX\_FLOAT[i]=WAVELENGTH →  
\_NXcrystal.wavelength WAVELENGTH

- reflectivity:NXdata →  
\_NXcrystal.NXdata\_id reflectivity
- transmission:NXdata →  
\_NXcrystal.NXdata\_id transmission
- geometry1:NXgeometry →  
\_NXcrystal.NXgeometry\_id geometry1
- temperature\_log:NXlog →  
\_NXcrystal.NXlog\_id temperature\_log
- shape:NXshape →  
\_NXcrystal.NXshape\_id shape

## 6.11 NXdata

(required) NXdata is a template of plottable data and their dimension scales. It is mandatory that there is at least one NXdata group in each NXentry group. Note that the variable and data can be defined with different names. The signal and axes attribute of the data item define which items are plottable data and which are dimension scales.

Each NXdata group will consist of only one data set containing plottable data and their standard deviations. This data set may be of arbitrary rank up to a maximum of NX\_MAXRANK=32. The plottable data will be identified by the attribute: signal=1 The plottable data will identify the dimension scale specification(s) in the axes attribute.

If available, the standard deviations of the data are to be stored in a data set of the same rank and dimensions, with the name errors.

For each data dimension, there should be a one-dimensional array of the same length. These one-dimensional arrays are the dimension scales of the data, i.e. the values of the independent variables at which the data is measured, such as scattering angle or energy transfer.

There are two methods of linking each data dimension to its respective dimension scale.

The preferred (and recommended) method uses the axes attribute to specify the names of each dimension scale.

The older method uses the axis attribute on each dimension scale to identify with an integer the axis whose value is the number of the dimension.

NXdata is used to implement one of the basic motivations in NeXus, to provide a default plot for the data of this NXentry. The actual data might be stored in another group and (hard) linked to the NXdata group.

NXdata (base class, version 1.0)

```
data:NX_NUMBER[n]
 @signal
 @axes
 @uncertainties
 @long_name
errors:NX_NUMBER[n]
offset:NX_FLOAT
scaling_factor:NX_FLOAT
variable:NX_NUMBER[n]
 @long_name
 @distribution
 @first_good
 @last_good
```

```

@axis
variable_errors:NX_NUMBER[n]
x:NX_FLOAT[nx]
y:NX_FLOAT[ny]
z:NX_FLOAT[nz]

```

- DATA:NXdata →  
    \_NXdata.NX\_tree\_path NEXUSTREEPATH  
    \_NXdata.NX\_id DATA  
    \_NXdata.NX\_scan\_id SCANID  
    \_NXdata.NX\_diffraction\_id DIFFRACTIONID  
    \_NXdata.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXdata\_DATA’ where DATA is the name of this group, typically ‘data’. The SCANID, DIFFRACTIONID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

```
_NXdata.NX_id DATA
```

is an alias for

```
_array_intensities.array_id DATA
```

and requires that we have a row in ARRAY\_DATA for which we have

```
_array_data.array_id DATAID
```

- data:NX\_NUMBER[n]=DATA →  
    \_NXdata.data DATA  
    is an alias of \_array\_data.data DATA for a row of ARRAY\_DATA for which we have \_array\_data.array\_id DATAID
- @signal=SIGNAL →  
    \_NXdata.data\_\_signal SIGNAL
- @axes=AXES →  
    \_NXdata.data\_\_axes AXES
- @uncertainties=UNCERTAINTIES →  
    \_NXdata.data\_\_uncertainties UNCERTAINTIES
- @long\_name=LONG\_NAME →  
    \_NXdata.data\_\_long\_name LONG\_NAME
- errors:NX\_NUMBER[n]=ERRORS →  
    \_NXdata.errors ERRORS
- offset:NX\_FLOAT=OFFSET →  
    \_NXdata.offset OFFSET is an alias for \_array\_intensities.offset OFFSET
- scaling\_factor:NX\_FLOAT=SCALING\_FACTOR →  
    \_NXdata.scaling\_factor SCALING\_FACTOR  
    is an alias for \_array\_intensities.scaling SCALING\_FACTOR
- variable:NX\_NUMBER[n]=VARIABLE →  
    \_NXdata.variable VARIABLE
- @long\_name=LONG\_NAME →  
    \_NXdata.variable\_\_long\_name LONG\_NAME
- @distribution=DISTRIBUTION →  
    \_NXdata.variable\_\_distribution DISTRIBUTION

- @first\_good=FIRST\_GOOD →  
\_NXdata.variable\_\_first\_good FIRST\_GOOD
- @last\_good=LAST\_GOOD →  
\_NXdata.variable\_\_last\_good LAST\_GOOD
- @axis=AXIS →  
\_NXdata.variable\_\_axis AXIS
- variable\_errors:NX\_NUMBER[n]=VARIABLE\_ERRORS →  
\_NXdata.variable\_errors VARIABLE\_ERRORS
- x:NX\_FLOAT[nx]=X →  
\_NXdata.x X
- y:NX\_FLOAT[ny]=Y →  
\_NXdata.y Y
- z:NX\_FLOAT[nz]=Z →  
\_NXdata.z Z

These items have a relationship with the CBF ARRAY\_INTENSITIES category, as noted above as aliases, Additional mappings into the various ARRAY categories should be established.

## 6.12 NXdetector

Template of a detector, detector bank, or multidetector.

```
NXdetector (base class, version 1.0)
 acquisition_mode:NX_CHAR
 angular_calibration:NX_FLOAT[i,j]
 angular_calibration_applied:NX_BOOLEAN
 azimuthal_angle:NX_FLOAT[np,i,j]
 beam_center_x:NX_FLOAT
 beam_center_y:NX_FLOAT
 bit_depth_readout:NX_INT
 calibration_date:NX_DATE_TIME
 count_time:NX_NUMBER[np]
 countrate_correction__applied:NX_BOOLEAN
 crate:NX_INT[i,j]
 @local_name
 data:NX_NUMBER[np,i,j,tof]
 @signal
 @axes
 @long_name
 @check_sum
 @link
 data_error:NX_NUMBER[np,i,j,tof]
 @units
 @link
 dead_time:NX_FLOAT[np,i,j]
 description:NX_CHAR
 detection_gas_path:NX_FLOAT
 detector_number:NX_INT[i,j]
 detector_readout_time:NX_FLOAT
```



```

diameter:NX_FLOAT
distance:NX_FLOAT[np,i,j]
flatfield:NX_FLOAT[i,j]
flatfield_applied:NX_BOOLEAN
flatfield_error:NX_FLOAT[i,j]
frame_start_number:NX_INT
frame_time:NX_FLOAT[NP]
gain_setting:NX_CHAR
gas_pressure:NX_FLOAT[i,j]
input:NX_INT[i,j]
 @local_name
layout:NX_CHAR
local_name:NX_CHAR
number_of_cycles:NX_INT
pixel_mask:NX_FLOAT[i,j]
pixel_mask_applied:NX_BOOLEAN
polar_angle:NX_FLOAT[np,i,j]
raw_time_of_flight:NX_INT[tof+1]
 @frequency
saturation_value:NX_INT
sensor_material:NX_CHAR
sensor_thickness:NX_FLOAT
sequence_number:NX_CHAR
slot:NX_INT[i,j]
 @local_name
solid_angle:NX_FLOAT[i,j]
threshold_energy:NX_FLOAT
time_of_flight:NX_FLOAT[tof+1]
 @axis
 @primary
 @long_name
 @link
trigger_dead_time:NX_FLOAT
trigger_delay_time:NX_FLOAT
type:NX_CHAR
x_pixel_offset:NX_FLOAT[i,j]
 @axis
 @primary
 @long_name
 @link
x_pixel_size:NX_FLOAT[i,j]
y_pixel_offset:NX_FLOAT[i,j]
 @axis
 @primary
 @long_name
y_pixel_size:NX_FLOAT[i,j]
NXcharacterization
efficiency:NXdata
 efficiency:NX_FLOAT[i,j,k]
 real_time:NX_NUMBER[i,j,k]

```

```

wavelength:NX_FLOAT[i,j,k]
NXgeometry
calibration_method:NXnote
data_file:NXnote

```

- DETECTOR:NXdetector →

```

_NXdetector.NX_tree_path NEXUSTREEPATH
_NXdetector.NX_id DETECTOR
_NXdetector.NX_scan_id SCANID
_NXdetector.NX_diffn_id DIFFRNID
_NXdetector.NX_entry_id ENTRYID

```

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXdetector\_\_DETECTOR’ where DETECTOR is the name of this group, typically ‘detector’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- acquisition\_mode:NX\_CHAR=ACQUISITION\_MODE →  
\_NXdetector.acquisition\_mode ACQUISITION\_MODE
- angular\_calibration:NX\_FLOAT[i,j]=ANGULAR\_CALIBRATION →  
\_NXdetector.angular\_calibration ANGULAR\_CALIBRATION
- angular\_calibration\_applied:NX\_BOOLEAN=ANGULAR\_CALIBRATION\_APPLIED →  
\_NXdetector.angular\_calibration\_applied ANGULAR\_CALIBRATION\_APPLIED
- azimuthal\_angle:NX\_FLOAT[np,i,j]=AZIMUTHAL\_ANGLE →  
\_NXdetector.azimuthal\_angle AZIMUTHAL\_ANGLE
- beam\_center\_x:NX\_FLOAT=BEAM\_CENTER\_X →  
\_NXdetector.beam\_center\_x BEAM\_CENTER\_X
- beam\_center\_y:NX\_FLOAT=BEAM\_CENTER\_Y →  
\_NXdetector.beam\_center\_y BEAM\_CENTER\_Y
- bit\_depth\_readout:NX\_INT=BIT\_DEPTH\_READOUT →  
\_NXdetector.bit\_depth\_readout BIT\_DEPTH\_READOUT
- calibration\_date:NX\_DATE.TIME=CALIBRATION\_DATE →  
\_NXdetector.calibration\_date CALIBRATION\_DATE
- count\_time:NX\_NUMBER[np]=COUNT\_TIME →  
\_NXdetector.count\_time COUNT\_TIME
- countrate\_correction\_applied:NX\_BOOLEAN=COUNTRATE\_CORRECTION\_\_APPLIED →  
\_NXdetector.countrate\_correction\_applied COUNTRATE\_CORRECTION\_\_APPLIED
- crate:NX\_INT[i,j]=CRATE →  
\_NXdetector.crate CRATE
- @local\_name=LOCAL\_NAME →  
\_NXdetector.crate\_\_local\_name LOCAL\_NAME
- data:NX\_NUMBER[np,i,j,tof]=DATA →  
\_NXdetector.data DATA
- @signal=SIGNAL →  
\_NXdetector.data\_\_signal SIGNAL
- @axes=AXES →  
\_NXdetector.data\_\_axes AXES

- @long\_name=LONG\_NAME →  
\_NXdetector.data\_\_long\_name LONG\_NAME
- @check\_sum=CHECK\_SUM →  
\_NXdetector.data\_\_check\_sum CHECK\_SUM
- @link=LINK →  
\_NXdetector.data\_\_link LINK
- data\_error:NX\_NUMBER[np,i,j,tof]=DATA\_ERROR →  
\_NXdetector.data\_error DATA\_ERROR
- @units=UNITS →  
\_NXdetector.data\_error\_\_units UNITS
- @link=LINK →  
\_NXdetector.data\_error\_\_link LINK
- dead\_time:NX\_FLOAT[np,i,j]=DEAD\_TIME →  
\_NXdetector.dead\_time DEAD\_TIME  
is an alias for  
\_diffrn\_detector.dtime DEAD\_TIME
- description:NX\_CHAR=DESCRIPTION →  
\_NXdetector.description DESCRIPTION
- detection\_gas\_path:NX\_FLOAT=DETECTION\_GAS\_PATH →  
\_NXdetector.detection\_gas\_path DETECTION\_GAS\_PATH
- detector\_number:NX\_INT[i,j]=DETECTOR\_NUMBER →  
\_NXdetector.detector\_number DETECTOR\_NUMBER
- detector\_readout\_time:NX\_FLOAT=DETECTOR\_READOUT\_TIME →  
\_NXdetector.detector\_readout\_time DETECTOR\_READOUT\_TIME
- diameter:NX\_FLOAT=DIAMETER →  
\_NXdetector.diameter DIAMETER
- distance:NX\_FLOAT[np,i,j]=DISTANCE →  
\_NXdetector.distance DISTANCE
- flatfield:NX\_FLOAT[i,j]=FLATFIELD →  
\_NXdetector.flatfield FLATFIELD
- flatfield\_applied:NX\_BOOLEAN=FLATFIELD\_APPLIED →  
\_NXdetector.flatfield\_applied FLATFIELD\_APPLIED
- flatfield\_error:NX\_FLOAT[i,j]=FLATFIELD\_ERROR →  
\_NXdetector.flatfield\_error FLATFIELD\_ERROR
- frame\_start\_number:NX\_INT=FRAME\_START\_NUMBER →  
\_NXdetector.frame\_start\_number FRAME\_START\_NUMBER
- frame\_time:NX\_FLOAT[NP]=FRAME\_TIME →  
\_NXdetector.frame\_time FRAME\_TIME
- gain\_setting:NX\_CHAR=GAIN\_SETTING →  
\_NXdetector.gain\_setting GAIN\_SETTING
- gas\_pressure:NX\_FLOAT[i,j]=GAS\_PRESSURE →  
\_NXdetector.gas\_pressure GAS\_PRESSURE
- input:NX\_INT[i,j]=INPUT →  
\_NXdetector.input INPUT

- @local\_name=LOCAL\_NAME →  
\_NXdetector.input\_\_local\_name LOCAL\_NAME
- layout:NX\_CHAR=LAYOUT →  
\_NXdetector.layout LAYOUT
- local\_name:NX\_CHAR=LOCAL\_NAME →  
\_NXdetector.local\_name LOCAL\_NAME
- number\_of\_cycles:NX\_INT=NUMBER\_OF\_CYCLES →  
\_NXdetector.number\_of\_cycles NUMBER\_OF\_CYCLES
- pixel\_mask:NX\_FLOAT[i,j]=PIXEL\_MASK →  
\_NXdetector.pixel\_mask PIXEL\_MASK
- pixel\_mask\_applied:NX\_BOOLEAN=PIXEL\_MASK\_APPLIED →  
\_NXdetector.pixel\_mask\_applied PIXEL\_MASK\_APPLIED
- polar\_angle:NX\_FLOAT[np,i,j]=POLAR\_ANGLE →  
\_NXdetector.polar\_angle POLAR\_ANGLE
- raw\_time\_of\_flight:NX\_INT[tof+1]=RAW\_TIME\_OF\_FLIGHT →  
\_NXdetector.raw\_time\_of\_flight RAW\_TIME\_OF\_FLIGHT
- @frequency=FREQUENCY →  
\_NXdetector.raw\_time\_of\_flight\_\_frequency FREQUENCY
- saturation\_value:NX\_INT=SATURATION\_VALUE →  
\_NXdetector.saturation\_value SATURATION\_VALUE
- sensor\_material:NX\_CHAR=SENSOR\_MATERIAL →  
\_NXdetector.sensor\_material SENSOR\_MATERIAL
- sensor\_thickness:NX\_FLOAT=SENSOR\_THICKNESS →  
\_NXdetector.sensor\_thickness SENSOR\_THICKNESS
- sequence\_number:NX\_CHAR=SEQUENCE\_NUMBER →  
\_NXdetector.sequence\_number SEQUENCE\_NUMBER
- slot:NX\_INT[i,j]=SLOT →  
\_NXdetector.slot SLOT
- @local\_name=LOCAL\_NAME →  
\_NXdetector.slot\_\_local\_name LOCAL\_NAME
- solid\_angle:NX\_FLOAT[i,j]=SOLID\_ANGLE →  
\_NXdetector.solid\_angle SOLID\_ANGLE
- threshold\_energy:NX\_FLOAT=THRESHOLD\_ENERGY →  
\_NXdetector.threshold\_energy THRESHOLD\_ENERGY
- time\_of\_flight:NX\_FLOAT[tof+1]=TIME\_OF\_FLIGHT →  
\_NXdetector.time\_of\_flight TIME\_OF\_FLIGHT
- @axis=AXIS →  
\_NXdetector.time\_of\_flight\_\_axis AXIS
- @primary=PRIMARY →  
\_NXdetector.time\_of\_flight\_\_primary PRIMARY
- @long\_name=LONG\_NAME →  
\_NXdetector.time\_of\_flight\_\_long\_name LONG\_NAME
- @link=LINK →  
\_NXdetector.time\_of\_flight\_\_link LINK

- trigger\_dead\_time:NX\_FLOAT=TRIGGER\_DEAD\_TIME →  
\_NXdetector.trigger\_dead\_time TRIGGER\_DEAD\_TIME
- trigger\_delay\_time:NX\_FLOAT=TRIGGER\_DELAY\_TIME →  
\_NXdetector.trigger\_delay\_time TRIGGER\_DELAY\_TIME
- type:NX\_CHAR=TYPE →  
\_NXdetector.type TYPE
- x\_pixel\_offset:NX\_FLOAT[i,j]=X\_PIXEL\_OFFSET →  
\_NXdetector.x\_pixel\_offset X\_PIXEL\_OFFSET
- @axis=AXIS →  
\_NXdetector.x\_pixel\_offset\_\_axis AXIS
- @primary=PRIMARY →  
\_NXdetector.x\_pixel\_offset\_\_primary PRIMARY
- @long\_name=LONG\_NAME →  
\_NXdetector.x\_pixel\_offset\_\_long\_name LONG\_NAME
- @link=LINK →  
\_NXdetector.x\_pixel\_offset\_\_link LINK
- x\_pixel\_size:NX\_FLOAT[i,j]=X\_PIXEL\_SIZE →  
\_NXdetector.x\_pixel\_size X\_PIXEL\_SIZE
- y\_pixel\_offset:NX\_FLOAT[i,j]=Y\_PIXEL\_OFFSET →  
\_NXdetector.y\_pixel\_offset Y\_PIXEL\_OFFSET
- @axis=AXIS →  
\_NXdetector.y\_pixel\_offset\_\_axis AXIS
- @primary=PRIMARY →  
\_NXdetector.y\_pixel\_offset\_\_primary PRIMARY
- @long\_name=LONG\_NAME →  
\_NXdetector.y\_pixel\_offset\_\_long\_name LONG\_NAME
- y\_pixel\_size:NX\_FLOAT[i,j]=Y\_PIXEL\_SIZE →  
\_NXdetector.y\_pixel\_size Y\_PIXEL\_SIZE
- characterization1:NXcharacterization →  
\_NXdetector.NXcharacterization\_id characterization1
- efficiency:NXdata →  
\_NXdetector.NXdata\_id efficiency
- efficiency:NX\_FLOAT[i,j,k]=EFFICIENCY →  
\_NXdetector.efficiency EFFICIENCY
- real\_time:NX\_NUMBER[i,j,k]=REAL\_TIME →  
\_NXdetector.real\_time REAL\_TIME
- wavelength:NX\_FLOAT[i,j,k]=WAVELENGTH →  
\_NXdetector.wavelength WAVELENGTH  
is an alias for  
\_diffrn\_radiation\_wavelength.wavelength WAVELENGTH
- geometry1:NXgeometry →  
\_NXdetector.NXgeometry\_id geometry1
- calibration\_method:NXnote →  
\_NXdetector.NXnote\_id calibration\_method
- data\_file:NXnote →  
\_NXdetector.NXnote\_id data\_file

## 6.13 NXdetector\_group

NXdetector\_group (base class, version 1.0)

```
group_index:NX_INT[i]
group_names:NX_CHAR
group_parent:NX_INT[]
group_type:NX_INT[]
```

- DETECTOR\_GROUP:NXdetector\_group →  
\_NXdetector\_group.NX\_tree\_path NEXUSTREEPATH  
\_NXdetector\_group.NX\_id DETECTOR\_GROUP  
\_NXdetector\_group.NX\_scan\_id SCANID  
\_NXdetector\_group.NX\_diffn\_id DIFFRNID  
\_NXdetector\_group.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘‘/NXdetector\_group\_DETECTOR\_GROUP’’ where DETECTOR\_GROUP is the name of this group, typically ‘‘detector\_group’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- group\_index:NX\_INT[i]=GROUP\_INDEX →  
\_NXdetector\_group.group\_index GROUP\_INDEX
- group\_names:NX\_CHAR=GROUP\_NAMES →  
\_NXdetector\_group.group\_names GROUP\_NAMES
- group\_parent:NX\_INT[]=GROUP\_PARENT →  
\_NXdetector\_group.group\_parent GROUP\_PARENT
- group\_type:NX\_INT[]=GROUP\_TYPE →  
\_NXdetector\_group.group\_type GROUP\_TYPE

This is closely related to the DIFFRN\_DETECTOR\_ELEMENT category and a method of integration needs to be found.

## 6.14 NXdisk\_chopper

NXdisk\_chopper (base class, version 1.0)

```
distance:NX_FLOAT
pair_separation:NX_FLOAT
phase:NX_FLOAT
radius:NX_FLOAT
ratio:NX_INT
rotation_speed:NX_FLOAT
slit_angle:NX_FLOAT
slit_height:NX_FLOAT
slits:NX_INT
type:NX_CHAR
wavelength_range:NX_FLOAT[2]
NXgeometry
```

- DISK\_CHOPPER:NXdisk\_chopper →  
\_NXdisk\_chopper.NX\_tree\_path NEXUSTREEPATH  
\_NXdisk\_chopper.NX\_id DISK\_CHOPPER

```

_NXdisk.chopper.NX_scan_id SCANID
_NXdisk.chopper.NX_diffn_id DIFFRNID
_NXdisk.chopper.NX_entry_id ENTRYID

```

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXdisk.chopper\_DISK\_CHOPPER’ where DISK\_CHOPPER is the name of this group, typically ‘disk\_chopper’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- distance:NX\_FLOAT=DISTANCE →  
\_NXdisk.chopper.distance DISTANCE
- pair\_separation:NX\_FLOAT=PAIR\_SEPARATION →  
\_NXdisk.chopper.pair\_separation PAIR\_SEPARATION
- phase:NX\_FLOAT=PHASE →  
\_NXdisk.chopper.phase PHASE
- radius:NX\_FLOAT=RADIUS →  
\_NXdisk.chopper.radius RADIUS
- ratio:NX\_INT=RATIO →  
\_NXdisk.chopper.ratio RATIO
- rotation\_speed:NX\_FLOAT=ROTATION\_SPEED →  
\_NXdisk.chopper.rotation\_speed ROTATION\_SPEED
- slit\_angle:NX\_FLOAT=SLIT\_ANGLE →  
\_NXdisk.chopper.slit\_angle SLIT\_ANGLE
- slit\_height:NX\_FLOAT=SLIT\_HEIGHT →  
\_NXdisk.chopper.slit\_height SLIT\_HEIGHT
- slits:NX\_INT=SLITS →  
\_NXdisk.chopper.slits SLITS
- type:NX\_CHAR=TYPE →  
\_NXdisk.chopper.type TYPE
- wavelength\_range:NX\_FLOAT[2]=WAVELENGTH\_RANGE →  
\_NXdisk.chopper.wavelength\_range WAVELENGTH\_RANGE
- geometry1:NXgeometry →  
\_NXdisk.chopper.NXgeometry\_id geometry1

## 6.15 NXentry

NXentry (base class, version 1.0)

```

@IDF_Version
collection_description:NX_CHAR
collection_identifier:NX_CHAR
collection_time:NX_FLOAT
definition:NX_CHAR
 @version
 @URL
definition_local:NX_CHAR
 @version
 @URL

```

```

duration:NX_INT
end_time:NX_DATE_TIME
entry_identifier:NX_CHAR
experiment_description:NX_CHAR
experiment_identifier:NX_CHAR
pre_sample_flightpath:NX_FLOAT
program_name:NX_CHAR
 @version
 @configuration
revision:NX_CHAR
 @comment
run_cycle:NX_CHAR
start_time:NX_DATE_TIME
title:NX_CHAR
NXcharacterization
NXdata
NXinstrument
NXmonitor
experiment_documentation:NXnote
notes:NXnote
thumbnail:NXnote
 @mime_type
NXprocess
NXsample
NXsubentry
NXuser

```

- ENTRY:NXentry →
  - \_NXentry.NX\_tree\_path NEXUSTREEPATH
  - \_NXentry.NX\_id ENTRY
  - \_NXentry.NX\_scan\_id SCANID
  - \_NXentry.NX\_diffn\_id DIFFRNID
  - \_NXentry.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXentry\_ENTRY’ where ENTRY is the name of this group, typically ‘entry’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- @IDF\_Version=IDF\_VERSION →
  - \_NXentry.NX\_class\_attribute\_\_IDF\_Version IDF\_VERSION
- collection\_description:NX\_CHAR=COLLECTION\_DESCRIPTION →
  - \_NXentry.collection\_description COLLECTION\_DESCRIPTION
- collection\_identifier:NX\_CHAR=COLLECTION\_IDENTIFIER →
  - \_NXentry.collection\_identifier COLLECTION\_IDENTIFIER
- collection\_time:NX\_FLOAT=COLLECTION\_TIME →
  - \_NXentry.collection\_time COLLECTION\_TIME
- definition:NX\_CHAR=DEFINITION →
  - \_NXentry.definition DEFINITION
- @version=VERSION →
  - \_NXentry.definition\_\_version VERSION



- @URL=URL →  
\_NXentry.definition\_\_URL URL
- definition\_local:NX\_CHAR=DEFINITION\_LOCAL →  
\_NXentry.definition\_local DEFINITION\_LOCAL
- @version=VERSION →  
\_NXentry.definition\_local\_\_version VERSION
- @URL=URL →  
\_NXentry.definition\_local\_\_URL URL
- duration:NX\_INT=DURATION →  
\_NXentry.duration DURATION
- end\_time:NX\_DATE\_TIME=END\_TIME →  
\_NXentry.end\_time END\_TIME
- entry\_identifier:NX\_CHAR=ENTRY\_IDENTIFIER →  
\_NXentry.entry\_identifier ENTRY\_IDENTIFIER
- experiment\_description:NX\_CHAR=EXPERIMENT\_DESCRIPTION →  
\_NXentry.experiment\_description EXPERIMENT\_DESCRIPTION
- experiment\_identifier:NX\_CHAR=EXPERIMENT\_IDENTIFIER →  
\_NXentry.experiment\_identifier EXPERIMENT\_IDENTIFIER
- pre\_sample\_flightpath:NX\_FLOAT=PRE\_SAMPLE\_FLIGHTPATH →  
\_NXentry.pre\_sample\_flightpath PRE\_SAMPLE\_FLIGHTPATH
- program\_name:NX\_CHAR=PROGRAM\_NAME →  
\_NXentry.program\_name PROGRAM\_NAME
- @version=VERSION →  
\_NXentry.program\_name\_\_version VERSION
- @configuration=CONFIGURATION →  
\_NXentry.program\_name\_\_configuration CONFIGURATION
- revision:NX\_CHAR=REVISION →  
\_NXentry.revision REVISION
- @comment=COMMENT →  
\_NXentry.revision\_\_comment COMMENT
- run\_cycle:NX\_CHAR=RUN\_CYCLE →  
\_NXentry.run\_cycle RUN\_CYCLE
- start\_time:NX\_DATE\_TIME=START\_TIME →  
\_NXentry.start\_time START\_TIME
- title:NX\_CHAR=TITLE →  
\_NXentry.title TITLE
- characterization1:NXcharacterization →  
\_NXentry.NXcharacterization\_id characterization1
- data1:NXdata →  
\_NXentry.NXdata\_id data1
- instrument1:NXinstrument →  
\_NXentry.NXinstrument\_id instrument1
- monitor1:NXmonitor →  
\_NXentry.NXmonitor\_id monitor1

- `experiment_documentation:NXnote` →  
`_NXentry.NXnote_id experiment_documentation`
- `notes:NXnote` →  
`_NXentry.NXnote_id notes`
- `thumbnail:NXnote` →  
`_NXentry.NXnote_id thumbnail`
- `@mime_type=MIME_TYPE` →  
`_NXentry.title_mime_type MIME_TYPE`
- `process1:NXprocess` →  
`_NXentry.NXprocess_id process1`
- `sample1:NXsample` →  
`_NXentry.NXsample_id sample1`
- `subentry1:NXsubentry` →  
`_NXentry.NXsubentry_id subentry1`
- `user1:NXuser` →  
`_NXentry.NXuser_id user1`

## 6.16 NXenvironment

NXenvironment (base class, version 1.0)

```
description:NX_CHAR
name:NX_CHAR
program:NX_CHAR
short_name:NX_CHAR
type:NX_CHAR
position:NXgeometry
NXnote
NXsensor
```

- `ENVIRONMENT:NXenvironment` →  
`_NXenvironment.NX_tree_path NEXUSTREEPATH`  
`_NXenvironment.NX_id ENVIRONMENT`  
`_NXenvironment.NX_scan_id SCANID`  
`_NXenvironment.NX_diffn_id DIFFRNID`  
`_NXenvironment.NX_entry_id ENTRYID`

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘`/NXenvironment_ENVIRONMENT`’ where ENVIRONMENT is the name of this group, typically ‘`environment`’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- `description:NX_CHAR=DESCRIPTION` →  
`_NXenvironment.description DESCRIPTION`
- `name:NX_CHAR=NAME` →  
`_NXenvironment.name NAME`
- `program:NX_CHAR=PROGRAM` →  
`_NXenvironment.program PROGRAM`

- short\_name:NX\_CHAR=SHORT\_NAME →  
\_NXenvironment.short\_name SHORT\_NAME
- type:NX\_CHAR=TYPE →  
\_NXenvironment.type TYPE
- position:NXgeometry →  
\_NXenvironment.NXgeometry\_id position
- note1:NXnote →  
\_NXenvironment.NXnote\_id note1
- sensor1:NXsensor →  
\_NXenvironment.NXsensor\_id sensor1

## 6.17 NXevent\_data

NXevent\_data (base class, version 1.0)

```
events_per_pulse:NX_INT[j]
pixel_number:NX_INT[i]
pulse_height:NX_FLOAT[i,k]
pulse_time:NX_INT[j]
@offset
time_of_flight:NX_INT[i]
```

- EVENT\_DATA:NXevent\_data →  
\_NXevent\_data.NX\_tree\_path NEXUSTREEPATH  
\_NXevent\_data.NX\_id EVENT\_DATA  
\_NXevent\_data.NX\_scan\_id SCANID  
\_NXevent\_data.NX\_diffn\_id DIFFRNID  
\_NXevent\_data.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXevent\_data\_EVENT\_DATA’ where EVENT\_DATA is the name of this group, typically ‘event\_data’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- events\_per\_pulse:NX\_INT[j]=EVENTS\_PER\_PULSE →  
\_NXevent\_data.events\_per\_pulse EVENTS\_PER\_PULSE
- pixel\_number:NX\_INT[i]=PIXEL\_NUMBER →  
\_NXevent\_data.pixel\_number PIXEL\_NUMBER
- pulse\_height:NX\_FLOAT[i,k]=PULSE\_HEIGHT →  
\_NXevent\_data.pulse\_height PULSE\_HEIGHT
- pulse\_time:NX\_INT[j]=PULSE\_TIME →  
\_NXevent\_data.pulse\_time PULSE\_TIME
- @offset=OFFSET →  
\_NXevent\_data.pulse\_time\_\_offset OFFSET
- time\_of\_flight:NX\_INT[i]=TIME\_OF\_FLIGHT →  
\_NXevent\_data.time\_of\_flight TIME\_OF\_FLIGHT

## 6.18 NXfermi\_chopper

NXfermi\_chopper (base class, version 1.0)

```
absorbing_material:NX_CHAR
distance:NX_FLOAT
energy:NX_FLOAT
height:NX_FLOAT
number:NX_INT
r_slit:NX_FLOAT
radius:NX_FLOAT
rotation_speed:NX_FLOAT
slit:NX_FLOAT
transmitting_material:NX_CHAR
type:NX_CHAR
wavelength:NX_FLOAT
width:NX_FLOAT
NXgeometry
```

- FERMI\_CHOPPER:NXfermi\_chopper →  
\_NXfermi\_chopper.NX\_tree\_path NEXUSTREEPATH  
\_NXfermi\_chopper.NX\_id FERMI\_CHOPPER  
\_NXfermi\_chopper.NX\_scan\_id SCANID  
\_NXfermi\_chopper.NX\_diffn\_id DIFFRNID  
\_NXfermi\_chopper.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXfermi\_chopper\_FERMI\_CHOPPER’ where FERMI\_CHOPPER is the name of this group, typically ‘fermi\_chopper’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- absorbing\_material:NX\_CHAR=ABSORBING\_MATERIAL →  
\_NXfermi\_chopper.absorbing\_material ABSORBING\_MATERIAL
- distance:NX\_FLOAT=DISTANCE →  
\_NXfermi\_chopper.distance DISTANCE
- energy:NX\_FLOAT=ENERGY →  
\_NXfermi\_chopper.energy ENERGY
- height:NX\_FLOAT=HEIGHT →  
\_NXfermi\_chopper.height HEIGHT
- number:NX\_INT=NUMBER →  
\_NXfermi\_chopper.number NUMBER
- r\_slit:NX\_FLOAT=R\_SLIT →  
\_NXfermi\_chopper.r\_slit R\_SLIT
- radius:NX\_FLOAT=RADIUS →  
\_NXfermi\_chopper.radius RADIUS
- rotation\_speed:NX\_FLOAT=ROTATION\_SPEED →  
\_NXfermi\_chopper.rotation\_speed ROTATION\_SPEED
- slit:NX\_FLOAT=SLIT →  
\_NXfermi\_chopper.slit SLIT

- transmitting\_material:NX\_CHAR=TRANSMITTING\_MATERIAL →  
\_NXfermi\_chopper.transmitting\_material TRANSMITTING\_MATERIAL
- type:NX\_CHAR=TYPE →  
\_NXfermi\_chopper.type TYPE
- wavelength:NX\_FLOAT=WAVELENGTH →  
\_NXfermi\_chopper.wavelength WAVELENGTH
- width:NX\_FLOAT=WIDTH →  
\_NXfermi\_chopper.width WIDTH
- geometry1:NXgeometry →  
\_NXfermi\_chopper.NXgeometry\_id geometry1

## 6.19 NXfilter

NXfilter (base class, version 1.0)

```
chemical_formula:NX_CHAR
coating_material:NX_CHAR
coating_roughness:NX_FLOAT[nsurf]
density:NX_NUMBER
description:NX_CHAR
m_value:NX_FLOAT
orientation_matrix:NX_FLOAT[n_comp,3,3]
status:NX_CHAR
substrate_material:NX_CHAR
substrate_roughness:NX_FLOAT
substrate_thickness:NX_FLOAT
temperature:NX_FLOAT
thickness:NX_FLOAT
unit_cell_a:NX_FLOAT
unit_cell_alpha:NX_FLOAT
unit_cell_b:NX_FLOAT
unit_cell_beta:NX_FLOAT
unit_cell_c:NX_FLOAT
unit_cell_gamma:NX_FLOAT
unit_cell_volume:NX_FLOAT[n_comp]
transmission:NXdata
NXgeometry
temperature_log:NXlog
sensor_type:NXsensor
```

- FILTER:NXfilter →  
\_NXfilter.NX.tree\_path NEXUSTREEPATH  
\_NXfilter.NX\_id FILTER  
\_NXfilter.NX.scan\_id SCANID  
\_NXfilter.NX.diffraction\_id DIFFRACTIONID  
\_NXfilter.NX.entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXfilter\_FILTER’ where FILTER is the name of this group, typically ‘filter’. The SCANID, DIFFRACTIONID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- chemical\_formula:NX\_CHAR=CHEMICAL\_FORMULA →  
\_NXfilter.chemical\_formula CHEMICAL\_FORMULA
- coating\_material:NX\_CHAR=COATING\_MATERIAL →  
\_NXfilter.coating\_material COATING\_MATERIAL
- coating\_roughness:NX\_FLOAT[nsurf]=COATING\_ROUGHNESS →  
\_NXfilter.coating\_roughness COATING\_ROUGHNESS
- density:NX\_NUMBER=DENSITY →  
\_NXfilter.density DENSITY
- description:NX\_CHAR=DESCRIPTION →  
\_NXfilter.description DESCRIPTION
- m\_value:NX\_FLOAT=M\_VALUE →  
\_NXfilter.m\_value M\_VALUE
- orientation\_matrix:NX\_FLOAT[n\_comp,3,3]=ORIENTATION\_MATRIX →  
\_NXfilter.orientation\_matrix ORIENTATION\_MATRIX
- status:NX\_CHAR=STATUS →  
\_NXfilter.status STATUS
- substrate\_material:NX\_CHAR=SUBSTRATE\_MATERIAL →  
\_NXfilter.substrate\_material SUBSTRATE\_MATERIAL
- substrate\_roughness:NX\_FLOAT=SUBSTRATE\_ROUGHNESS →  
\_NXfilter.substrate\_roughness SUBSTRATE\_ROUGHNESS
- substrate\_thickness:NX\_FLOAT=SUBSTRATE\_THICKNESS →  
\_NXfilter.substrate\_thickness SUBSTRATE\_THICKNESS
- temperature:NX\_FLOAT=TEMPERATURE →  
\_NXfilter.temperature TEMPERATURE
- thickness:NX\_FLOAT=THICKNESS →  
\_NXfilter.thickness THICKNESS
- unit\_cell\_a:NX\_FLOAT=UNIT\_CELL\_A →  
\_NXfilter.unit\_cell\_a UNIT\_CELL\_A
- unit\_cell\_alpha:NX\_FLOAT=UNIT\_CELL\_ALPHA →  
\_NXfilter.unit\_cell\_alpha UNIT\_CELL\_ALPHA
- unit\_cell\_b:NX\_FLOAT=UNIT\_CELL\_B →  
\_NXfilter.unit\_cell\_b UNIT\_CELL\_B
- unit\_cell\_beta:NX\_FLOAT=UNIT\_CELL\_BETA →  
\_NXfilter.unit\_cell\_beta UNIT\_CELL\_BETA
- unit\_cell\_c:NX\_FLOAT=UNIT\_CELL\_C →  
\_NXfilter.unit\_cell\_c UNIT\_CELL\_C
- unit\_cell\_gamma:NX\_FLOAT=UNIT\_CELL\_GAMMA →  
\_NXfilter.unit\_cell\_gamma UNIT\_CELL\_GAMMA
- unit\_cell\_volume:NX\_FLOAT[n\_comp]=UNIT\_CELL\_VOLUME →  
\_NXfilter.unit\_cell\_volume UNIT\_CELL\_VOLUME
- transmission:NXdata →  
\_NXfilter.NXdata\_id transmission
- geometry1:NXgeometry →  
\_NXfilter.NXgeometry\_id geometry1

- temperature\_log:NXlog →  
\_NXfilter.NXlog\_id temperature\_log
- sensor\_type:NXsensor →  
\_NXfilter.NXsensor\_id sensor\_type

## 6.20 NXflipper

NXflipper (base class, version 1.0)

```
comp_current:NX_FLOAT
comp_turns:NX_FLOAT
flip_current:NX_FLOAT
flip_turns:NX_FLOAT
guide_current:NX_FLOAT
guide_turns:NX_FLOAT
thickness:NX_FLOAT
type:NX_CHAR
```

- FLIPPER:NXflipper →  
\_NXflipper.NX\_tree\_path NEXUSTREEPATH  
\_NXflipper.NX\_id FLIPPER  
\_NXflipper.NX\_scan\_id SCANID  
\_NXflipper.NX\_diffn\_id DIFFRNID  
\_NXflipper.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXflipper\_FLIPPER’ where FLIPPER is the name of this group, typically ‘flipper’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- comp\_current:NX\_FLOAT=COMP\_CURRENT →  
\_NXflipper.comp\_current COMP\_CURRENT
- comp\_turns:NX\_FLOAT=COMP\_TURNS →  
\_NXflipper.comp\_turns COMP\_TURNS
- flip\_current:NX\_FLOAT=FLIP\_CURRENT →  
\_NXflipper.flip\_current FLIP\_CURRENT
- flip\_turns:NX\_FLOAT=FLIP\_TURNS →  
\_NXflipper.flip\_turns FLIP\_TURNS
- guide\_current:NX\_FLOAT=GUIDE\_CURRENT →  
\_NXflipper.guide\_current GUIDE\_CURRENT
- guide\_turns:NX\_FLOAT=GUIDE\_TURNS →  
\_NXflipper.guide\_turns GUIDE\_TURNS
- thickness:NX\_FLOAT=THICKNESS →  
\_NXflipper.thickness THICKNESS
- type:NX\_CHAR=TYPE →  
\_NXflipper.type TYPE

## 6.21 NXgeometry

NXgeometry (base class, version 1.0)

```

component_index:NX_INT
description:NX_CHAR
NXorientation
NXshape
NXtranslation

```

- GEOMETRY:NXgeometry →  
     NXgeometry.NX\_tree\_path NEXUSTREEPATH  
     NXgeometry.NX\_id GEOMETRY  
     NXgeometry.NX\_scan\_id SCANID  
     NXgeometry.NX\_diffn\_id DIFFRNID  
     NXgeometry.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXgeometry\_GEOMETRY’ where GEOMETRY is the name of this group, typically ‘geometry’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- component\_index:NX\_INT=COMPONENT\_INDEX →  
     NXgeometry.component\_index COMPONENT\_INDEX
- description:NX\_CHAR=DESCRIPTION →  
     NXgeometry.description DESCRIPTION
- orientation1:NXorientation →  
     NXgeometry.NXorientation\_id orientation1
- shape1:NXshape →  
     NXgeometry.NXshape\_id shape1
- translation1:NXtranslation →  
     NXgeometry.NXtranslation\_id translation1

## 6.22 NXguide

```

NXguide (base class, version 1.0)
bend_angle_x:NX_FLOAT
bend_angle_y:NX_FLOAT
coating_material:NX_FLOAT[nsurf]
coating_roughness:NX_FLOAT[nsurf]
description:NX_CHAR
external_material:NX_CHAR
incident_angle:NX_FLOAT
interior_atmosphere:NX_CHAR
m_value:NX_FLOAT[nsurf]
number_sections:NX_INT
substrate_material:NX_FLOAT[nsurf]
substrate_roughness:NX_FLOAT[nsurf]
substrate_thickness:NX_FLOAT[nsurf]
reflectivity:NXdata
 data:NX_NUMBER[nsurf,nw1]
 @signal
 @axes
 surface:NX_NUMBER[nsurf]
 wavelength:NX_NUMBER[nw1]

```



## NXgeometry

- GUIDE:NXguide →  
\_NXguide.NX\_tree\_path NEXUSTREEPATH  
\_NXguide.NX\_id GUIDE  
\_NXguide.NX\_scan\_id SCANID  
\_NXguide.NX\_diffn\_id DIFFRNID  
\_NXguide.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXguide\_GUIDE’ where GUIDE is the name of this group, typically ‘guide’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- bend\_angle\_x:NX\_FLOAT=BEND\_ANGLE\_X →  
\_NXguide.bend\_angle\_x BEND\_ANGLE\_X
- bend\_angle\_y:NX\_FLOAT=BEND\_ANGLE\_Y →  
\_NXguide.bend\_angle\_y BEND\_ANGLE\_Y
- coating\_material:NX\_FLOAT[nsurf]=COATING\_MATERIAL →  
\_NXguide.coating\_material COATING\_MATERIAL
- coating\_roughness:NX\_FLOAT[nsurf]=COATING\_ROUGHNESS →  
\_NXguide.coating\_roughness COATING\_ROUGHNESS
- description:NX\_CHAR=DESCRIPTION →  
\_NXguide.description DESCRIPTION
- external\_material:NX\_CHAR=EXTERNAL\_MATERIAL →  
\_NXguide.external\_material EXTERNAL\_MATERIAL
- incident\_angle:NX\_FLOAT=INCIDENT\_ANGLE →  
\_NXguide.incident\_angle INCIDENT\_ANGLE
- interior\_atmosphere:NX\_CHAR=INTERIOR\_ATMOSPHERE →  
\_NXguide.interior\_atmosphere INTERIOR\_ATMOSPHERE
- m\_value:NX\_FLOAT[nsurf]=M\_VALUE →  
\_NXguide.m\_value M\_VALUE
- number\_sections:NX\_INT=NUMBER\_SECTIONS →  
\_NXguide.number\_sections NUMBER\_SECTIONS
- substrate\_material:NX\_FLOAT[nsurf]=SUBSTRATE\_MATERIAL →  
\_NXguide.substrate\_material SUBSTRATE\_MATERIAL
- substrate\_roughness:NX\_FLOAT[nsurf]=SUBSTRATE\_ROUGHNESS →  
\_NXguide.substrate\_roughness SUBSTRATE\_ROUGHNESS
- substrate\_thickness:NX\_FLOAT[nsurf]=SUBSTRATE\_THICKNESS →  
\_NXguide.substrate\_thickness SUBSTRATE\_THICKNESS
- reflectivity:NXdata →  
\_NXguide.NXdata\_id reflectivity
- data:NX\_NUMBER[nsurf,nwl]=DATA →  
\_NXguide.data DATA
- @signal=SIGNAL →  
\_NXguide.data\_\_signal SIGNAL

- @axes=AXES →  
\_NXguide.data.axes AXES
- surface:NX\_NUMBER[nsurf]=SURFACE →  
\_NXguide.surface SURFACE
- wavelength:NX\_NUMBER[nwl]=WAVELENGTH →  
\_NXguide.wavelength WAVELENGTH
- geometry1:NXgeometry →  
\_NXguide.NXgeometry\_id geometry1

## 6.23 NXinsertion\_device

NXinsertion\_device (base class, version 1.0)

```
bandwidth:NX_FLOAT
energy:NX_FLOAT
gap:NX_FLOAT
harmonic:NX_INT
k:NX_FLOAT
length:NX_FLOAT
magnetic_wavelength:NX_FLOAT
phase:NX_FLOAT
poles:NX_INT
power:NX_FLOAT
taper:NX_FLOAT
type:NX_CHAR
spectrum:NXdata
NXgeometry
```

- INSERTION\_DEVICE:NXinsertion\_device →  
\_NXinsertion\_device.NX\_tree\_path NEXUSTREEPATH  
\_NXinsertion\_device.NX\_id INSERTION\_DEVICE  
\_NXinsertion\_device.NX\_scan\_id SCANID  
\_NXinsertion\_device.NX\_diffn\_id DIFFRNID  
\_NXinsertion\_device.NX\_entry\_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXinsertion\_device\_\_INSERTION\_DEVICE’ where INSERTION\_DEVICE is the name of this group, typically ‘insertion\_device’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- bandwidth:NX\_FLOAT=BANDWIDTH →  
\_NXinsertion\_device.bandwidth BANDWIDTH
- energy:NX\_FLOAT=ENERGY →  
\_NXinsertion\_device.energy ENERGY
- gap:NX\_FLOAT=GAP →  
\_NXinsertion\_device.gap GAP
- harmonic:NX\_INT=HARMONIC →  
\_NXinsertion\_device.harmonic HARMONIC
- k:NX\_FLOAT=K →  
\_NXinsertion\_device.k K

- `length:NX_FLOAT=LENGTH →`  
`_NXinsertion_device.length LENGTH`
- `magnetic_wavelength:NX_FLOAT=MAGNETIC_WAVELENGTH →`  
`_NXinsertion_device.magnetic_wavelength MAGNETIC_WAVELENGTH`
- `phase:NX_FLOAT=PHASE →`  
`_NXinsertion_device.phase PHASE`
- `poles:NX_INT=POLES →`  
`_NXinsertion_device.poles POLES`
- `power:NX_FLOAT=POWER →`  
`_NXinsertion_device.power POWER`
- `taper:NX_FLOAT=TAPER →`  
`_NXinsertion_device.taper TAPER`
- `type:NX_CHAR=TYPE →`  
`_NXinsertion_device.type TYPE`
- `spectrum:NXdata →`  
`_NXinsertion_device.NXdata_id spectrum`
- `geometry1:NXgeometry →`  
`_NXinsertion_device.NXgeometry_id geometry1`

## 6.24 NXinstrument

NXinstrument (base class, version 1.0)

`name:NX_CHAR`

`@short_name`

`NXaperture`

`NXattenuator`

`NXbeam`

`NXbeam_stop`

`NXbending_magnet`

`NXcollimator`

`NXcrystal`

`NXdetector`

`NXdisk_chopper`

`NXfermi_chopper`

`NXfilter`

`NXflipper`

`NXguide`

`NXinsertion_device`

`NXmirror`

`NXmoderator`

`NXpolarizer`

`NXsource`

`NXvelocity_selector`

- `INSTRUMENT:NXinstrument →`  
`_NXinstrument.NX_tree_path NEXUSTREEPATH`  
`_NXinstrument.NX_id INSTRUMENT`  
`_NXinstrument.NX_scan_id SCANID`

- ```

_NXinstrument.NX.diffrn_id DIFFRNID
_NXinstrument.NX.entry_id ENTRYID

```
- where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘‘/NXinstrument__INSTRUMENT’’ where INSTRUMENT is the name of this group, typically ‘‘instrument’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.
- name:NX_CHAR=NAME →
_NXinstrument.name NAME
 - @short_name=SHORT_NAME →
_NXinstrument.name__short_name SHORT_NAME
 - aperture1:NXaperture →
_NXinstrument.NXaperture_id aperture1
 - attenuator1:NXattenuator →
_NXinstrument.NXattenuator_id attenuator1
 - beam1:NXbeam →
_NXinstrument.NXbeam_id beam1
 - beam_stop1:NXbeam_stop →
_NXinstrument.NXbeam_stop_id beam_stop1
 - bending_magnet1:NXbending_magnet →
_NXinstrument.NXbending_magnet_id bending_magnet1
 - collimator1:NXcollimator →
_NXinstrument.NXcollimator_id collimator1
 - crystal1:NXcrystal →
_NXinstrument.NXcrystal_id crystal1
 - detector1:NXdetector →
_NXinstrument.NXdetector_id detector1
 - disk_chopper1:NXdisk_chopper →
_NXinstrument.NXdisk_chopper_id disk_chopper1
 - fermi_chopper1:NXfermi_chopper →
_NXinstrument.NXfermi_chopper_id fermi_chopper1
 - filter1:NXfilter →
_NXinstrument.NXfilter_id filter1
 - flipper1:NXflipper →
_NXinstrument.NXflipper_id flipper1
 - guide1:NXguide →
_NXinstrument.NXguide_id guide1
 - insertion_device1:NXinsertion_device →
_NXinstrument.NXinsertion_device_id insertion_device1
 - mirror1:NXmirror →
_NXinstrument.NXmirror_id mirror1
 - moderator1:NXmoderator →
_NXinstrument.NXmoderator_id moderator1

- polarizer1:NXpolarizer →
_NXinstrument.NXpolarizer_id polarizer1
- source1:NXsource →
_NXinstrument.NXsource_id source1
- velocity_selector1:NXvelocity_selector →
_NXinstrument.NXvelocity_selector_id velocity_selector1

6.25 NXlog

NXlog (base class, version 1.0)

```
average_value:NX_FLOAT
average_value_error:NX_FLOAT
description:NX_CHAR
duration:NX_FLOAT
maximum_value:NX_FLOAT
minimum_value:NX_FLOAT
raw_value:NX_NUMBER
time:NX_FLOAT
  @start
value:NX_NUMBER
```

- LOG:NXlog →
_NXlog.NX_tree_path NEXUSTREEPATH
_NXlog.NX_id LOG
_NXlog.NX_scan_id SCANID
_NXlog.NX_diffn_id DIFFRNID
_NXlog.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXlog_LOG’ where LOG is the name of this group, typically ‘log’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- average_value:NX_FLOAT=AVERAGE_VALUE →
_NXlog.average_value AVERAGE_VALUE
- average_value_error:NX_FLOAT=AVERAGE_VALUE_ERROR →
_NXlog.average_value_error AVERAGE_VALUE_ERROR
- description:NX_CHAR=DESCRIPTION →
_NXlog.description DESCRIPTION
- duration:NX_FLOAT=DURATION →
_NXlog.duration DURATION
- maximum_value:NX_FLOAT=MAXIMUM_VALUE →
_NXlog.maximum_value MAXIMUM_VALUE
- minimum_value:NX_FLOAT=MINIMUM_VALUE →
_NXlog.minimum_value MINIMUM_VALUE
- raw_value:NX_NUMBER=RAW_VALUE →
_NXlog.raw_value RAW_VALUE
- time:NX_FLOAT=TIME →
_NXlog.time TIME

- @start=START →
_NXlog.time_start START
- value:NX_NUMBER=VALUE →
_NXlog.value VALUE

6.26 NXmirror

NXmirror (base class, version 1.0)

```
bend_angle_x:NX_FLOAT
bend_angle_y:NX_FLOAT
coating_material:NX_CHAR
coating_roughness:NX_FLOAT
description:NX_CHAR
even_layer_density:NX_FLOAT
even_layer_material:NX_CHAR
external_material:NX_CHAR
incident_angle:NX_FLOAT
interior_atmosphere:NX_CHAR
layer_thickness:NX_FLOAT
m_value:NX_FLOAT
odd_layer_density:NX_FLOAT
odd_layer_material:NX_CHAR
substrate_density:NX_FLOAT
substrate_material:NX_CHAR
substrate_roughness:NX_FLOAT
substrate_thickness:NX_FLOAT
type:NX_CHAR
reflectivity:NXdata
figure_data:NXdata
NXgeometry
shape:NXshape
```

- MIRROR:NXmirror →
_NXmirror.NX_tree_path NEXUSTREEPATH
_NXmirror.NX_id MIRROR
_NXmirror.NX_scan_id SCANID
_NXmirror.NX_diffn_id DIFFRNID
_NXmirror.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXmirror_MIRROR’ where MIRROR is the name of this group, typically ‘mirror’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- bend_angle_x:NX_FLOAT=BEND_ANGLE_X →
_NXmirror.bend_angle_x BEND_ANGLE_X
- bend_angle_y:NX_FLOAT=BEND_ANGLE_Y →
_NXmirror.bend_angle_y BEND_ANGLE_Y
- coating_material:NX_CHAR=COATING_MATERIAL →
_NXmirror.coating_material COATING_MATERIAL

- coating_roughness:NX_FLOAT=COATING_ROUGHNESS →
_NXmirror.coating_roughness COATING_ROUGHNESS
- description:NX_CHAR=DESCRIPTION →
_NXmirror.description DESCRIPTION
- even_layer_density:NX_FLOAT=EVEN_LAYER_DENSITY →
_NXmirror.even_layer_density EVEN_LAYER_DENSITY
- even_layer_material:NX_CHAR=EVEN_LAYER_MATERIAL →
_NXmirror.even_layer_material EVEN_LAYER_MATERIAL
- external_material:NX_CHAR=EXTERNAL_MATERIAL →
_NXmirror.external_material EXTERNAL_MATERIAL
- incident_angle:NX_FLOAT=INCIDENT_ANGLE →
_NXmirror.incident_angle INCIDENT_ANGLE
- interior_atmosphere:NX_CHAR=INTERIOR_ATMOSPHERE →
_NXmirror.interior_atmosphere INTERIOR_ATMOSPHERE
- layer_thickness:NX_FLOAT=LAYER_THICKNESS →
_NXmirror.layer_thickness LAYER_THICKNESS
- m_value:NX_FLOAT=M_VALUE →
_NXmirror.m_value M_VALUE
- odd_layer_density:NX_FLOAT=ODD_LAYER_DENSITY →
_NXmirror.odd_layer_density ODD_LAYER_DENSITY
- odd_layer_material:NX_CHAR=ODD_LAYER_MATERIAL →
_NXmirror.odd_layer_material ODD_LAYER_MATERIAL
- substrate_density:NX_FLOAT=SUBSTRATE_DENSITY →
_NXmirror.substrate_density SUBSTRATE_DENSITY
- substrate_material:NX_CHAR=SUBSTRATE_MATERIAL →
_NXmirror.substrate_material SUBSTRATE_MATERIAL
- substrate_roughness:NX_FLOAT=SUBSTRATE_ROUGHNESS →
_NXmirror.substrate_roughness SUBSTRATE_ROUGHNESS
- substrate_thickness:NX_FLOAT=SUBSTRATE_THICKNESS →
_NXmirror.substrate_thickness SUBSTRATE_THICKNESS
- type:NX_CHAR=TYPE →
_NXmirror.type TYPE
- reflectivity:NXdata →
_NXmirror.NXdata_id reflectivity
- figure_data:NXdata →
_NXmirror.NXdata_id figure_data
- geometry1:NXgeometry →
_NXmirror.NXgeometry_id geometry1
- shape:NXshape →
_NXmirror.NXshape_id shape

6.27 NXmoderator

NXmoderator (base class, version 1.0)

```
coupled:NX_BOOLEAN
coupling_material:NX_CHAR
distance:NX_FLOAT
poison_depth:NX_FLOAT
poison_material:NX_CHAR
temperature:NX_FLOAT
type:NX_CHAR
pulse_shape:NXdata
NXgeometry
temperature_log:NXlog
```

- MODERATOR:NXmoderator →
_NXmoderator.NX_tree_path NEXUSTREEPATH
_NXmoderator.NX_id MODERATOR
_NXmoderator.NX_scan_id SCANID
_NXmoderator.NX_diffn_id DIFFRNID
_NXmoderator.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXmoderator_MODERATOR’ where MODERATOR is the name of this group, typically ‘moderator’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- coupled:NX_BOOLEAN=COUPLED →
_NXmoderator.coupled COUPLED
- coupling_material:NX_CHAR=COUPLING_MATERIAL →
_NXmoderator.coupling_material COUPLING_MATERIAL
- distance:NX_FLOAT=DISTANCE →
_NXmoderator.distance DISTANCE
- poison_depth:NX_FLOAT=POISON_DEPTH →
_NXmoderator.poison_depth POISON_DEPTH
- poison_material:NX_CHAR=POISON_MATERIAL →
_NXmoderator.poison_material POISON_MATERIAL
- temperature:NX_FLOAT=TEMPERATURE →
_NXmoderator.temperature TEMPERATURE
- type:NX_CHAR=TYPE →
_NXmoderator.type TYPE
- pulse_shape:NXdata →
_NXmoderator.NXdata_id pulse_shape
- geometry1:NXgeometry →
_NXmoderator.NXgeometry_id geometry1
- temperature_log:NXlog →
_NXmoderator.NXlog_id temperature_log

6.28 NXmonitor

NXmonitor (base class, version 1.0)

```
count_time:NX_FLOAT
data:NX_NUMBER[n]
  @signal
  @axes
distance:NX_FLOAT
efficiency:NX_NUMBER[]
end_time:NX_DATE_TIME
integral:NX_NUMBER
mode:NX_CHAR
preset:NX_NUMBER
range:NX_FLOAT[2]
sampled_fraction:NX_FLOAT
start_time:NX_DATE_TIME
time_of_flight:NX_FLOAT[]
type:NX_CHAR
NXgeometry
integral_log:NXlog
```

- MONITOR:NXmonitor →
_NXmonitor.NX_tree_path NEXUSTREEPATH
_NXmonitor.NX_id MONITOR
_NXmonitor.NX_scan_id SCANID
_NXmonitor.NX_diffn_id DIFFRNID
_NXmonitor.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXmonitor_MONITOR’ where MONITOR is the name of this group, typically ‘monitor’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- count_time:NX_FLOAT=COUNT_TIME →
_NXmonitor.count_time COUNT_TIME
- data:NX_NUMBER[n]=DATA →
_NXmonitor.data DATA
- @signal=SIGNAL →
_NXmonitor.data__signal SIGNAL
- @axes=AXES →
_NXmonitor.data__axes AXES
- distance:NX_FLOAT=DISTANCE →
_NXmonitor.distance DISTANCE
- efficiency:NX_NUMBER[]=EFFICIENCY →
_NXmonitor.efficiency EFFICIENCY
- end_time:NX_DATE_TIME=END_TIME →
_NXmonitor.end_time END_TIME
- integral:NX_NUMBER=INTEGRAL →
_NXmonitor.integral INTEGRAL

- `mode:NX_CHAR=MODE →`
`_NXmonitor.mode MODE`
- `preset:NX_NUMBER=PRESET →`
`_NXmonitor.preset PRESET`
- `range:NX_FLOAT[2]=RANGE →`
`_NXmonitor.range RANGE`
- `sampled_fraction:NX_FLOAT=SAMPLED_FRACTION →`
`_NXmonitor.sampled_fraction SAMPLED_FRACTION`
- `start_time:NX_DATE_TIME=START_TIME →`
`_NXmonitor.start_time START_TIME`
- `time_of_flight:NX_FLOAT[]=TIME_OF_FLIGHT →`
`_NXmonitor.time_of_flight TIME_OF_FLIGHT`
- `type:NX_CHAR=TYPE →`
`_NXmonitor.type TYPE`
- `geometry1:NXgeometry →`
`_NXmonitor.NXgeometry_id geometry1`
- `integral_log:NXlog →`
`_NXmonitor.NXlog_id integral_log`

6.29 NXmonochromator

NXmonochromator (base class, version 1.0)

`energy:NX_FLOAT`
`energy_error:NX_FLOAT`
`wavelength:NX_FLOAT`
`wavelength_error:NX_FLOAT`
`NXcrystal`
`distribution:NXdata`
`geometry:NXgeometry`
`NXvelocity_selector`

- `MONOCHROMATOR:NXmonochromator →`
`_NXmonochromator.NX_tree_path NEXUSTREEPATH`
`_NXmonochromator.NX_id MONOCHROMATOR`
`_NXmonochromator.NX_scan_id SCANID`
`_NXmonochromator.NX_diffn_id DIFFRNID`
`_NXmonochromator.NX_entry_id ENTRYID`

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with `‘/NXmonochromator_MONOCHROMATOR’` where MONOCHROMATOR is the name of this group, typically `‘monochromator’`. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- `energy:NX_FLOAT=ENERGY →`
`_NXmonochromator.energy ENERGY`
- `energy_error:NX_FLOAT=ENERGY_ERROR →`
`_NXmonochromator.energy_error ENERGY_ERROR`

- wavelength:NX_FLOAT=WAVELENGTH →
_NXmonochromator.wavelength WAVELENGTH
- wavelength_error:NX_FLOAT=WAVELENGTH_ERROR →
_NXmonochromator.wavelength_error WAVELENGTH_ERROR
- crystal1:NXcrystal →
_NXmonochromator.NXcrystal_id crystal1
- distribution:NXdata →
_NXmonochromator.NXdata_id distribution
- geometry:NXgeometry →
_NXmonochromator.NXgeometry_id geometry
- velocity_selector1:NXvelocity_selector →
_NXmonochromator.NXvelocity_selector_id velocity_selector1

6.30 NXnote

NXnote (base class, version 1.0)

author:NX_CHAR
data:NX_BINARY
date:NX_DATE_TIME
description:NX_CHAR
file_name:NX_CHAR
type:NX_CHAR

- NOTE:NXnote →
_NXnote.NX_tree_path NEXUSTREEPATH
_NXnote.NX_id NOTE
_NXnote.NX_scan_id SCANID
_NXnote.NX_diffn_id DIFFRNID
_NXnote.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXnote_’NOTE’ where NOTE is the name of this group, typically ‘note’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- author:NX_CHAR=AUTHOR →
_NXnote.author AUTHOR
- data:NX_BINARY=DATA →
_NXnote.data DATA
- date:NX_DATE_TIME=DATE →
_NXnote.date DATE
- description:NX_CHAR=DESCRIPTION →
_NXnote.description DESCRIPTION
- file_name:NX_CHAR=FILE_NAME →
_NXnote.file_name FILE_NAME
- type:NX_CHAR=TYPE →
_NXnote.type TYPE

6.31 NXorientation

NXorientation (base class, version 1.0)

value:NX_FLOAT[numobj,6]

NXgeometry

- ORIENTATION:NXorientation →
_NXorientation.NX_tree_path NEXUSTREEPATH
_NXorientation.NX_id ORIENTATION
_NXorientation.NX_scan_id SCANID
_NXorientation.NX_diffn_id DIFFRNID
_NXorientation.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXorientation_ORIENTATION’ where ORIENTATION is the name of this group, typically ‘orientation’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- value:NX_FLOAT[numobj,6]=VALUE →
_NXorientation.value VALUE
- geometry1:NXgeometry →
_NXorientation.NXgeometry_id geometry1

6.32 NXparameters

NXparameters (base class, version 1.0)

term:NX_CHAR

@units

- PARAMETERS:NXparameters →
_NXparameters.NX_tree_path NEXUSTREEPATH
_NXparameters.NX_id PARAMETERS
_NXparameters.NX_scan_id SCANID
_NXparameters.NX_diffn_id DIFFRNID
_NXparameters.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXparameters_PARAMETERS’ where PARAMETERS is the name of this group, typically ‘parameters’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- term:NX_CHAR=TERM →
_NXparameters.term TERM
- @units=UNITS →
_NXparameters.term__units UNITS

6.33 NXpolarizer

NXpolarizer (base class, version 1.0)

composition:NX_CHAR

efficiency:NX_FLOAT

reflection:NX_INT[3]

type:NX_CHAR

- POLARIZER:NXpolarizer →
 NXpolarizer.NX_tree_path NEXUSTREEPATH
 NXpolarizer.NX_id POLARIZER
 NXpolarizer.NX_scan_id SCANID
 NXpolarizer.NX_diffn_id DIFFRNID
 NXpolarizer.NX_entry_id ENTRYID
 where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘‘/NXpolarizer_POLARIZER’’ where POLARIZER is the name of this group, typically ‘‘polarizer’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.
- composition:NX_CHAR=COMPOSITION →
 NXpolarizer.composition COMPOSITION
- efficiency:NX_FLOAT=EFFICIENCY →
 NXpolarizer.efficiency EFFICIENCY
- reflection:NX_INT[3]=REFLECTION →
 NXpolarizer.reflection REFLECTION
- type:NX_CHAR=TYPE →
 NXpolarizer.type TYPE

6.34 NXpositioner

NXpositioner (base class, version 1.0)

```
acceleration_time:NX_NUMBER
controller_record:NX_CHAR
description:NX_CHAR
name:NX_CHAR
raw_value:NX_NUMBER[n]
soft_limit_max:NX_NUMBER
soft_limit_min:NX_NUMBER
target_value:NX_NUMBER[n]
tolerance:NX_NUMBER[n]
value:NX_NUMBER[n]
velocity:NX_NUMBER
```

- POSITIONER:NXpositioner →
 NXpositioner.NX_tree_path NEXUSTREEPATH
 NXpositioner.NX_id POSITIONER
 NXpositioner.NX_scan_id SCANID
 NXpositioner.NX_diffn_id DIFFRNID
 NXpositioner.NX_entry_id ENTRYID
 where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘‘/NXpositioner_POSITIONER’’ where POSITIONER is the name of this group, typically ‘‘positioner’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.
- acceleration_time:NX_NUMBER=ACCELERATION_TIME →
 NXpositioner.acceleration_time ACCELERATION_TIME

- controller_record:NX_CHAR=CONTROLLER_RECORD →
_NXpositioner.controller_record CONTROLLER_RECORD
- description:NX_CHAR=DESCRIPTION →
_NXpositioner.description DESCRIPTION
- name:NX_CHAR=NAME →
_NXpositioner.name NAME
- raw_value:NX_NUMBER[n]=RAW_VALUE →
_NXpositioner.raw_value RAW_VALUE
- soft_limit_max:NX_NUMBER=SOFT_LIMIT_MAX →
_NXpositioner.soft_limit_max SOFT_LIMIT_MAX
- soft_limit_min:NX_NUMBER=SOFT_LIMIT_MIN →
_NXpositioner.soft_limit_min SOFT_LIMIT_MIN
- target_value:NX_NUMBER[n]=TARGET_VALUE →
_NXpositioner.target_value TARGET_VALUE
- tolerance:NX_NUMBER[n]=TOLERANCE →
_NXpositioner.tolerance TOLERANCE
- value:NX_NUMBER[n]=VALUE →
_NXpositioner.value VALUE
- velocity:NX_NUMBER=VELOCITY →
_NXpositioner.velocity VELOCITY

6.35 NXprocess

NXprocess (base class, version 1.0)

date:NX_DATE_TIME
program:NX_CHAR
version:NX_CHAR
NXnote

- PROCESS:NXprocess →
_NXprocess.NX_tree_path NEXUSTREEPATH
_NXprocess.NX_id PROCESS
_NXprocess.NX_scan_id SCANID
_NXprocess.NX_diffn_id DIFFRNID
_NXprocess.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXprocess_PROCESS’ where PROCESS is the name of this group, typically ‘process’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- date:NX_DATE_TIME=DATE →
_NXprocess.date DATE
- program:NX_CHAR=PROGRAM →
_NXprocess.program PROGRAM
- version:NX_CHAR=VERSION →
_NXprocess.version VERSION
- note1:NXnote →
_NXprocess.NXnote_id note1

6.36 NXroot

NXroot (base class, version 1.0)

```
@NX_class
@file_time
@file_name
@file_update_time
@NeXus_version
@HDF_version
@HDF5_Version
@XML_version
@creator
NXentry
```

- ROOT:NXroot →
_NXroot.NX_tree_path NEXUSTREEPATH
_NXroot.NX_id ROOT
_NXroot.NX_scan_id SCANID
_NXroot.NX_diffn_id DIFFRNID
_NXroot.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the group, ending with ‘/NXroot_ROOT’ where ROOT is the name of this group, typically ‘root’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- @NX_class=NX.CLASS →
_NXroot.NX_class_attribute__NX_class NX.CLASS
- @file_time=FILE.TIME →
_NXroot.NX_class_attribute__file_time FILE.TIME
- @file_name=FILE.NAME →
_NXroot.NX_class_attribute__file_name FILE.NAME
- @file_update_time=FILE.UPDATE.TIME →
_NXroot.NX_class_attribute__file_update_time FILE.UPDATE.TIME
- @NeXus_version=NEXUS.VERSION →
_NXroot.NX_class_attribute__NeXus_version NEXUS.VERSION
- @HDF_version=HDF.VERSION →
_NXroot.NX_class_attribute__HDF_version HDF.VERSION
- @HDF5_Version=HDF5.VERSION →
_NXroot.NX_class_attribute__HDF5.Version HDF5.VERSION
- @XML_version=XML.VERSION →
_NXroot.NX_class_attribute__XML_version XML.VERSION
- @creator=CREATOR →
_NXroot.NX_class_attribute__creator CREATOR
- entry1:NXentry →
_NXroot.NXentry_id entry1

6.37 NXsample

```
NXsample (base class, version 1.0)
  changer_position:NX_INT
  chemical_formula:NX_CHAR
  component:NX_CHAR
  concentration:NX_FLOAT[n_comp]
  density:NX_FLOAT[n_comp]
  description:NX_CHAR
  distance:NX_FLOAT
  electric_field:NX_FLOAT[n_eField]
    @direction
  external_DAC:NX_FLOAT
  magnetic_field:NX_FLOAT[n_mField]
    @direction
  mass:NX_FLOAT[n_comp]
  name:NX_CHAR
  orientation_matrix:NX_FLOAT[n_comp,3,3]
  path_length:NX_FLOAT
  path_length_window:NX_FLOAT
  preparation_date:NX_DATE_TIME
  pressure:NX_FLOAT[n_pField]
  relative_molecular_mass:NX_FLOAT[n_comp]
  rotation_angle:NX_FLOAT
  sample_component:NX_CHAR
  sample_orientation:NX_FLOAT[3]
  scattering_length_density:NX_FLOAT[n_comp]
  short_title:NX_CHAR
  situation:NX_CHAR
  stress_field:NX_FLOAT[n_sField]
    @direction
  temperature:NX_FLOAT[n_Temp]
  thickness:NX_FLOAT
  type:NX_CHAR
  unit_cell:NX_FLOAT[n_comp,6]
  unit_cell_class:NX_CHAR
  unit_cell_group:NX_CHAR
  unit_cell_volume:NX_FLOAT[n_comp]
  volume_fraction:NX_FLOAT[n_comp]
  x_translation:NX_FLOAT
NXbeam
  transmission:NXdata
  temperature_env:NXenvironment
  magnetic_field_env:NXenvironment
  geometry:NXgeometry
  temperature_log:NXlog
  magnetic_field_log:NXlog
  external_ADC:NXlog
```

- SAMPLE:NXsample →
_NXsample.NX_tree_path NEXUSTREEPATH


```

_NXsample.NX_id SAMPLE
_NXsample.NX_scan_id SCANID
_NXsample.NX_diffn_id DIFFRNID
_NXsample.NX_entry_id ENTRYID

```

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with ‘‘/NXsample_SAMPLE’’ where SAMPLE is the name of this NeXus class instance, typically ‘‘sample’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- changer_position:NX_INT=CHANGER_POSITION →
_NXsample.changer_position CHANGER_POSITION
- chemical_formula:NX_CHAR=CHEMICAL_FORMULA →
_NXsample.chemical_formula CHEMICAL_FORMULA
- component:NX_CHAR=COMPONENT →
_NXsample.component COMPONENT
- concentration:NX_FLOAT[n_comp]=CONCENTRATION →
_NXsample.concentration CONCENTRATION
- density:NX_FLOAT[n_comp]=DENSITY →
_NXsample.density DENSITY
- description:NX_CHAR=DESCRIPTION →
_NXsample.description DESCRIPTION
- distance:NX_FLOAT=DISTANCE →
_NXsample.distance DISTANCE
- electric_field:NX_FLOAT[n_eField]=ELECTRIC_FIELD →
_NXsample.electric_field ELECTRIC_FIELD
- @direction=DIRECTION →
_NXsample.electric_field_direction DIRECTION
- external_DAC:NX_FLOAT=EXTERNAL_DAC →
_NXsample.external_DAC EXTERNAL_DAC
- magnetic_field:NX_FLOAT[n_mField]=MAGNETIC_FIELD →
_NXsample.magnetic_field MAGNETIC_FIELD
- @direction=DIRECTION →
_NXsample.magnetic_field_direction DIRECTION
- mass:NX_FLOAT[n_comp]=MASS →
_NXsample.mass MASS
- name:NX_CHAR=NAME →
_NXsample.name NAME
- orientation_matrix:NX_FLOAT[n_comp,3,3]=ORIENTATION_MATRIX →
_NXsample.orientation_matrix ORIENTATION_MATRIX
- path_length:NX_FLOAT=PATH_LENGTH →
_NXsample.path_length PATH_LENGTH
- path_length_window:NX_FLOAT=PATH_LENGTH_WINDOW →
_NXsample.path_length_window PATH_LENGTH_WINDOW
- preparation_date:NX_DATE.TIME=PREPARATION_DATE →
_NXsample.preparation_date PREPARATION_DATE

- `pressure:NX_FLOAT[n_pField]=PRESSURE →`
`_NXsample.pressure PRESSURE`
- `relative_molecular_mass:NX_FLOAT[n_comp]=RELATIVE_MOLECULAR_MASS →`
`_NXsample.relative_molecular_mass RELATIVE_MOLECULAR_MASS`
- `rotation_angle:NX_FLOAT=ROTATION_ANGLE →`
`_NXsample.rotation_angle ROTATION_ANGLE`
- `sample_component:NX_CHAR=SAMPLE_COMPONENT →`
`_NXsample.sample_component SAMPLE_COMPONENT`
- `sample_orientation:NX_FLOAT[3]=SAMPLE_ORIENTATION →`
`_NXsample.sample_orientation SAMPLE_ORIENTATION`
- `scattering_length_density:NX_FLOAT[n_comp]=SCATTERING_LENGTH_DENSITY →`
`_NXsample.scattering_length_density SCATTERING_LENGTH_DENSITY`
- `short_title:NX_CHAR=SHORT_TITLE →`
`_NXsample.short_title SHORT_TITLE`
- `situation:NX_CHAR=SITUATION →`
`_NXsample.situation SITUATION`
- `stress_field:NX_FLOAT[n_sField]=STRESS_FIELD →`
`_NXsample.stress_field STRESS_FIELD`
- `@direction=DIRECTION →`
`_NXsample.stress_field_direction DIRECTION`
- `temperature:NX_FLOAT[n_Temp]=TEMPERATURE →`
`_NXsample.temperature TEMPERATURE`
- `thickness:NX_FLOAT=THICKNESS →`
`_NXsample.thickness THICKNESS`
- `type:NX_CHAR=TYPE →`
`_NXsample.type TYPE`
- `unit_cell:NX_FLOAT[n_comp,6]=UNIT_CELL →`
`_NXsample.unit_cell UNIT_CELL`
- `unit_cell_class:NX_CHAR=UNIT_CELL_CLASS →`
`_NXsample.unit_cell_class UNIT_CELL_CLASS`
- `unit_cell_group:NX_CHAR=UNIT_CELL_GROUP →`
`_NXsample.unit_cell_group UNIT_CELL_GROUP`
- `unit_cell_volume:NX_FLOAT[n_comp]=UNIT_CELL_VOLUME →`
`_NXsample.unit_cell_volume UNIT_CELL_VOLUME`
- `volume_fraction:NX_FLOAT[n_comp]=VOLUME_FRACTION →`
`_NXsample.volume_fraction VOLUME_FRACTION`
- `x_translation:NX_FLOAT=X_TRANSLATION →`
`_NXsample.x_translation X_TRANSLATION`
- `beam1:NXbeam →`
`_NXsample.NXbeam_id beam1`
- `transmission:NXdata →`
`_NXsample.NXdata_id transmission`
- `temperature_env:NXenvironment →`
`_NXsample.NXenvironment_id temperature.env`

- magnetic_field.env:NXenvironment →
_NXsample.NXenvironment_id magnetic_field.env
- geometry:NXgeometry →
_NXsample.NXgeometry_id geometry
- temperature_log:NXlog →
_NXsample.NXlog_id temperature_log
- magnetic_field_log:NXlog →
_NXsample.NXlog_id magnetic_field_log
- external_ADC:NXlog →
_NXsample.NXlog_id external_ADC

6.38 NXsensor

NXsensor (base class, version 1.0)

```
attached_to:NX_CHAR
external_field_brief:NX_CHAR
high_trip_value:NX_FLOAT
low_trip_value:NX_FLOAT
measurement:NX_CHAR
model:NX_CHAR
name:NX_CHAR
run_control:NX_BOOLEAN
short_name:NX_CHAR
type:NX_CHAR
value:NX_FLOAT[n]
value_deriv1:NX_FLOAT[]
value_deriv2:NX_FLOAT[]
geometry:NXgeometry
value_log:NXlog
value_deriv1_log:NXlog
value_deriv2_log:NXlog
external_field_full:NXorientation
```

- SENSOR:NXsensor →
_NXsensor.NX_tree_path NEXUSTREEPATH
_NXsensor.NX_id SENSOR
_NXsensor.NX_scan_id SCANID
_NXsensor.NX_diffn_id DIFFRNID
_NXsensor.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with ‘‘/NXsensor_SENSOR’’ where SENSOR is the name of this NeXus class instance, typically ‘‘sensor’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- attached_to:NX_CHAR=ATTACHED_TO →
_NXsensor.attached_to ATTACHED_TO
- external_field_brief:NX_CHAR=EXTERNAL_FIELD_BRIEF →
_NXsensor.external_field_brief EXTERNAL_FIELD_BRIEF

- `high_trip_value:NX_FLOAT=HIGH_TRIP_VALUE →`
`_NXsensor.high_trip_value HIGH_TRIP_VALUE`
- `low_trip_value:NX_FLOAT=LOW_TRIP_VALUE →`
`_NXsensor.low_trip_value LOW_TRIP_VALUE`
- `measurement:NX_CHAR=MEASUREMENT →`
`_NXsensor.measurement MEASUREMENT`
- `model:NX_CHAR=MODEL →`
`_NXsensor.model MODEL`
- `name:NX_CHAR=NAME →`
`_NXsensor.name NAME`
- `run_control:NX_BOOLEAN=RUN_CONTROL →`
`_NXsensor.run_control RUN_CONTROL`
- `short_name:NX_CHAR=SHORT_NAME →`
`_NXsensor.short_name SHORT_NAME`
- `type:NX_CHAR=TYPE →`
`_NXsensor.type TYPE`
- `value:NX_FLOAT[n]=VALUE →`
`_NXsensor.value VALUE`
- `value_deriv1:NX_FLOAT[]=VALUE_DERIV1 →`
`_NXsensor.value_deriv1 VALUE_DERIV1`
- `value_deriv2:NX_FLOAT[]=VALUE_DERIV2 →`
`_NXsensor.value_deriv2 VALUE_DERIV2`
- `geometry:NXgeometry →`
`_NXsensor.NXgeometry_id geometry`
- `value_log:NXlog →`
`_NXsensor.NXlog_id value_log`
- `value_deriv1_log:NXlog →`
`_NXsensor.NXlog_id value_deriv1_log`
- `value_deriv2_log:NXlog →`
`_NXsensor.NXlog_id value_deriv2_log`
- `external_field_full:NXorientation →`
`_NXsensor.NXorientation_id external_field_full`

6.39 NXshape

NXshape (base class, version 1.0)

`direction:NX_CHAR`

`shape:NX_CHAR`

`size:NX_FLOAT[numobj,nshapepar]`

- `SHAPE:NXshape →`
`_NXshape.NX_tree_path NEXUSTREEPATH`
`_NXshape.NX_id SHAPE`
`_NXshape.NX_scan_id SCANID`
`_NXshape.NX_diffn_id DIFFRNID`

_NXshape.NX_entry_id ENTRYID
 where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with ‘‘/NXshape__SHAPE’’ where SHAPE is the name of this NeXus class instance, typically ‘‘shape’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- direction:NX_CHAR=DIRECTION →
 _NXshape.direction DIRECTION
- shape:NX_CHAR=SHAPE →
 _NXshape.shape SHAPE
- size:NX_FLOAT[numobj,nshapepar]=SIZE →
 _NXshape.size SIZE

6.40 NXsource

NXsource (base class, version 1.0)

bunch_distance:NX_FLOAT
 bunch_length:NX_FLOAT
 current:NX_FLOAT
 distance:NX_FLOAT
 emittance_x:NX_FLOAT
 emittance_y:NX_FLOAT
 energy:NX_FLOAT
 flux:NX_FLOAT
 frequency:NX_FLOAT
 last_fill:NX_NUMBER
 @time
 mode:NX_CHAR
 name:NX_CHAR
 @short_name
 number_of_bunches:NX_INT
 period:NX_FLOAT
 power:NX_FLOAT
 probe:NX_CHAR
 pulse_width:NX_FLOAT
 sigma_x:NX_FLOAT
 sigma_y:NX_FLOAT
 target_material:NX_CHAR
 top_up:NX_BOOLEAN
 type:NX_CHAR
 voltage:NX_FLOAT
 bunch_pattern:NXdata
 title:NX_CHAR
 pulse_shape:NXdata
 distribution:NXdata
 geometry:NXgeometry
 notes:NXnote

- SOURCE:NXsource →
 _NXsource.NX_tree_path NEXUSTREEPATH

```

_NXsource.NX_id SOURCE
_NXsource.NX_scan_id SCANID
_NXsource.NX_diffn_id DIFFRNID
_NXsource.NX_entry_id ENTRYID

```

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with ‘‘/NXsource__SOURCE’’ where SOURCE is the name of this NeXus class instance, typically ‘‘source’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- bunch_distance:NX_FLOAT=BUNCH_DISTANCE →
_NXsource.bunch_distance BUNCH_DISTANCE
- bunch_length:NX_FLOAT=BUNCH_LENGTH →
_NXsource.bunch_length BUNCH_LENGTH
- current:NX_FLOAT=CURRENT →
_NXsource.current CURRENT
- distance:NX_FLOAT=DISTANCE →
_NXsource.distance DISTANCE
- emittance_x:NX_FLOAT=EMITTANCE_X →
_NXsource.emittance_x EMITTANCE_X
- emittance_y:NX_FLOAT=EMITTANCE_Y →
_NXsource.emittance_y EMITTANCE_Y
- energy:NX_FLOAT=ENERGY →
_NXsource.energy ENERGY
- flux:NX_FLOAT=FLUX →
_NXsource.flux FLUX
- frequency:NX_FLOAT=FREQUENCY →
_NXsource.frequency FREQUENCY
- last_fill:NX_NUMBER=LAST_FILL →
_NXsource.last_fill LAST_FILL
- @time=TIME →
_NXsource.last_fill__time TIME
- mode:NX_CHAR=MODE →
_NXsource.mode MODE
- name:NX_CHAR=NAME →
_NXsource.name NAME
- @short_name=SHORT_NAME →
_NXsource.name__short_name SHORT_NAME
- number_of_bunches:NX_INT=NUMBER_OF_BUNCHES →
_NXsource.number_of_bunches NUMBER_OF_BUNCHES
- period:NX_FLOAT=PERIOD →
_NXsource.period PERIOD
- power:NX_FLOAT=POWER →
_NXsource.power POWER
- probe:NX_CHAR=PROBE →
_NXsource.probe PROBE

- pulse_width:NX_FLOAT=PULSE_WIDTH →
_NXsource.pulse_width PULSE_WIDTH
- sigma_x:NX_FLOAT=SIGMA_X →
_NXsource.sigma_x SIGMA_X
- sigma_y:NX_FLOAT=SIGMA_Y →
_NXsource.sigma_y SIGMA_Y
- target_material:NX_CHAR=TARGET_MATERIAL →
_NXsource.target_material TARGET_MATERIAL
- top_up:NX_BOOLEAN=TOP_UP →
_NXsource.top_up TOP_UP
- type:NX_CHAR=TYPE →
_NXsource.type TYPE
- voltage:NX_FLOAT=VOLTAGE →
_NXsource.voltage VOLTAGE
- bunch_pattern:NXdata →
_NXsource.NXdata_id bunch_pattern
- title:NX_CHAR=TITLE →
_NXsource.title TITLE
- pulse_shape:NXdata →
_NXsource.NXdata_id pulse_shape
- distribution:NXdata →
_NXsource.NXdata_id distribution
- geometry:NXgeometry →
_NXsource.NXgeometry_id geometry
- notes:NXnote →
_NXsource.NXnote_id notes

6.41 NXsubentry

NXsubentry (base class, version 1.0)

```
@IDF_Version
collection_description:NX_CHAR
collection_identifier:NX_CHAR
collection_time:NX_FLOAT
definition:NX_CHAR
    @version
    @URL
definition_local:NX_CHAR
    @version
    @URL
duration:NX_INT
end_time:NX_DATE_TIME
entry_identifier:NX_CHAR
experiment_description:NX_CHAR
experiment_identifier:NX_CHAR
pre_sample_flightpath:NX_FLOAT
```

```

program_name:NX_CHAR
  @version
  @configuration
revision:NX_CHAR
  @comment
run_cycle:NX_CHAR
start_time:NX_DATE_TIME
title:NX_CHAR
NXcharacterization
NXdata
NXinstrument
NXmonitor
experiment_documentation:NXnote
notes:NXnote
thumbnail:NXnote
  @mime_type
NXprocess
NXsample
NXuser

```

- SUBENTRY:NXsubentry →

```

  _NXsubentry.NX_tree_path NEXUSTREEPATH
  _NXsubentry.NX_id SUBENTRY
  _NXsubentry.NX_scan_id SCANID
  _NXsubentry.NX_diffn_id DIFFRNID
  _NXsubentry.NX_entry_id ENTRYID

```

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with ‘‘/NXsubentry__SUBENTRY’’ where SUBENTRY is the name of this NeXus class instance, typically ‘‘subentry’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- @IDF_Version=IDF.VERSION →
_NXsubentry.NX_class_attribute__IDF_Version IDF.VERSION
- collection_description:NX_CHAR=COLLECTION_DESCRIPTION →
_NXsubentry.collection_description COLLECTION_DESCRIPTION
- collection_identfier:NX_CHAR=COLLECTION_IDENTIFIER →
_NXsubentry.collection_identfier COLLECTION_IDENTIFIER
- collection_time:NX_FLOAT=COLLECTION_TIME →
_NXsubentry.collection_time COLLECTION.TIME
- definition:NX_CHAR=DEFINITION →
_NXsubentry.definition DEFINITION
- @version=VERSION →
_NXsubentry.definition__version VERSION
- @URL=URL →
_NXsubentry.definition__URL URL
- definition_local:NX_CHAR=DEFINITION_LOCAL →
_NXsubentry.definition_local DEFINITION_LOCAL

- @version=VERSION →
_NXsubentry.definition.local__version VERSION
- @URL=URL →
_NXsubentry.definition.local__URL URL
- duration:NX_INT=DURATION →
_NXsubentry.duration DURATION
- end.time:NX_DATE_TIME=END_TIME →
_NXsubentry.end.time END.TIME
- entry_identifier:NX_CHAR=ENTRY_IDENTIFIER →
_NXsubentry.entry_identifier ENTRY_IDENTIFIER
- experiment_description:NX_CHAR=EXPERIMENT_DESCRIPTION →
_NXsubentry.experiment_description EXPERIMENT_DESCRIPTION
- experiment_identifier:NX_CHAR=EXPERIMENT_IDENTIFIER →
_NXsubentry.experiment_identifier EXPERIMENT_IDENTIFIER
- pre_sample_flightpath:NX_FLOAT=PRE_SAMPLE_FLIGHTPATH →
_NXsubentry.pre_sample_flightpath PRE.SAMPLE.FLIGHTPATH
- program_name:NX_CHAR=PROGRAM_NAME →
_NXsubentry.program_name PROGRAM_NAME
- @version=VERSION →
_NXsubentry.program_name__version VERSION
- @configuration=CONFIGURATION →
_NXsubentry.program_name__configuration CONFIGURATION
- revision:NX_CHAR=REVISION →
_NXsubentry.revision REVISION
- @comment=COMMENT →
_NXsubentry.revision__comment COMMENT
- run_cycle:NX_CHAR=RUN_CYCLE →
_NXsubentry.run_cycle RUN_CYCLE
- start.time:NX_DATE_TIME=START_TIME →
_NXsubentry.start.time START.TIME
- title:NX_CHAR=TITLE →
_NXsubentry.title TITLE
- characterization1:NXcharacterization →
_NXsubentry.NXcharacterization_id characterization1
- data1:NXdata →
_NXsubentry.NXdata_id data1
- instrument1:NXinstrument →
_NXsubentry.NXinstrument_id instrument1
- monitor1:NXmonitor →
_NXsubentry.NXmonitor_id monitor1
- experiment_documentation:NXnote →
_NXsubentry.NXnote_id experiment_documentation
- notes:NXnote →
_NXsubentry.NXnote_id notes

- thumbnail:NXnote →
_NXsubentry.NXnote_id thumbnail
- @mime_type=MIME.TYPE →
_NXsubentry.title_mime_type MIME.TYPE
- process1:NXprocess →
_NXsubentry.NXprocess_id process1
- sample1:NXsample →
_NXsubentry.NXsample_id sample1
- user1:NXuser →
_NXsubentry.NXuser_id user1

6.42 NXtransformations

NXtransformations (base class, version 1.0)

@default:NX_CHAR

AXISNAME:

@transformation_type: (optional) NX_CHAR

@vector: (required) NX_NUMBER[3]

@offset: (optional) NX_NUMBER[3]

@offset_units: (optional) NX_CHAR

@depends_on: (optional) NX_CHAR

@equipment_component: (optional) NX_CHAR

AXISNAME_end: (optional) NX_CHAR

AXISNAME_increment_set: (optional) NX_CHAR

- AXISNAME:NXtransformations →
_NXtransformations_transformation_type XFORMTYPE
_NXtransformations_vector XFORMVECTOR
_NXtransformations_offset XFORMOFFSET
_NXtransformations_offset_units XFORMOFFSTUN
_NXtransformations_depends_on XFORMDEPON
_NXtransformations_equipment_component XFORMCOMPON

for which, in the CIF version, AXISNAME XFORMTYPE is mapped to `_axis.id`, XFORMVECTOR and XFORMOFFSET need to be converted from NeXus McStas conventions to appropriate CBF/imgCIF conventions and then mapped to `_axis.vector[1]`, `_axis.vector[2]`, and `_axis.vector[3]`, and `_axis.offset[1]`, `_axis.offset[2]`, and `_axis.offset[3]`, XFORMDEPON is mapped to `_axis.depends_on`, and XFORMCOMPON is mapped to `_axis.equipment` and `_axis.equipment_component`.

6.43 NXtranslation

NXtranslation (base class, version 1.0)

distances:NX_FLOAT[numobj,3]

geometry:NXgeometry

- TRANSLATION:NXtranslation →
_NXtranslation.NX_tree_path NEXUSTREEPATH
_NXtranslation.NX_id TRANSLATION
_NXtranslation.NX_scan_id SCANID
_NXtranslation.NX_diffn_id DIFFRNID

`_NXtranslation.NX_entry_id ENTRYID`
 where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with `“/NXtranslation_TRANSLATION”` where TRANSLATION is the name of this NeXus class instance, typically `“translation”`.
 The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- `distances:NX_FLOAT[numobj,3]=DISTANCES →`
`_NXtranslation.distances DISTANCES`
- `geometry:NXgeometry →`
`_NXtranslation.NXgeometry_id geometry`

6.44 NXuser

NXuser (base class, version 1.0)

`address:NX_CHAR`
`affiliation:NX_CHAR`
`email:NX_CHAR`
`facility_user_id:NX_CHAR`
`fax_number:NX_CHAR`
`name:NX_CHAR`
`role:NX_CHAR`
`telephone_number:NX_CHAR`

- `USER:NXuser →`
`_NXuser.NX_tree_path NEXUSTREEPATH`
`_NXuser.NX_id USER`
`_NXuser.NX_scan_id SCANID`
`_NXuser.NX_diffn_id DIFFRNID`
`_NXuser.NX_entry_id ENTRYID`

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with `“/NXuser_USER”` where USER is the name of this NeXus class instance, typically `“user”`. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- `address:NX_CHAR=ADDRESS →`
`_NXuser.address ADDRESS`
- `affiliation:NX_CHAR=AFFILIATION →`
`_NXuser.affiliation AFFILIATION`
- `email:NX_CHAR=EMAIL →`
`_NXuser.email EMAIL`
- `facility_user_id:NX_CHAR=FACILITY_USER_ID →`
`_NXuser.facility_user_id FACILITY_USER_ID`
- `fax_number:NX_CHAR=FAX_NUMBER →`
`_NXuser.fax_number FAX_NUMBER`
- `name:NX_CHAR=NAME →`
`_NXuser.name NAME`
- `role:NX_CHAR=ROLE →`
`_NXuser.role ROLE`

- `telephone_number:NX_CHAR=TELEPHONE_NUMBER →`
`_NXuser.telephone_number TELEPHONE_NUMBER`

6.45 NXvelocity_selector

NXvelocity_selector (base class, version 1.0)

```
height:NX_FLOAT
length:NX_FLOAT
num:NX_INT
radius:NX_FLOAT
rotation_speed:NX_FLOAT
spwidth:NX_FLOAT
table:NX_FLOAT
twist:NX_FLOAT
type:NX_CHAR
wavelength:NX_FLOAT
wavelength_spread:NX_FLOAT
width:NX_FLOAT
geometry:NXgeometry
```

- `VELOCITY_SELECTOR:NXvelocity_selector →`
`_NXvelocity_selector.NX_tree_path NEXUSTREEPATH`
`_NXvelocity_selector.NX_id VELOCITY_SELECTOR`
`_NXvelocity_selector.NX_scan_id SCANID`
`_NXvelocity_selector.NX_diffn_id DIFFRNID`
`_NXvelocity_selector.NX_entry_id ENTRYID`

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with `“/NXvelocity_selector__VELOCITY_SELECTOR”` where VELOCITY_SELECTOR is the name of this NeXus class instance, typically `“velocity_selector”`. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- `height:NX_FLOAT=HEIGHT →`
`_NXvelocity_selector.height HEIGHT`
- `length:NX_FLOAT=LENGTH →`
`_NXvelocity_selector.length LENGTH`
- `num:NX_INT=NUM →`
`_NXvelocity_selector.num NUM`
- `radius:NX_FLOAT=RADIUS →`
`_NXvelocity_selector.radius RADIUS`
- `rotation_speed:NX_FLOAT=ROTATION_SPEED →`
`_NXvelocity_selector.rotation_speed ROTATION_SPEED`
- `spwidth:NX_FLOAT=SPWIDTH →`
`_NXvelocity_selector.spwidth SPWIDTH`
- `table:NX_FLOAT=TABLE →`
`_NXvelocity_selector.table TABLE`
- `twist:NX_FLOAT=TWIST →`
`_NXvelocity_selector.twist TWIST`

- `type:NX_CHAR=TYPE →`
`_NXvelocity_selector.type TYPE`
- `wavelength:NX_FLOAT=WAVELENGTH →`
`_NXvelocity_selector.wavelength WAVELENGTH`
- `wavelength_spread:NX_FLOAT=WAVELENGTH_SPREAD →`
`_NXvelocity_selector.wavelength_spread WAVELENGTH_SPREAD`
- `width:NX_FLOAT=WIDTH →`
`_NXvelocity_selector.width WIDTH`
- `geometry:NXgeometry →`
`_NXvelocity_selector.NXgeometry_id geometry`

6.46 NXxraylens

NXxraylens (base class, version 1.0)

```
aperture:NX_FLOAT
curvature:NX_FLOAT
cylindrical:NX_BOOLEAN
focus_type:NX_CHAR
gas:NX_CHAR
gas_pressure:NX_FLOAT
lens_geometry:NX_CHAR
lens_length:NX_FLOAT
lens_material:NX_CHAR
lens_thickness:NX_FLOAT
number_of_lenses:NX_INT
symmetric:NX_BOOLEAN
cylinder_orientation:NXnote
```

- `XRAYLENS:NXxraylens →`
`_NXxraylens.NX_tree_path NEXUSTREEPATH`
`_NXxraylens.NX_id XRAYLENS`
`_NXxraylens.NX_scan_id SCANID`
`_NXxraylens.NX_diffn_id DIFFRNID`
`_NXxraylens.NX_entry_id ENTRYID`

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with `‘/NXxraylens_XRAYLENS’` where XRAYLENS is the name of this NeXus class instance, typically `‘xraylens’`. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- `aperture:NX_FLOAT=APERTURE →`
`_NXxraylens.aperture APERTURE`
- `curvature:NX_FLOAT=CURVATURE →`
`_NXxraylens.curvature CURVATURE`
- `cylindrical:NX_BOOLEAN=CYLINDRICAL →`
`_NXxraylens.cylindrical CYLINDRICAL`
- `focus_type:NX_CHAR=FOCUS_TYPE →`
`_NXxraylens.focus_type FOCUS_TYPE`

- `gas:NX_CHAR=GAS →
_NXxraylens.gas GAS`
- `gas_pressure:NX_FLOAT=GAS_PRESSURE →
_NXxraylens.gas_pressure GAS_PRESSURE`
- `lens_geometry:NX_CHAR=LENS_GEOMETRY →
_NXxraylens.lens_geometry LENS_GEOMETRY`
- `lens_length:NX_FLOAT=LENS_LENGTH →
_NXxraylens.lens_length LENS_LENGTH`
- `lens_material:NX_CHAR=LENS_MATERIAL →
_NXxraylens.lens_material LENS_MATERIAL`
- `lens_thickness:NX_FLOAT=LENS_THICKNESS →
_NXxraylens.lens_thickness LENS_THICKNESS`
- `number_of_lenses:NX_INT=NUMBER_OF_LENSES →
_NXxraylens.number_of_lenses NUMBER_OF_LENSES`
- `symmetric:NX_BOOLEAN=SYMMETRIC →
_NXxraylens.symmetric SYMMETRIC`
- `cylinder_orientation:NXnote →
_NXxraylens.NXnote_id cylinder_orientation`

7 Proposed Pixel Array Detector Application Definitions

The following has been derived from the current Dectris Eiger test data and presentations and the Dectris web site:

https://www.dectris.com/nexus.html#main_head_navigation

In this version we have noted the multi-NXDATA use of `NXentry`, and note the extensions to `NXdectector` by including what Dectris has called `detectorSpecific:DetectorSpecific` `detectorSpecific:DECTRIS_detectorSpecific` as a recommended placeholders for such information in an update to the `NXdectector` base class. We have also called `detectorModule_xxx:DetectorModule` `detectorModule_xxx:DECTRIS_detector_module` and `detectorChip_xxx:DetectorChip` `detectorChip_xxx:DECTRIS_detector_chip`. These changes in class names have no impact on the relevant HDF5 paths.

7.1 NXentry

The base class `NXentry` currently has the following structure:

```

NXentry (base class, version 1.0)
  @IDF_Version
  collection_description:NX_CHAR
  collection_identifier:NX_CHAR
  collection_time:NX_FLOAT
  definition:NX_CHAR
    @version
    @URL
  definition_local:NX_CHAR
    @version
    @URL
  duration:NX_INT
  end_time:NX_DATE_TIME
  entry_identifier:NX_CHAR
  experiment_description:NX_CHAR
  experiment_identifier:NX_CHAR
  pre_sample_flightpath:NX_FLOAT
  program_name:NX_CHAR
    @version
    @configuration
  revision:NX_CHAR
    @comment
  run_cycle:NX_CHAR
  start_time:NX_DATE_TIME
  title:NX_CHAR
  NXcharacterization
  NXdata
  NXinstrument
  NXmonitor
  experiment_documentation:NXnote
  notes:NXnote
  thumbnail:NXnote
    @mime_type
  NXprocess
  NXsample
  NXsubentry
  NXuser

```

There is no conflict between the Dectris proposal and this class. We recommend that NIAC formally adopt the use of multiple NXdata NeXus class instances in NXentry and note that in the NXentry base class definition, not just imply that possibility in the NXdata base class definition as is now the case, so that the Dectris proposal

```

NXentry_pad (application definition, version 0.1) (overlays NXentry)
  data_000001:NXDATA
  data_000002:NXDATA
  ...
  data_nnnnnn:NXDATA
  instrument:NXinstrument
  ...

```

will not raise any questions. The CIF mapping of the Dectris proposal would then be:

- ENTRY:NXentry →
 - _NXentry.NX_tree_path NEXUSTREEPATH
 - _NXentry.NX_id ENTRY
 - _NXentry.NX_scan_id SCANID
 - _NXentry.NX_diffraction_id DIFFRACTIONID
 - _NXentry.NX_entry_id ENTRYID
 where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with ‘‘/NXentry_ENTRY’’ where ENTRY is the name of this NeXus class instance, typically ‘‘entry’’. The SCANID, DIFFRACTIONID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.
- data_000001:NXDATA →
 - _NXentry.NXDATA_id data_000001
- data_000002:NXDATA →
 - _NXentry.NXDATA_id data_000002
- ...
- data_nnnnnn:NXDATA →
 - _NXentry.NXDATA_id data_nnnnnn
- instrument:NXinstrument →
 - _NXentry.NXinstrument_id instrument

7.2 NXinstrument

There is no conflict between the Dectris proposal and this class. The Dectris use is

NXinstrument

dectector:NXdectector

which is a standard use and would map in CIF to

- INSTRUMENT:NXinstrument →
 - _NXinstrument.NX_tree_path NEXUSTREEPATH
 - _NXinstrument.NX_id INSTRUMENT
 - _NXinstrument.NX_scan_id SCANID
 - _NXinstrument.NX_diffraction_id DIFFRACTIONID
 - _NXinstrument.NX_entry_id ENTRYID
 where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with ‘‘/NXinstrument_INSTRUMENT’’ where INSTRUMENT is the name of this NeXus class instance, typically ‘‘instrument’’. The SCANID, DIFFRACTIONID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.
- dectector:NXdectector →
 - _NXinstrument.NXdectector_id dectector

7.3 NXdetector

There are minor conflicts between the Dectris proposal and this class. Dectris has used the detector_number field as a character string, rather than as a number. We have added the ‘Dectris_’ prefix for the character string version and changes the name to ‘Dectris_detector_serial_number’. The use of character string serial ‘numbers’ is a common practice for electronic equipment, so we recommend that a ‘detector_serial_number:NX_CHAR’ field be added to the NXdetector base class to carry this information. The field detector_distance is not currently in this class. We have added the ‘Dectris_’ prefix as a placeholder until a field with this information is added to the base class. Mark Koennecke has suggested ‘sample_detector_distance’ as the new field name, so the name used here is ‘Dectris_sample_detector_distance’

Dectris has added one subclass, which is likely to be needed for almost all pixel array detectors. However, at the moment, what is proposed is a Dectris-specific class definition, so rather than proposing the subclass immediately as a new NeXus base class, we propose the subclass with a ‘DECTRIS_’ prefix. The revised NXdetector base class would then be as follows. The fields explicitly used by DECTRIS are noted with [**]. Fields proposed by Dectris that are not currently in the NXdetector base class are prefixed with the ‘DECTRIS_’ prefix

```
NXdetector (base class, version 1.1)
  acquisition_mode:NX_CHAR                      [**]
  angular_calibration:NX_FLOAT[i,j]             [**]
  angular_calibration_applied:NX_BOOLEAN         [**]
  azimuthal_angle:NX_FLOAT[np,i,j]
  beam_center_x:NX_FLOAT                        [**]
  beam_center_y:NX_FLOAT                        [**]
  bit_depth_readout:NX_UINT                     [**]
  calibration_date:NX_DATE_TIME
  countrate_correction_applied:NX_BOOLEAN       [**]
  count_time:NX_NUMBER[np]                     [**]
  crate:NX_INT[i,j]
  @local_name
  data:NX_NUMBER[np,i,j,tof]
  @signal
  @axes
  @long_name
  @check_sum
  @link
  data_error:NX_NUMBER[np,i,j,tof]
  @units
  @link
  dead_time:NX_FLOAT[np,i,j]
  description:NX_CHAR                          [**]
  DECTRIS_sample_detector_distance:NX_FLOAT      [**]
  detection_gas_path:NX_FLOAT
  detector_number:NX_INT[i,j]
  DECTRIS_detector_serial_number:NX_CHAR        [**]
  detector_readout_time:NX_FLOAT                [**]
  DECTRIS_efficiency_correction_applied:NX_BOOLEAN [**]
  diameter:NX_FLOAT
  distance:NX_FLOAT[np,i,j]
  flatfield:NX_FLOAT[i,j]
```

```

flatfield_applied:NX_BOOLEAN          [**]
flatfield_error:NX_FLOAT[i,j]
frame_start_number:NX_INT
frame_time:NX_FLOAT[NP]              [**]
gain_setting:NX_CHAR                 [**]
gas_pressure:NX_FLOAT[i,j]
input:NX_INT[i,j]
    @local_name
layout:NX_CHAR
local_name:NX_CHAR
number_of_cycles:NX_INT              [**]
pixel_mask:NX_FLOAT[i,j]
pixel_mask_applied:NX_BOOLEAN        [**]
polar_angle:NX_FLOAT[np,i,j]
raw_time_of_flight:NX_INT[tof+1]
    @frequency
saturation_value:NX_INT
sensor_material:NX_CHAR              [**]
sensor_thickness:NX_FLOAT            [**]
sequence_number:NX_CHAR
slot:NX_INT[i,j]
    @local_name
solid_angle:NX_FLOAT[i,j]
threshold_energy:NX_FLOAT            [**]
time_of_flight:NX_FLOAT[tof+1]
    @axis
    @primary
    @long_name
    @link
trigger_dead_time:NX_FLOAT
trigger_delay_time:NX_FLOAT
type:NX_CHAR
DECTRIS_virtual_pixel_correction_applied:NX_BOOLEAN [**]
x_pixel_offset:NX_FLOAT[i,j]
    @axis
    @primary
    @long_name
    @link
x_pixel_size:NX_FLOAT[i,j]          [**]
y_pixel_offset:NX_FLOAT[i,j]
    @axis
    @primary
    @long_name
y_pixel_size:NX_FLOAT[i,j]          [**]
characterization:NXcharacterization
detectorSpecific:DECTRIS\_detector\_specific [**]
efficiency:NXdata
    efficiency:NX_FLOAT[i,j,k]
    real_time:NX_NUMBER[i,j,k]
    wavelength:NX_FLOAT[i,j,k]

```

```

geometry:NXgeometry
calibration_method:NXnote
data_file:NXnote

    adding

DECTRIS_sample_detector_distance:NX_FLOAT
DECTRIS_detector_serial_number:NX_CHAR
DECTRIS_efficiency_correction_applied:NX_BOOLEAN
DECTRIS_virtual_pixel_correction_applied:NX_BOOLEAN

    and giving explicit suggested names to the NXcharacterization and NXgeometry instances.

NXdetector_detriss (application definition, version 0.1)
    (overlays NXDetector)
NXdetector
    acquisition_mode:NX_CHAR
    angular_calibration_applied:NX_BOOLEAN
    beam_center_x:NX_FLOAT
        @units
    beam_center_y:NX_FLOAT
        @units
    bit_depth_readout:NX_UINT
    count_time:NX_FLOAT[np]
        @units
    countrate_correction_applied:NX_BOOLEAN
    description:NX_CHAR
    DECTRIS_sample_detector_distance:NX_FLOAT
        @units
    DECTRIS_detector_serial_number:NX_CHAR
    detector_number:NX_CHAR
    detector_readout_time:NX_FLOAT[np]
        @units
    detectorSpecific:DECTRIS_detector_specific
    DECTRIS_efficiency_correction_applied:NX_BOOLEAN
    flatfield_correction_applied:NX_BOOL
    frame_time:NX_FLOAT[np]
        @units
    gain_setting:NX_CHAR
    number_of_cycles:NX_UINT
    pixel_mask_applied:NX_BOOL
    sensor_material:NX_STRING
    sensor_thickness:NX_FLOAT
        @units
    threshold_energy:NX_FLOAT
        @units
    DECTRIS_virtual_pixel_correction_applied:NX_BOOL
    x_pixel_size:NX_FLOAT
        @units
    y_pixel_size:NX_FLOAT
        @units

```

- DETECTOR:NXdetector →
 _NXdetector.NX_tree_path NEXUSTREEPATH
 _NXdetector.NX_id DETECTOR
 _NXdetector.NX_scan_id SCANID
 _NXdetector.NX_diffraction_id DIFFRACTIONID
 _NXdetector.NX_entry_id ENTRYID
 where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with ‘‘/NXdetector_DETECTOR’’ where DETECTOR is the name of this NeXus class instance, typically ‘‘detector’’. The SCANID, DIFFRACTIONID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.
- acquisition_mode:NX_CHAR=ACQUISITION_MODE →
 _NXdetector.acquisition_mode ACQUISITION_MODE
- angular_calibration_applied:NX_BOOLEAN=ANGULAR_CALIBRATION_APPLIED →
 _NXdetector.angular_calibration_applied ANGULAR_CALIBRATION_APPLIED
- beam_center_x:NX_FLOAT=BEAM_CENTER_X →
 _NXdetector.beam_center_x BEAM_CENTER_X
- @units=UNITS →
 _NXdetector.beam_center_x__units UNITS
- beam_center_y:NX_FLOAT=BEAM_CENTER_Y →
 _NXdetector.beam_center_y BEAM_CENTER_Y
- @units=UNITS →
 _NXdetector.beam_center_y__units UNITS
- bit_depth_readout:NX_UINT=BIT_DEPTH_READOUT →
 _NXdetector.bit_depth_readout BIT_DEPTH_READOUT
- count_time:NX_FLOAT[np]=COUNT_TIME →
 _NXdetector.count_time COUNT_TIME
- @units=UNITS →
 _NXdetector.count_time__units UNITS
- countrate_correction_applied:NX_BOOLEAN=COUNTRATE_CORRECTION_APPLIED →
 _NXdetector.countrate_correction_applied COUNTRATE_CORRECTION_APPLIED
- description:NX_CHAR=DESCRIPTION →
 _NXdetector.description DESCRIPTION
- detector_number:NX_CHAR=DETECTOR_NUMBER →
 _NXdetector.detector_number DETECTOR_NUMBER
- detectorSpecific:DECTRIS_detector_specific →
 _NXdetector.DECTRIS_detector_specific_id detectorSpecific
- DETRIS_detector_number:NX_CHAR=DETECTOR_NUMBER →
 _NXdetector.DECTRIS_detector_number DETECTOR_NUMBER
- detector_readout_time:NX_FLOAT[np]=DETECTOR_READOUT_TIME →
 _NXdetector.detector_readout_time DETECTOR_READOUT_TIME
- @units=UNITS →
 _NXdetector.detector_readout_time__units UNITS
- efficiency_correction_applied:NX_BOOL=EFFICIENCY_CORRECTION_APPLIED →
 _NXdetector.efficiency_correction_applied EFFICIENCY_CORRECTION_APPLIED

- flatfield_correction_applied:NX_BOOL=FLATFIELD.CORRECTION.APPLIED →
_NXdetector.flatfield_correction_applied FLATFIELD.CORRECTION.APPLIED
- frame_time:NX_FLOAT[np]=FRAME.TIME →
_NXdetector.frame_time FRAME.TIME
- @units=UNITS →
_NXdetector.frame_time__units UNITS
- gain_setting:NX_CHAR=GAIN.SETTING →
_NXdetector.gain_setting GAIN.SETTING
- number_of_cycles:NX_UINT=NUMBER.OF.CYCLES →
_NXdetector.number_of_cycles NUMBER.OF.CYCLES
- pixel_mask_applied:NX_BOOL=PIXEL.MASK.APPLIED →
_NXdetector.pixel_mask_applied PIXEL.MASK.APPLIED
- sensor_material:NX_STRING=SENSOR.MATERIAL →
_NXdetector.sensor_material SENSOR.MATERIAL
- sensor_thickness:NX_FLOAT=SENSOR.THICKNESS →
_NXdetector.sensor_thickness SENSOR.THICKNESS
- @units=UNITS →
_NXdetector.sensor_thickness__units UNITS
- threshold_energy:NX_FLOAT=THRESHOLD.ENERGY →
_NXdetector.threshold_energy THRESHOLD.ENERGY
- @units=UNITS →
_NXdetector.threshold_energy__units UNITS
- DECTRIS.virtual_pixel_correction_applied:NX_BOOL=VIRTUAL.PIXEL.CORRECTION.APPLIED →
_NXdetector.DECTRIS.virtual_pixel_correction_applied VIRTUAL.PIXEL.CORRECTION.APPLIED
- x_pixel_size:NX_FLOAT=X.PIXEL.SIZE →
_NXdetector.x_pixel_size X.PIXEL.SIZE
- @units=UNITS →
_NXdetector.x_pixel_size__units UNITS
- y_pixel_size:NX_FLOAT=Y.PIXEL.SIZE →
_NXdetector.y_pixel_size Y.PIXEL.SIZE
- @units=UNITS →
_NXdetector.y_pixel_size__units UNITS

7.4 detectorSpecific:DECTRIS_detector_specific

These are some additions added by Dectris, that they designated as a new detectorSpecific:DetectorSpecific NeXus class instance.

```
DECTRIS\_detector\_specific (application definition, version 0.1)
  countrate_correction_bunch_mode:NX_CHAR
  countrate_correction_count_cutoff:NX_UINT
  countrate_correction_lookup_table:NX_FLOAT[1000000]
  data_collection_date:NX_CHAR
  detectorModule_000:DECTRIS_detector_module
  detectorModule_001:DECTRIS_detector_module
```

[illegible]

```

detectorModule_052:DECTRIS_detector_module
detectorModule_053:DECTRIS_detector_module
detectorModule_054:DECTRIS_detector_module
detectorModule_055:DECTRIS_detector_module
detectorModule_056:DECTRIS_detector_module
detectorModule_057:DECTRIS_detector_module
detectorModule_058:DECTRIS_detector_module
detectorModule_059:DECTRIS_detector_module
detector_origin:NX_FLOAT
    @depends_on
    @transformation
    @units
    @vector
transformLabToDetector:NX_FLOAT
    @rotation:NX_FLOAT[9]
    @translation:NX_FLOAT[3]
    @units
flat field:NX_FLOAT[number of x pixels,number of y pixels]
mode_register:NX_UINT
nimages:NX_UINT
number_of_excluded_pixels:NX_UINT
photon_energy:NX_FLOAT
pixel_mask:NX_UINT[number of x pixels,number of y pixels]
readout_mode:NX_CHAR
software_version:NX_CHAR
sub_image_exposure_time:NX_FLOAT
summation_mode:NX_CHAR
summation_nimages:NX_UINT
trigger_mode:NX_CHAR
x_pixels_in_detector:NX_UINT
y_pixels_in_detector:NX_UINT

```

- DETECTOR_SPECIFIC:DECTRIS_detector_specific →
 _DECTRIS_detector_specific.NX_tree_path NEXUSTREEPATH
 _DECTRIS_detector_specific.NX_id DETECTOR_SPECIFIC
 _DECTRIS_detector_specific.NX_scan_id SCANID
 _DECTRIS_detector_specific.NX_diffraction_id DIFFRACTIONID
 _DECTRIS_detector_specific.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with ‘‘/NXcollection-__DETECTOR_SPECIFIC’’ where _DETECTOR_SPECIFIC is the name of this NeXus class instance, typically ‘‘collection’’. The SCANID, DIFFRACTIONID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- countrate_correction_bunch_mode:NX_CHAR=COUNTRATE_CORRECTION_BUNCH_MODE →
 _DECTRIS_detector_specific.countrate_correction_bunch_mode
 COUNTRATE_CORRECTION_BUNCH_MODE
- countrate_correction_count_cutoff:NX_UINT=
 COUNTRATE_CORRECTION_COUNT_CUTOFF →
 _DECTRIS_detector_specific.countrate_correction_count_cutoff
 COUNTRATE_CORRECTION_COUNT_CUTOFF

- `countrate_correction_lookup_table:NX_FLOAT[1000000]`
`=COUNTRATE_CORRECTION_LOOKUP_TABLE →`
`_DECTRIS_detector_specific.countrate_correction_lookup_table`
`COUNTRATE_CORRECTION_LOOKUP_TABLE`
- `data_collection_date:NX_CHAR=DATA_COLLECTION_DATE →`
`_DECTRIS_detector_specific.data_collection_date DATA_COLLECTION_DATE`
- `detectorModule.000:DECTRIS_detector_module`
`detectorModule.001:DECTRIS_detector_module`
`...`
`detectorModule.059:DECTRIS_detector_module →`
`_DECTRIS_detector_specific.DECTRIS_detector_module_id`
`["detectorModule.000","detectorModule.001",...,"detectorModule.059"]`
- `detector_origin:NX_FLOAT64=DETECTOR_ORIGIN →`
`_DECTRIS_detector_specific.detector_origin DETECTOR_ORIGIN`
- `@depends_on=DEPENDS_ON →`
`_DECTRIS_detector_specific.detector_origin__depends_on DEPENDS_ON`
- `@transformation=TRANSFORMATION →`
`_DECTRIS_detector_specific.detector_origin__transformation TRANSFORMATION`
- `@units=UNITS →`
`_DECTRIS_detector_specific.detector_origin__units UNITS`
- `@vector=VECTOR →`
`_DECTRIS_detector_specific.detector_origin__vector VECTOR`
- `transformLabToDetector:NX_FLOAT64=TRANSFORMLABTODETECTOR →`
`_DECTRIS_detector_specific.transformLabToDetector TRANSFORMLABTODETECTOR`
- `@rotation=ROTATION →`
`_DECTRIS_detector_specific.transformLabToDetector__rotation ROTATION`
- `@translation=TRANSLATION →`
`_DECTRIS_detector_specific.transformLabToDetector__translation TRANSLATION`
- `@units=UNITS →`
`_DECTRIS_detector_specific.transformLabToDetector__units UNITS`
- `flatfield:NX_FLOAT[numberofxpixels,numberofypixels]=FLATFIELD →`
`_DECTRIS_detector_specific.flatfield FLATFIELD`
- `mode_register:NX_UINT=MODE_REGISTER →`
`_DECTRIS_detector_specific.mode_register MODE_REGISTER`
- `nimages:NX_UINT=NIMAGES →`
`_DECTRIS_detector_specific.nimages NIMAGES`
- `number_of_excluded_pixels:NX_UINT=NUMBER_OF_EXCLUDED_PIXELS →`
`_DECTRIS_detector_specific.number_of_excluded_pixels NUMBER_OF_EXCLUDED_PIXELS`
- `photon_energy:NX_FLOAT=PHOTON_ENERGY →`
`_DECTRIS_detector_specific.photon_energy PHOTON_ENERGY`
- `pixel_mask:NX_UINT[numberofxpixels,numberofypixels]=PIXEL_MASK →`
`_DECTRIS_detector_specific.pixel_mask PIXEL_MASK`
- `readout_mode:NX_CHAR=READOUT_MODE →`
`_DECTRIS_detector_specific.readout_mode READOUT_MODE`

- software_version:NX_CHAR=SOFTWARE_VERSION →
_DECTRIS_detector_specific.software_version SOFTWARE_VERSION
- sub_image_exposure_time:NX_FLOAT=SUB_IMAGE_EXPOSURE_TIME →
_DECTRIS_detector_specific.sub_image_exposure_time SUB_IMAGE_EXPOSURE_TIME
- summation_mode:NX_CHAR=SUMMATION_MODE →
_DECTRIS_detector_specific.summation_mode SUMMATION_MODE
- summation_nimages:NX_UINT=SUMMATION_NIMAGES →
_DECTRIS_detector_specific.summation_nimages SUMMATION_NIMAGES
- trigger_mode:NX_CHAR=TRIGGER_MODE →
_DECTRIS_detector_specific.trigger_mode TRIGGER_MODE
- x_pixels_in_detector:NX_UINT=X_PIXELS_IN_DETECTOR →
_DECTRIS_detector_specific.x_pixels_in_detector X_PIXELS_IN_DETECTOR
- y_pixels_in_detector:NX_UINT=Y_PIXELS_IN_DETECTOR →
_DECTRIS_detector_specific.y_pixels_in_detector Y_PIXELS_IN_DETECTOR

The specific number of modules may vary.

7.5 detectorModule_nnn:DECTRIS_detector_module

DECTRIS_detector_module (application definition, version 0.1)

```

dac_names:NX_CHAR[6]
dac_values:NX_UINT[6]
data_origin:NX_UINT[2]
data_size:NX_UINT[2]
detectorChip_000:DECTRIS_detector_chip
detectorChip_001:DECTRIS_detector_chip
detectorChip_002:DECTRIS_detector_chip
detectorChip_003:DECTRIS_detector_chip
detectorChip_004:DECTRIS_detector_chip
detectorChip_005:DECTRIS_detector_chip
detectorChip_006:DECTRIS_detector_chip
detectorChip_007:DECTRIS_detector_chip
detectorChip_008:DECTRIS_detector_chip
detectorChip_009:DECTRIS_detector_chip
detectorChip_010:DECTRIS_detector_chip
detectorChip_011:DECTRIS_detector_chip
detectorChip_012:DECTRIS_detector_chip
detectorChip_013:DECTRIS_detector_chip
detectorChip_014:DECTRIS_detector_chip
detectorChip_015:DECTRIS_detector_chip
fast_pixel_direction:NX_FLOAT64
  @depends_on
  @transformation
  @units
  @vector
firmware_version:NX_CHAR
module_offset:NX_FLOAT64
  @depends_on
  @transformation

```

```

    @units
    @vector
nbits:NX_UINT
nchips:NX_UINT
readout_frequency:NX_FLOAT
    @units
region_of_interest:NX_UINT[4]
slow_pixel_direction:NX_FLOAT64
    @depends_on
    @transformation
    @units
    @vector
x_pixels_in_module:NX_UINT
y_pixels_in_module:NX_UINT

```

- DETECTORMODULE:DECTRIS_detector_module →
 _DECTRIS_detector_module.NX_tree_path NEXUSTREEPATH
 _DECTRIS_detector_module.NX_id DETECTORMODULE
 _DECTRIS_detector_module.NX_scan_id SCANID
 _DECTRIS_detector_module.NX_diffn_id DIFFRNID
 _DECTRIS_detector_module.NX_entry_id ENTRYID
 where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under
 score and, finally, the name of the NeXus class instance, ending with ‘‘/DECTRIS_detector_module-
 _DETECTORMODULE’’ where DETECTORMODULE is the name of this NeXus class instance, typically
 ‘‘detectorModule_nnn’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when
 multiple scans, *etc.* are aggregated in the same CBF.
- dac_names:NX_CHAR[6]=DAC_NAMES →
 _DECTRIS_detector_module.dac_names DAC_NAMES
- dac_values:NX_UINT[7]=DAC_VALUES →
 _DECTRIS_detector_module.dac_values DAC_VALUES
- data_origin:NX_UINT[2]=DATA_ORIGIN →
 _DECTRIS_detector_module.data_origin DATA_ORIGIN
- data_size:NX_UINT[2]=DATA_SIZE →
 _DECTRIS_detector_module.data_size DATA_SIZE
- detectorChip_000:DECTRIS_detector_chip
 detectorChip_001:DECTRIS_detector_chip
 ...
 detectorChip_015:DECTRIS_detector_chip →
 _DECTRIS_detector_module.NXcollection_id
 ["detectorChip_000","detectorChip_001",...,"detectorChip_015"]
- fast_pixel_direction:NX_FLOAT64=FAST_PIXEL_DIRECTION →
 _DECTRIS_detector_module.fast_pixel_direction FAST_PIXEL_DIRECTION
- @depends_on=DEPENDS_ON →
 _DECTRIS_detector_module.fast_pixel_direction__depends_on DEPENDS_ON
- @transformation=TRANSFORMATION →
 _DECTRIS_detector_module.fast_pixel_direction__transformation TRANSFORMATION
- @units=UNITS →
 _DECTRIS_detector_module.fast_pixel_direction__units UNITS

- @vector=VECTOR →
_DECTRIS_detector_module.fast_pixel_direction__vector VECTOR
- firmware_version:NX_CHAR=FIRMWARE_VERSION →
_DECTRIS_detector_module.firmware_version FIRMWARE_VERSION
- module_offset:NX_FLOAT64=MODULE_OFFSET →
_DECTRIS_detector_module.module_offset MODULE_OFFSET
- @depends_on=DEPENDS_ON →
_DECTRIS_detector_module.module_offset__depends_on DEPENDS_ON
- @transformation=TRANSFORMATION →
_DECTRIS_detector_module.module_offset__transformation TRANSFORMATION
- @units=UNITS →
_DECTRIS_detector_module.module_offset__units UNITS
- @vector=VECTOR →
_DECTRIS_detector_module.module_offset__vector VECTOR
- nbits:NX_UINT=NBITS →
_DECTRIS_detector_module.nbits NBITS
- nchips:NX_UINT=NCHIPS →
_DECTRIS_detector_module.nchips NCHIPS
- readout_frequency:NX_FLOAT=READOUT_FREQUENCY →
_DECTRIS_detector_module.readout_frequency READOUT_FREQUENCY
- @units=UNITS →
_DECTRIS_detector_module.readout_frequency__units UNITS
- region_of_interest:NX_UINT[4]=REGION_OF_INTEREST →
_DECTRIS_detector_module.region_of_interest REGION_OF_INTEREST
- slow_pixel_direction:NX_FLOAT64=SLOW_PIXEL_DIRECTION →
_DECTRIS_detector_module.slow_pixel_direction SLOW_PIXEL_DIRECTION
- @depends_on=DEPENDS_ON →
_DECTRIS_detector_module.slow_pixel_direction__depends_on DEPENDS_ON
- @transformation=TRANSFORMATION →
_DECTRIS_detector_module.slow_pixel_direction__transformation TRANSFORMATION
- @units=UNITS →
_DECTRIS_detector_module.slow_pixel_direction__units UNITS
- @vector=VECTOR →
_DECTRIS_detector_module.slow_pixel_direction__vector VECTOR
- x_pixels_in_module:NX_UINT=X_PIXELS_IN_MODULE →
_DECTRIS_detector_module.x_pixels_in_module X_PIXELS_IN_MODULE
- y_pixels_in_module:NX_UINT=Y_PIXELS_IN_MODULE →
_DECTRIS_detector_module.y_pixels_in_module Y_PIXELS_IN_MODULE

The specific number of chips may vary.

7.6 detectorChip_nn:DECTRIS_detector_chip

```
DECTRIS_detector_chip
  chip_type:NX_CHAR
  x_pixels_in_chip:NX_UINT
  x_position:NX_UINT
  y_pixels_in_chip:NX_UINT
  y_position:NX_UINT
```

- DETECTORCHIP:DECTRIS_detector_chip →
_DECTRIS_detector_chip.NX_tree_path NEXUSTREEPATH
_DECTRIS_detector_chip.NX_id DETECTORCHIP
_DECTRIS_detector_chip.NX_scan_id SCANID
_DECTRIS_detector_chip.NX_diffn_id DIFFRNID
_DECTRIS_detector_chip.NX_entry_id ENTRYID

where components of NEXUSTREEPATH are composed of the relevant NeXus class, a double under score and, finally, the name of the NeXus class instance, ending with ‘‘/DECTRIS_detector_chip-_DETECTORCHIP’’ where DETECTORCHIP is the name of this NeXus class instance, typically ‘‘detectorChip_nnn’’. The SCANID, DIFFRNID and ENTRYID are optional keys for use when multiple scans, *etc.* are aggregated in the same CBF.

- chip_type:NX_CHAR=CHIP_TYPE →
_NXcollection.chip_type CHIP_TYPE
- x_pixels_in_chip:NX_UINT=X_PIXELS_IN_CHIP →
_NXcollection.x_pixels_in_chip X_PIXELS_IN_CHIP
- x_position:NX_UINT=X_POSITION →
_NXcollection.x_position X_POSITION
- y_pixels_in_chip:NX_UINT=Y_PIXELS_IN_CHIP →
_NXcollection.y_pixels_in_chip Y_PIXELS_IN_CHIP
- y_position:NX_UINT=Y_POSITION →
_NXcollection.y_position Y_POSITION

7.7 Consolidated Dectris Eiger Application Definition

In addition to the above classes, the Dectris application definition draws on NXsample for the rotation_angle_step and NXmonochromator for the wavelength. The Dectris Eiger NeXus format as proposed by Dectris with the proposed changes in this document marked with ‘‘**[...]**’’ is as follows. However, rotation_angle_step is not in the NeXus base class dictionary. Therefore, we flag it with the DETRIS prefix until it is formally adopted as part of the NXsample base class.

```
NXdectris_eiger (application definition, version 0.1)
  (overlays NXentry)
  entry:NXentry
    data_000001:NXDATA
      data:NXINT[np_000001,number of x pixels,number of y pixels]
      @image_nr_low
      @image_nr_high
    data_000002:NXDATA
      data:NXINT[np_000001,number of x pixels,number of y pixels]
      @image_nr_low
```

```

    @image_nr_high
...
data_nnnnnn:NXDATA
    data:NXINT[np_nnnnnn,number of x pixels,number of y pixels]
        @image_nr_low
        @image_nr_high
instrument:NXinstrument
    detector:NXdetector
        acquisition_mode:NX_CHAR
        angular_calibration_applied:NX_FLOAT32[number of x pixels,number of y pixels]
        beam_center_x:NX_FLOAT
            @units
        beam_center_y:NX_FLOAT
            @units
        bit_depth_readout:NX_UINT
        countrate_correction_applied:NX_BOOLEAN
        count_time:NX_FLOAT[np]
            @units
        description:NX_CHAR
        detector_number:NX_CHAR
        detectorSpecific:DetectorSpecific                **[detectorSpecific:DECTRIS_detector_specific]**
            countrate_correction_bunch_mode:NX_CHAR
            countrate_correction_count_cutoff:NX_UINT
            countrate_correction_lookup_table:NX_FLOAT[1000000]
            data_collection_date:NX_CHAR
            detectorModule_000:detectorModule            **[detectorModule_000:DECTRIS_detector_module]**
                dac_names:NX_CHAR[6]
                dac_values:NX_UINT[7]
                data_origin:NX_UINT[2]
                data_size:NX_UINT[2]
                detectorChip_000:DetectorChip            **[detectorChip_000:DECTRIS_detector_chip]**
                    chip_type:NX_CHAR
                    x_pixels_in_chip:NX_UINT
                    x_position:NX_UINT
                    y_pixels_in_chip:NX_UINT
                    y_position:NX_UINT
                detectorChip_001:DetectorChip            **[detectorChip_001:DECTRIS_detector_chip]**
                    chip_type:NX_CHAR
                    x_pixels_in_chip:NX_UINT
                    x_position:NX_UINT
                    y_pixels_in_chip:NX_UINT
                    y_position:NX_UINT
            ...
            detectorChip_015:DetectorChip                **[detectorChip_015:DECTRIS_detector_chip]**
                chip_type:NX_CHAR
                x_pixels_in_chip:NX_UINT
                x_position:NX_UINT
                y_pixels_in_chip:NX_UINT
                y_position:NX_UINT
            fast_pixel_direction:NX_FLOAT

```

```

        @transformation
        @vector:NX_FLOAT[3]
        @units
        @depends_on
firmware_version:NX_CHAR
module_offset:NX_FLOAT
        @transformation
        @vector:NX_FLOAT[3]
        @units
        @depends_on
nbits:NX_UINT
nchips:NX_UINT
readout_frequency:NX_FLOAT
        @units
region_of_interest:NX_UINT[4]
slow_pixel_direction:NX_FLOAT
        @transformation
        @vector:NX_FLOAT[3]
        @units
        @depends_on
x_pixels_in_module:NX_UINT
y_pixels_in_module:NX_UINT
detectorModule_001:detectorModule
dac_names:NX_CHAR[6]
dac_values:NX_UINT[7]
data_origin:NX_UINT[2]
data_size:NX_UINT[2]
detectorChip_000:DetectorChip
        chip_type:NX_CHAR
        x_pixels_in_chip:NX_UINT
        x_position:NX_UINT
        y_pixels_in_chip:NX_UINT
        y_position:NX_UINT
detectorChip_001:DetectorChip
        chip_type:NX_CHAR
        x_pixels_in_chip:NX_UINT
        x_position:NX_UINT
        y_pixels_in_chip:NX_UINT
        y_position:NX_UINT
...
detectorChip_015:DetectorChip
        chip_type:NX_CHAR
        x_pixels_in_chip:NX_UINT
        x_position:NX_UINT
        y_pixels_in_chip:NX_UINT
        y_position:NX_UINT
fast_pixel_direction:NX_FLOAT
        @transformation
        @vector:NX_FLOAT[3]
        @units

```

****[detectorModule_001:DECTRIS_detector_module]****

****[detectorChip_000:DECTRIS_detector_chip]****

****[detectorChip_001:DECTRIS_detector_chip]****

****[detectorChip_015:DECTRIS_detector_chip]****

```

        @depends_on
        firmware_version:NX_CHAR
        module_offset:NX_FLOAT
        @transformation
        @vector:NX_FLOAT[3]
        @units
        @depends_on
        nbits:NX_UINT
        nchips:NX_UINT
        readout_frequency:NX_FLOAT
        @units
        region_of_interest:NX_UINT[4]
        slow_pixel_direction:NX_FLOAT
        @transformation
        @vector:NX_FLOAT[3]
        @units
        @depends_on
        x_pixels_in_module:NX_UINT
        y_pixels_in_module:NX_UINT
        ...
detectorModule_059:detectorModule                                **[detectorModule_059:DECTRIS_detector_module]**
    dac_names:NX_CHAR[6]
    dac_values:NX_UINT[7]
    data_origin:NX_UINT[2]
    data_size:NX_UINT[2]
    detectorChip_000:DetectorChip                                **[detectorChip_000:DECTRIS_detector_chip]**
        chip_type:NX_CHAR
        x_pixels_in_chip:NX_UINT
        x_position:NX_UINT
        y_pixels_in_chip:NX_UINT
        y_position:NX_UINT
    detectorChip_001:DetectorChip                                **[detectorChip_001:DECTRIS_detector_chip]**
        chip_type:NX_CHAR
        x_pixels_in_chip:NX_UINT
        x_position:NX_UINT
        y_pixels_in_chip:NX_UINT
        y_position:NX_UINT
    ...
    detectorChip_015:DetectorChip                                **[detectorChip_015:DECTRIS_detector_chip]**
        chip_type:NX_CHAR
        x_pixels_in_chip:NX_UINT
        x_position:NX_UINT
        y_pixels_in_chip:NX_UINT
        y_position:NX_UINT
    fast_pixel_direction:NX_FLOAT
        @transformation
        @vector:NX_FLOAT[3]
        @units
        @depends_on
        firmware_version:NX_CHAR

```

```

module_offset:NX_FLOAT
    @transformation
    @vector:NX_FLOAT[3]
    @units
    @depends_on
nbits:NX_UINT
nchips:NX_UINT
readout_frequency:NX_FLOAT
    @units
region_of_interest:NX_UINT[4]
slow_pixel_direction:NX_FLOAT
    @transformation
    @vector:NX_FLOAT[3]
    @units
    @depends_on
x_pixels_in_module:NX_UINT
y_pixels_in_module:NX_UINT
detector_origin:NX_FLOAT
    @transformation
    @vector:NX_FLOAT[3]
    @units
    @depends_on
transformLabToDetector:NX_FLOAT
    @rotation:NX_FLOAT[9]
    @translation:NX_FLOAT[3]
    @units
flatfield:NX_FLOAT[number of x pixels,number of y pixels]
mode_register:NX_UINT
nimages:NX_UINT
number_of_excluded_pixels:NX_UINT
photon_energy:NX_FLOAT
pixel_mask:NX_UINT[number of x pixels,number of y pixels]
readout_mode:NX_CHAR
software_version:NX_CHAR
sub_image_exposure_time:NX_FLOAT
summation_mode:NX_CHAR
summation_nimages:NX_UINT
trigger_mode:NX_CHAR
x_pixels_in_detector:NX_UINT
y_pixels_in_detector:NX_UINT
detector_number:NX_CHAR
detector_readout_time:NX_FLOAT[np]
    @units
efficiency_correction_applied:NX_BOOL
flatfield_correction_applied:NX_BOOL
frame_time:NX_FLOAT[np]
    @units
gain_setting:NX_CHAR
number_of_cycles:NX_UINT
pixel_mask_applied:NX_BOOL

```



```

sensor_material:NX_STRING
sensor_thickness:NX_FLOAT
    @units
threshold_energy:NX_FLOAT
    @units
virtual_pixel_correction_applied:NX_BOOLEAN
x_pixel_size:NX_FLOAT
    @units
y_pixel_size:NX_FLOAT
    @units
sample:NXsample
    rotation_angle_step:NX_FLOAT[np]                **[DECTRIS_rotation_angle_step]**
    @units
monochromator:NXmonochromator
    wavelength:NX_FLOAT32
    @units

```

References

- [1] Herbert J Bernstein, Andreas Förster, Asmit Bhowmick, Aaron S Brewster, Sandor Brockhauser, Luca Gelisio, David R Hall, Filip Leonarski, Valerio Mariani, Gianluca Santoni, Clemens Vonnrhein, and Graeme Winter. Gold standard for macromolecular crystallography diffraction data. *IUCrJ*, 7(5):784 -- 792, 2020.
- [2] HJ Bernstein and AP Hammersley. *International tables for crystallography, volume G: Definition and exchange of crystallographic data*, chapter 2.3. Specification of the Crystallographic Binary File (CBF/imgCIF), page 37. Springer, 2006.
- [3] HJ Bernstein, JM Sloan, G Winter, and TS Richter. Managing crystallographic data in facilities using integrated cif, hdf5 and nexus. *Work*, 25:29, 2013.
- [4] J Drenth. chapter Introduction to basic crystallography, pages 45 -- 63. Wiley Online Library, 2012.
- [5] Mark Könnecke, Frederick A Akeroyd, Herbert J Bernstein, Aaron S Brewster, Stuart I. Campbell, Björn Clausen, Stephen Cottrell, Jens Uwe Hoffmann, Pete R. Jemian, David Männicke, Raymond Osborn, Peter F. Peterson, Tobias Richter, Jiro Suzuki, Benjamin Watts Watts, Eugen Wintersberger, and Joachim Wuttke Wuttke. The nexus data format. *J. Appl. Cryst.*, 48(1):301 -- 305, 2015.
- [6] Andrew GW Leslie and Harry R Powell. Processing diffraction data with mosflm. In *Evolving Methods for Macromolecular Crystallography*, volume 245 of *NATO Science Series*, pages 41 -- 51. Springer, 2007.
- [7] Protein Data Bank. Protein data bank contents guide and dictionary. Format specification, Brookhaven National Laboratory, 1992. Available from: <https://www.wwpdb.org/documentation/file-format>.
- [8] George H Stout and Lyle H Jensen. *X-ray structure determination: a practical guide*. John Wiley & Sons, 1989.
- [9] John D. Westbrook, Jasmine Y. Young, Chenghua Shao, Zukang Feng, Vladimir Guranovic, Catherine L. Lawson, Brinda Vallat, Paul D. Adams, John M. Berrisford, Gerard Bricogne, Robbie P. Diederichs, Kay and Joosten, Peter Keller, Nigel W. Moriarty, Oleg V. Sobolev, Sameer Velankar, Clemens Vonnrhein, David G. Waterman, Genji Kurisu,

Helen M. Berman, Stephen K. Burley, and Ezra Peisach.
foundational semantic tools for structural biology.
2022.

Pdbx/mmCIF ecosystem:
J. Mol. Bio., 434(11):167599,